



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 6/19/2025

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 06/18/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 06/19/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB136193

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q2353-01	BC271315-1-1	1	1.00	1.00	2.00	2.00	100.0	PILC
Q2353-02	BC271315-1-2	2	1.00	1.00	2.00	2.00	100.0	PILC
Q2354-01	3321	3	1.15	10.12	11.27	10.35	90.9	
Q2354-03	3323	4	1.00	1.00	2.00	2.00	100.0	oil sample
Q2355-01	3897	5	1.00	1.00	2.00	2.00	100.0	DEBRIS
Q2355-03	20071	6	1.14	10.46	11.6	10.42	88.7	
Q2355-04	20071-E2	7	1.14	10.56	11.7	10.54	89.0	
Q2356-01	Y2308-0051-1-1	8	1.00	1.00	2.00	2.00	100.0	PILC
Q2356-02	Y2308-0051-1-2	9	1.00	1.00	2.00	2.00	100.0	PILC
Q2356-03	348-9010-1-1	10	1.00	1.00	2.00	2.00	100.0	PILC
Q2356-04	348-9010-1-2	11	1.00	1.00	2.00	2.00	100.0	PILC
Q2357-01	TR-06-06182025	12	1.18	10.45	11.63	10.66	90.7	
Q2357-02	TR-06-06182025-E2	13	1.12	10.86	11.98	11.71	97.5	
Q2358-01	NB-07-06182025	14	1.19	10.06	11.25	10.00	87.6	
Q2358-02	NB-07-06182025-E2	15	1.16	10.57	11.73	10.85	91.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$