



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 7/29/2025

OVENTEMP IN Celsius(°C): 108
Time IN: 16:45
In Date: 07/28/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:14
Out Date: 07/29/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID-OVEN

QC:LB136637

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
Q2709-01	7211935-295	9	1.12	10.77	11.89	11.53	96.7	
Q2710-01	JC-03-07282025	1	1.18	10.25	11.43	10.85	94.3	
Q2710-02	JC-03-07282025-E2	2	1.19	10.44	11.63	10.77	91.8	
Q2710-04	JC-03-07282025-E3	3	1.18	10.52	11.7	10.7	90.5	
Q2711-01	PCB15-33	4	1.00	1.00	2.00	2.00	100.0	Caulk sample
Q2711-02	PCB15-31	5	1.00	1.00	2.00	2.00	100.0	Caulk sample
Q2711-03	PCB15-28	6	1.00	1.00	2.00	2.00	100.0	Caulk sample
Q2711-04	PCB14-3	7	1.00	1.00	2.00	2.00	100.0	Caulk sample
Q2711-05	PCB14-5	8	1.00	1.00	2.00	2.00	100.0	Caulk sample

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$