



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 7/30/2025

OVENTEMP IN Celsius(°C): 108
Time IN: 16:45
In Date: 07/29/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:14
Out Date: 07/30/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID-OVEN

QC:LB136647

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q2720-01	M3A	1	1.00	1.00	2.00	2.00	100.0	PILC
Q2720-02	M3B	2	1.00	1.00	2.00	2.00	100.0	PILC
Q2720-03	M4A	3	1.00	1.00	2.00	2.00	100.0	PILC
Q2721-01	FDH025P-1-1	4	1.00	1.00	2.00	2.00	100.0	PILC
Q2721-02	FDH025P-1-2	5	1.00	1.00	2.00	2.00	100.0	PILC
Q2721-03	BC278720-1-1	6	1.00	1.00	2.00	2.00	100.0	PILC
Q2721-04	BC278720-1-2	7	1.00	1.00	2.00	2.00	100.0	PILC
Q2722-01	LAW-25-0107	8	1.00	1.00	2.00	2.00	100.0	oil sample
Q2722-03	LAW-25-0108	9	1.00	1.00	2.00	2.00	100.0	oil sample
Q2723-01	L11A	10	1.00	1.00	2.00	2.00	100.0	PILC
Q2723-02	L11B	11	1.00	1.00	2.00	2.00	100.0	PILC
Q2723-03	L12A	12	1.00	1.00	2.00	2.00	100.0	PILC
Q2723-04	L12B	13	1.00	1.00	2.00	2.00	100.0	PILC
Q2723-05	L13A	14	1.00	1.00	2.00	2.00	100.0	PILC

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$