



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 8/5/2025

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 08/04/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:14
Out Date: 08/05/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID-OVEN

QC:LB136687

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q2752-03	PIPE-LINE-LIQUIDS	1	1.00	1.00	2.00	2.00	100.0	oil sample
Q2757-01	WASTE	2	1.14	9.99	11.13	9.28	81.5	
Q2757-02	VOC	3	1.13	10.85	11.98	9.96	81.4	
Q2757-03	1	4	1.19	10.50	11.69	9.83	82.3	
Q2757-04	2	5	1.18	10.81	11.99	10.1	82.5	
Q2757-05	3	6	1.18	10.81	11.99	9.98	81.4	
Q2757-06	4	7	1.13	10.59	11.72	9.66	80.5	
Q2757-07	5	8	1.13	10.86	11.99	9.95	81.2	
Q2758-01	WASTE	9	1.12	10.57	11.69	10.94	92.9	
Q2759-01	SOLID-DRUMS	10	1.15	10.84	11.99	5.75	42.4	
Q2760-01	NB-07-08042025	11	1.13	10.76	11.89	11.36	95.1	
Q2760-02	NB-07-08042025-E2	12	1.13	10.37	11.5	10.44	89.8	
Q2763-01	TP-2	13	1.16	10.02	11.18	9.36	81.8	
Q2763-02	TP-2 EPH	14	1.18	10.12	11.3	9.75	84.7	
Q2763-03	TP-2 VOC	15	1.18	10.48	11.66	9.45	78.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$