



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 9/3/2025

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 09/02/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:33
Out Date: 09/03/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID-OVEN

QC:LB137037

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q2999-01	HIH071C-1-1	8	1.00	1.00	2.00	2.00	100.0	PILC
Q2999-02	HIH071C-1-2	9	1.00	1.00	2.00	2.00	100.0	PILC
Q2999-03	HIH638W-1-1	10	1.00	1.00	2.00	2.00	100.0	PILC
Q2999-04	HIH638W-1-2	11	1.00	1.00	2.00	2.00	100.0	PILC
Q2999-05	HEH766W-1-1	12	1.00	1.00	2.00	2.00	100.0	PILC
Q2999-06	HEH766W-1-2	13	1.00	1.00	2.00	2.00	100.0	PILC
Q2999-07	HEH766W-2-1	14	1.00	1.00	2.00	2.00	100.0	PILC
Q2999-08	HEH766W-2-2	15	1.00	1.00	2.00	2.00	100.0	PILC
Q3001-01	S-1	1	1.15	10.07	11.22	10.59	93.7	
Q3001-02	S-2	2	1.19	10.47	11.66	11.1	94.7	
Q3001-03	S-3	3	1.19	10.27	11.46	10.85	94.1	
Q3001-04	S-4	4	1.18	10.00	11.18	10.19	90.1	
Q3001-05	S-5	5	1.19	10.43	11.62	10.7	91.2	
Q3001-06	S-6	6	1.18	10.29	11.47	10.56	91.2	
Q3001-06D	S-6DUP	7	1.14	10.36	11.5	10.6	91.3	
Q3002-01	FDH119M-1-1	16	1.00	1.00	2.00	2.00	100.0	PILC
Q3002-02	FDH119M-1-2	17	1.00	1.00	2.00	2.00	100.0	PILC
Q3002-03	BC271332-1-1	18	1.00	1.00	2.00	2.00	100.0	PILC
Q3002-04	BC271332-1-2	19	1.00	1.00	2.00	2.00	100.0	PILC
Q3003-01	VNJ-231	20	1.14	10.45	11.59	11.14	95.7	
Q3003-02	VNJ-231	21	1.15	10.38	11.53	11.01	95.0	
Q3004-02	Drum Speedy Dry/RAGS	22	1.16	10.64	11.8	10.95	92.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$