



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 9/11/2025

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 09/10/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:37
Out Date: 09/11/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID-OVEN

QC:LB137131

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q3064-01	905-A	1	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
Q3064-02	905-B-FRONT	2	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
Q3064-03	905-B-BACK	3	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
Q3065-01	VNJ-204	4	1.18	10.85	12.03	10.24	83.5	
Q3065-02	VNJ-204-E2	5	1.16	10.61	11.77	10.15	84.7	
Q3066-01	354	6	1.13	10.53	11.66	11.18	95.4	
Q3066-02	354-EPH-2	7	1.14	10.16	11.3	10.8	95.1	
Q3066-03	VNJ-242	8	1.13	10.84	11.97	11.05	91.5	
Q3066-04	VNJ-242-EPH-2	9	1.16	10.56	11.72	10.71	90.4	
Q3068-01	AMBRIC-TECH-9-9	10	1.18	10.42	11.6	10.45	89.0	
Q3072-01	30-GARFIELD-PL-COMP	11	1.14	10.30	11.44	10.62	92.0	
Q3072-01D	30-GARFIELD-PL-COMPDUP	12	1.16	10.30	11.46	10.65	92.1	
Q3073-01	BEL-25-0022	13	1.00	1.00	2.00	2.00	100.0	oil sample
Q3073-02	BEL-25-0023	14	1.00	1.00	2.00	2.00	100.0	oily-debris
Q3074-01	MOO-25-0141	15	1.00	1.00	2.00	2.00	100.0	oily-debris
Q3075-01	MOO-25-0254	16	1.12	10.76	11.88	9.93	81.9	
Q3075-02	MOO-25-0254-0256	17	1.19	10.70	11.89	9.77	80.2	
Q3075-04	MOO-25-0258	18	1.00	1.00	2.00	2.00	100.0	debris
Q3075-06	MOO-25-0259	19	1.00	1.00	2.00	2.00	100.0	debris
Q3075-08	MOO-25-0260	20	1.13	10.83	11.96	11.56	96.3	
Q3075-09	MOO-25-0260	21	1.13	10.83	11.96	11.56	96.3	
Q3075-10	MOO-25-02563	22	1.14	10.51	11.65	10.36	87.7	
Q3075-11	MOO-25-02561-0263	23	1.15	10.35	11.5	10.22	87.6	
Q3076-01	L18A	24	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-02	L18B	25	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-03	L19A	26	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-04	L20A	27	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-05	L20B	28	1.00	1.00	2.00	2.00	100.0	PILC



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 9/11/2025

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 09/10/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:37
Out Date: 09/11/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID-OVEN

QC:LB137131

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q3076-06	L21A	29	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-07	L21B	30	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-08	L22A	31	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-09	L22B	32	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-10	245F75-1-1	33	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-11	245F75-1-2	34	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-12	JEB978E-1-1	35	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-13	JEB978E-1-2	36	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-14	JEB978E-2-1	37	1.00	1.00	2.00	2.00	100.0	PILC
Q3076-15	JEB978E-2-2	38	1.00	1.00	2.00	2.00	100.0	PILC
Q3077-01	LAW-25-0137	39	1.00	1.00	2.00	2.00	100.0	debris
Q3077-02	LAW-25-0137	40	1.16	10.46	11.62	10.2	86.4	
Q3077-03	LAW-25-0139	41	1.16	10.46	11.62	10.2	86.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$