



PERCENT SOLID

Supervisor: rubina
Analyst: jignesh
Date: 11/6/2025

OVENTEMP IN Celsius(°C): 104
Time IN: 17:30
In Date: 11/05/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:37
Out Date: 11/06/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID-OVEN

QC:LB137782

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q3544-01	ONYX-1	1	1.13	10.65	11.78	10.54	88.4	
Q3544-02	ONYX-1-EPH	2	1.18	10.40	11.58	10.28	87.5	
Q3544-03	ONYX-1-VOC	3	1.16	10.61	11.77	10.51	88.1	
Q3545-01	PCB-110325KC-01 Hard Doorframe Caulk	15	1.00	1.00	2.00	2.00	100.0	caulk sample
Q3546-01	TP-COMP	4	1.13	10.53	11.66	10.66	90.5	
Q3546-02	TP-1-EPH-1	5	1.13	10.48	11.61	10.38	88.3	
Q3546-03	TP-1-EPH-2	6	1.16	10.50	11.66	9.95	83.7	
Q3546-04	TP-1-EPH-3	7	1.18	10.18	11.36	10.37	90.3	
Q3546-05	TP-1-EPH-4	8	1.17	10.60	11.77	10.5	88.0	
Q3546-06	TP-2-EPH-5	9	1.15	10.76	11.91	11.1	92.5	
Q3546-07	TP-2-EPH-6	10	1.15	10.52	11.67	10.86	92.3	
Q3546-08	TP-2-EPH-7	11	1.14	10.32	11.46	10.93	94.9	
Q3546-09	TP-2-EPH-8	12	1.19	10.23	11.42	10.75	93.5	
Q3549-01	RT4922	13	1.11	10.45	11.56	1.76	6.2	GEL Matrix
Q3549-02	RT4922	14	1.16	10.64	11.8	3.29	20.0	GEL Matrix

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$