



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 12/3/2025

OVENTEMP IN Celsius(°C): 106
Time IN: 17:00
In Date: 12/02/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:25
Out Date: 12/03/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID-OVEN

QC:LB138082

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q3751-01	120125HU01-P	1	1.00	1.00	2.00	2.00	100.0	Caulk sample
Q3751-02	120125HU02-P	2	1.00	1.00	2.00	2.00	100.0	Caulk sample
Q3753-01	MOO-25-345-356	3	1.11	10.53	11.64	11.3	96.8	
Q3753-02	MOO-25-334-337	4	1.15	10.17	11.32	10.52	92.1	
Q3753-02D	MOO-25-334-337DUP	5	1.19	10.15	11.34	10.55	92.2	
Q3753-03	MOO-25-334-337	6	1.15	10.17	11.32	10.52	92.1	
Q3753-04	AR520-0002	7	1.15	10.61	11.76	10.59	89.0	
Q3753-05	AR520-0002	8	1.15	10.61	11.76	10.59	89.0	
Q3754-01	L73A	9	1.00	1.00	2.00	2.00	100.0	PILC
Q3754-02	L73B	10	1.00	1.00	2.00	2.00	100.0	PILC
Q3754-03	L74A	11	1.00	1.00	2.00	2.00	100.0	PILC
Q3754-04	L74B	12	1.00	1.00	2.00	2.00	100.0	PILC
Q3754-05	L75A	13	1.00	1.00	2.00	2.00	100.0	PILC
Q3754-06	L75B	14	1.00	1.00	2.00	2.00	100.0	PILC

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$