

Clear table

| ident      | concentrati<br>on F- | concentrati<br>on CL- | concentrati<br>on NO2 | concentrati<br>on BR- | concentrati<br>on NO3 | concentrati<br>on HPO4 | concentrati<br>on SO4 | file name    | date time      |
|------------|----------------------|-----------------------|-----------------------|-----------------------|-----------------------|------------------------|-----------------------|--------------|----------------|
| STD1       | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | _2016-12-19_ | 12/19/16 9:23  |
| STD2       | 0.3080               | 0.3830                | 0.3860                | 0.9820                | 0.3310                | 0.2850                 | 2.4720                | _2016-12-19_ | 12/19/16 9:52  |
| STD3       | 0.4910               | 0.7700                | 0.7620                | 2.3260                | 0.6160                | 1.1570                 | 4.2010                | _2016-12-19_ | 12/19/16 10:21 |
| STD4       | 0.8470               | 1.2660                | 1.2740                | 4.1810                | 0.9710                | 2.0570                 | 6.4890                | _2016-12-19_ | 12/19/16 10:50 |
| STD5       | 1.9400               | 2.8970                | 2.9010                | 9.7340                | 2.1420                | 4.8210                 | 14.0750               | _2016-12-19_ | 12/19/16 11:19 |
| STD6       | 3.8850               | 5.7000                | 5.8050                | 19.4350               | 4.2510                | 9.9100                 | 27.9370               | _2016-12-19_ | 12/19/16 11:48 |
| STD7       | 4.6150               | 7.1670                | 7.3290                | 24.2570               | 5.3470                | 11.2610                | 35.3110               | _2016-12-19_ | 12/19/16 12:16 |
| STD8       | 10.2150              | 15.2680               | 15.4550               | 50.5850               | 11.5410               | 24.4880                | 76.7650               | _2016-12-19_ | 12/19/16 12:45 |
| ICV        | 4.6490               | 7.2820                | 7.3520                | 24.5240               | 6.0510                | 11.5010                | 36.2640               | _2016-12-19_ | 12/19/16 14:14 |
| ICB        | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | _2016-12-19_ | 12/19/16 14:43 |
| CCV        | 5.1550               | 7.2630                | 7.4700                | 24.7900               | 5.4600                | 12.5390                | 36.3570               | _2017-01-04_ | 1/4/17 17:08   |
| CCB        | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | _2017-01-04_ | 1/4/17 17:36   |
| LB85241BLW | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | _2017-01-04_ | 1/4/17 18:05   |
| LB85241BSW | 4.9980               | 7.4400                | 7.6010                | 24.7790               | 6.0750                | 12.7200                | 36.8510               | _2017-01-04_ | 1/4/17 18:34   |
| I1002-08   | 0.1320               | 0.1580                | 0.1890                | 0.6430                | 0.1330                | 0.1510                 | 1.0200                | _2017-01-04_ | 1/4/17 19:03   |
| I1002-09   | 0.1310               | 0.1270                | 0.2000                | 0.6570                | 0.1290                | 0.1310                 | 1.0170                | _2017-01-04_ | 1/4/17 19:32   |
| I1002-10   | 0.1270               | 0.1330                | 0.1920                | 0.6380                | 0.1310                | 0.1060                 | 1.0360                | _2017-01-04_ | 1/4/17 20:01   |
| I1002-11   | 0.1340               | 0.1290                | 0.2040                | 0.6470                | 0.1330                | 0.1270                 | 1.0140                | _2017-01-04_ | 1/4/17 20:30   |
| I1002-12   | 0.1480               | 0.1520                | 0.1860                | 0.6370                | 0.1360                | 0.1040                 | 1.0140                | _2017-01-04_ | 1/4/17 20:59   |
| I1002-13   | 0.1370               | 0.1260                | 0.1900                | 0.6060                | 0.1250                | 0.1120                 | 0.9910                | _2017-01-04_ | 1/4/17 21:28   |
| I1002-14   | 0.1290               | 0.1310                | 0.1960                | 0.6370                | 0.1310                | 0.1530                 | 1.0600                | _2017-01-04_ | 1/4/17 21:57   |
| I1002-14   | 0.1360               | 0.1300                | 0.1980                | 0.6470                | 0.1290                | 0.1510                 | 1.0360                | _2017-01-04_ | 1/4/17 22:26   |
| CCV        | 5.2290               | 7.5010                | 7.7330                | 24.9550               | 6.3950                | 13.0800                | 36.6530               | _2017-01-04_ | 1/4/17 22:55   |
| CCB        | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | _2017-01-04_ | 1/4/17 23:24   |


M300-0

18d

1/4/17

CALIBRATION OF COMPONENT F-

Method: ANIONS 12-19-16.mtw  
 Equation:  $Q = 0.787523 \cdot A + 4.54623$   
 RSD: 7.224 %  
 Correlation coefficient: 0.998228




K3 = 0      K2 = 0      K1 = 0.787523      K0 = 4.54623

Base: Area  
 Ref.channel: Cond  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | No   |
| 2     | 0.0649 | 2.042 | 0.1   | 20      | 5.749     |      | Yes  |
| 3     | 0.3055 | 6.696 | 0.4   | 20      | 5.749     |      | Yes  |
| 4     | 0.6642 | 15.74 | 0.8   | 20      | 5.749     |      | Yes  |
| 5     | 1.926  | 43.49 | 2     | 20      | 5.749     |      | Yes  |
| 6     | 4.217  | 92.89 | 4     | 20      | 5.749     |      | Yes  |
| 7     | 5.206  | 111.4 | 5     | 20      | 5.749     |      | Yes  |
| 8     | 11.54  | 253.6 | 10    | 20      | 5.749     |      | Yes  |

CALIBRATION OF COMPONENT CL-

Method: ANIONS 12-19-16.mtw  
 Equation:  $Q = 1.17567 \cdot A + 4.48845$   
 RSD: 5.701 %  
 Correlation coefficient: 0.998897




K3 = 0      K2 = 0      K1 = 1.17567      K0 = 4.48845

Base: Area  
 Ref.channel: Cond  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | No   |
| 2     | 0.1287 | 2.693 | 0.15  | 20      | 8.479     |      | Yes  |
| 3     | 0.4753 | 9.278 | 0.6   | 20      | 8.479     |      | Yes  |
| 4     | 0.9236 | 17.71 | 1.2   | 20      | 8.479     |      | Yes  |
| 5     | 2.405  | 45.46 | 3     | 20      | 8.479     |      | Yes  |
| 6     | 5.013  | 93.15 | 6     | 20      | 8.479     |      | Yes  |
| 7     | 6.356  | 118.1 | 7.5   | 20      | 8.479     |      | Yes  |
| 8     | 13.97  | 255.9 | 15    | 20      | 8.479     |      | Yes  |

CALIBRATION OF COMPONENT NO2

Method: ANIONS 12-19-16.mtw  
 Equation:  $Q = 0.507925 \cdot A + 5.31582$   
 RSD: 5.119 %  
 Correlation coefficient: 0.999111




K3 = 0      K2 = 0      K1 = 0.507925      K0 = 5.31582  
 Base: Area  
 Ref.channel: Cond  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc.  | Vol/Dil | Retention | Used | File |
|-------|--------|-------|--------|---------|-----------|------|------|
| 1     | 0      | 0     | 0      | 0       | 0         |      | No   |
| 2     | 0.2222 | 4.727 | 0.152  | 20      | 10.05     |      | Yes  |
| 3     | 0.9375 | 19.53 | 0.608  | 20      | 10.05     |      | Yes  |
| 4     | 1.888  | 39.69 | 1.216  | 20      | 10.05     |      | Yes  |
| 5     | 5.008  | 103.8 | 3.04   | 20      | 10.05     |      | Yes  |
| 6     | 10.59  | 218.1 | 6.08   | 20      | 10.05     |      | Yes  |
| 7     | 13.46  | 278.1 | 7.6    | 20      | 10.05     |      | Yes  |
| 8     | 28.27  | 598.1 | 15.215 | 20      | 10.05     |      | Yes  |

CALIBRATION OF COMPONENT BR-

Method: ANIONS 12-19-16.mtw  
 Equation:  $Q = 2.82735 \cdot A - 0.058448$   
 RSD: 3.611 %  
 Correlation coefficient: 0.999557




K3 = 0      K2 = 0      K1 = 2.82735      K0 = -0.058448

Base: Area  
 Ref.channel: Cond  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |     |
|-------|--------|-------|-------|---------|-----------|------|------|-----|
| 1     | 0      | 0     | 0     | 0       | 0         |      |      | No  |
| 2     | 0.1804 | 6.964 | 0.5   | 20      | 12.84     |      |      | Yes |
| 3     | 0.4772 | 16.47 | 2     | 20      | 12.84     |      |      | Yes |
| 4     | 0.8432 | 29.6  | 4     | 20      | 12.84     |      |      | Yes |
| 5     | 1.853  | 68.87 | 10    | 20      | 12.84     |      |      | Yes |
| 6     | 3.323  | 137.5 | 20    | 20      | 12.84     |      |      | Yes |
| 7     | 3.853  | 171.6 | 25    | 20      | 12.84     |      |      | Yes |
| 8     | 6.801  | 357.8 | 50    | 20      | 12.84     |      |      | Yes |

CALIBRATION OF COMPONENT NO3

Method: ANIONS 12-19-16.mtw  
 Equation:  $Q = 0.459551 \cdot A + 4.85182$   
 RSD: 6.967 %  
 Correlation coefficient: 0.998351




K3 = 0      K2 = 0      K1 = 0.459551      K0 = 4.85182  
 Base: Area  
 Ref.channel: Cond  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | No   |
| 2     | 0.1506 | 3.847 | 0.113 | 20      | 14.38     |      | Yes  |
| 3     | 0.6345 | 16.24 | 0.452 | 20      | 14.38     |      | Yes  |
| 4     | 1.236  | 31.7  | 0.904 | 20      | 14.38     |      | Yes  |
| 5     | 3.27   | 82.68 | 2.26  | 20      | 14.38     |      | Yes  |
| 6     | 6.934  | 174.5 | 4.52  | 20      | 14.38     |      | Yes  |
| 7     | 8.792  | 222.1 | 5.65  | 20      | 14.38     |      | Yes  |
| 8     | 19.15  | 491.7 | 11.3  | 20      | 14.38     |      | Yes  |

CALIBRATION OF COMPONENT HPO4

Method: ANIONS 12-19-16.mtw  
 Equation:  $Q = 1.29463 \cdot A + 0.433908$   
 RSD: 5.455 %  
 Correlation coefficient: 0.998990




K3 = 0      K2 = 0      K1 = 1.29463      K0 = 0.433908

Base: Area  
 Ref.channel: Cond  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height  | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|---------|-------|-------|---------|-----------|------|------|
| 1     | 0       | 0     | 0     | 0       | 0         |      | No   |
| 2     | 0.06533 | 4.069 | 0.242 | 20      | 20.62     | Yes  |      |
| 3     | 0.3076  | 17.53 | 0.968 | 20      | 20.62     | Yes  |      |
| 4     | 0.5813  | 31.44 | 1.936 | 20      | 20.62     | Yes  |      |
| 5     | 1.439   | 74.14 | 4.841 | 20      | 20.62     | Yes  |      |
| 6     | 3.183   | 152.8 | 9.682 | 20      | 20.62     | Yes  |      |
| 7     | 3.734   | 173.6 | 12.1  | 20      | 20.62     | Yes  |      |
| 8     | 8.676   | 378   | 24.21 | 20      | 20.62     | Yes  |      |

CALIBRATION OF COMPONENT SO4

Method: ANIONS 12-19-16.mtw  
 Equation:  $Q = 1.46163 \cdot A + 36.8806$   
 RSD: 7.867 %  
 Correlation coefficient: 0.997898



K3 = 0      K2 = 0      K1 = 1.46163      K0 = 36.8806

Base: Area  
 Ref.channel: Cond  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |     |
|-------|--------|-------|-------|---------|-----------|------|------|-----|
| 1     | 0      | 0     | 0     | 0       | 0         |      |      | No  |
| 2     | 0.2002 | 8.596 | 0.75  | 20      | 23.34     |      |      | Yes |
| 3     | 0.8026 | 32.25 | 3     | 20      | 23.34     |      |      | Yes |
| 4     | 1.621  | 63.56 | 6     | 20      | 23.34     |      |      | Yes |
| 5     | 4.37   | 167.4 | 15    | 20      | 23.34     |      |      | Yes |
| 6     | 9.597  | 357   | 30    | 20      | 23.34     |      |      | Yes |
| 7     | 12.36  | 457.9 | 37.5  | 20      | 23.34     |      |      | Yes |
| 8     | 28.23  | 1025  | 75    | 20      | 23.34     |      |      | Yes |

Report date: 12/20/2016 12:30:23 PM  
Printed by: wet

Ident: STD1  
Analysis from: 12/19/2016 9:23:13 AM  
File: \_2016-12-19\_


Last save: 12/19/2016 2:11:28 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56225

Last save: 12/19/2016 8:39:58 AM

SAMPLE:

: AP/AK  
Vial number: 2  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 12/20/2016 12:30:37 PM  
Printed by: wet


Ident: STD2  
Analysis from: 12/19/2016 9:52:11 AM  
File: \_2016-12-19\_

Last save: 12/19/2016 2:11:28 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56226

Last save: 12/19/2016 8:39:58 AM

SAMPLE:  
: AP/AK  
Vial number: 3  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.75             | 0.380          | 0.06            | 6.36        | 2.042             | 6.20      | 0. |
| 2  | 8.43             | 0.293          | 0.13            | 12.65       | 2.693             | 8.18      | 0. |
| 3  | 9.98             | 0.317          | 0.22            | 21.90       | 4.727             | 14.35     | 0. |
| 4  | 12.65            | 0.570          | 0.18            | 17.74       | 6.964             | 21.14     | 0. |
| 5  | 14.23            | 0.398          | 0.15            | 14.77       | 3.847             | 11.68     | 0. |
| 6  | 20.53            | 0.942          | 0.06            | 6.36        | 4.069             | 12.35     | 0. |
| 7  | 23.08            | 0.614          | 0.20            | 19.70       | 8.596             | 26.10     | 0. |
| 7  | 26.00            | 0.502          | 1.01            | 99.48       | 32.938            | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 12/20/2016 12:30:49 PM  
Printed by: wet


Ident: STD3  
Analysis from: 12/19/2016 10:21:07 AM  
File: \_2016-12-19\_

Last save: 12/19/2016 2:11:28 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56227

Last save: 12/19/2016 8:39:58 AM

SAMPLE:  
: AP/AK  
Vial number: 4  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.75             | 0.309          | 0.31            | 7.74        | 6.696             | 5.67      | 0. |
| 2  | 8.44             | 0.287          | 0.47            | 12.05       | 9.278             | 7.86      | 0. |
| 3  | 9.99             | 0.311          | 0.94            | 23.77       | 19.530            | 16.55     | 0. |
| 4  | 12.62            | 0.525          | 0.48            | 12.10       | 16.474            | 13.96     | 0. |
| 5  | 14.18            | 0.392          | 0.63            | 16.08       | 16.236            | 13.76     | 0. |
| 6  | 20.49            | 0.836          | 0.31            | 7.78        | 17.533            | 14.86     | 0. |
| 7  | 23.11            | 0.596          | 0.80            | 20.35       | 32.252            | 27.33     | 0. |
| 7  | 26.00            | 0.465          | 3.94            | 99.87       | 117.999           | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 12/20/2016 12:31:01 PM  
Printed by: wet


Ident: STD4  
Analysis from: 12/19/2016 10:50:06 AM  
File: \_2016-12-19\_

Last save: 12/19/2016 2:11:28 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56228

Last save: 12/19/2016 8:39:58 AM

SAMPLE:  
: AP/AK  
Vial number: 5  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.74             | 0.305          | 0.66            | 8.56        | 15.739            | 6.86      | 0. |
| 2  | 8.41             | 0.283          | 0.92            | 11.90       | 17.714            | 7.72      | 0. |
| 3  | 9.96             | 0.312          | 1.89            | 24.34       | 39.691            | 17.30     | 0. |
| 4  | 12.58            | 0.546          | 0.84            | 10.86       | 29.598            | 12.90     | 0. |
| 5  | 14.14            | 0.393          | 1.24            | 15.93       | 31.705            | 13.82     | 0. |
| 6  | 20.48            | 0.787          | 0.58            | 7.48        | 31.439            | 13.70     | 0. |
| 7  | 23.12            | 0.586          | 1.62            | 20.89       | 63.556            | 27.70     | 0. |
| 7  | 26.00            | 0.459          | 7.75            | 99.96       | 229.443           | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 12/20/2016 12:31:12 PM  
Printed by: wet


Ident: STD5  
Analysis from: 12/19/2016 11:19:04 AM  
File: \_2016-12-19\_

Last save: 12/19/2016 2:11:28 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56229

Last save: 12/19/2016 8:39:58 AM

SAMPLE:  
: AP/AK  
Vial number: 6  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.77             | 0.300          | 1.92            | 9.49        | 43.490            | 7.42      | 0. |
| 2  | 8.43             | 0.281          | 2.40            | 11.86       | 45.465            | 7.76      | 0. |
| 3  | 9.97             | 0.308          | 5.01            | 24.70       | 103.751           | 17.71     | 0. |
| 4  | 12.61            | 0.605          | 1.85            | 9.14        | 68.874            | 11.76     | 0. |
| 5  | 14.12            | 0.385          | 3.27            | 16.13       | 82.678            | 14.11     | 0. |
| 6  | 20.49            | 0.738          | 1.44            | 7.10        | 74.145            | 12.66     | 0. |
| 7  | 23.16            | 0.574          | 4.37            | 21.55       | 167.362           | 28.57     | 0. |
| 7  | 26.00            | 0.456          | 20.27           | 99.98       | 585.764           | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 12/20/2016 12:31:23 PM  
Printed by: wet


Ident: STD6  
Analysis from: 12/19/2016 11:48:02 AM  
File: \_2016-12-19\_

Last save: 12/19/2016 2:11:30 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56230

Last save: 12/19/2016 8:39:58 AM

SAMPLE:  
: AP/AK  
Vial number: 7  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.75             | 0.296          | 4.22            | 9.84        | 92.887            | 7.58      | 0. |
| 2  | 8.41             | 0.277          | 5.01            | 11.69       | 93.149            | 7.60      | 0. |
| 3  | 9.94             | 0.307          | 10.59           | 24.71       | 218.121           | 17.79     | 0. |
| 4  | 12.62            | 0.721          | 3.32            | 7.75        | 137.500           | 11.22     | 0. |
| 5  | 14.08            | 0.383          | 6.93            | 16.18       | 174.469           | 14.23     | 0. |
| 6  | 20.39            | 0.678          | 3.18            | 7.42        | 152.763           | 12.46     | 0. |
| 7  | 23.11            | 0.560          | 9.60            | 22.39       | 357.038           | 29.12     | 0. |
| 7  | 26.00            | 0.460          | 42.86           | 100.00      | 1225.926          | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 12/20/2016 12:31:35 PM  
Printed by: wet


Ident: STD7  
Analysis from: 12/19/2016 12:16:58 PM  
File: \_2016-12-19\_

Last save: 12/19/2016 2:11:30 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56231

Last save: 12/19/2016 8:39:58 AM

SAMPLE:  
: AP/AK  
Vial number: 8  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.75             | 0.292          | 5.21            | 9.68        | 111.431           | 7.27      | 0. |
| 2  | 8.42             | 0.275          | 6.36            | 11.82       | 118.096           | 7.70      | 0. |
| 3  | 9.96             | 0.307          | 13.46           | 25.04       | 278.117           | 18.14     | 0. |
| 4  | 12.64            | 0.801          | 3.85            | 7.16        | 171.612           | 11.19     | 0. |
| 5  | 14.09            | 0.385          | 8.79            | 16.35       | 222.144           | 14.49     | 0. |
| 6  | 20.41            | 0.659          | 3.73            | 6.94        | 173.635           | 11.33     | 0. |
| 7  | 23.14            | 0.557          | 12.36           | 22.99       | 457.940           | 29.87     | 0. |
| 7  | 26.00            | 0.468          | 53.76           | 100.00      | 1532.975          | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 12/20/2016 12:31:46 PM  
Printed by: wet


Ident: STD8  
Analysis from: 12/19/2016 12:45:57 PM  
File: \_2016-12-19\_

Last save: 12/19/2016 2:12:28 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56232

Last save: 12/19/2016 8:39:58 AM

SAMPLE:  
: AP/AK  
Vial number: 9  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.76             | 0.307          | 11.54           | 9.90        | 253.639           | 7.55      | 0. |
| 2  | 8.42             | 0.272          | 13.98           | 11.98       | 255.919           | 7.62      | 0. |
| 3  | 9.94             | 0.318          | 28.27           | 24.23       | 598.078           | 17.80     | 0. |
| 4  | 12.72            | 1.032          | 6.80            | 5.83        | 357.847           | 10.65     | 0. |
| 5  | 14.09            | 0.394          | 19.15           | 16.42       | 491.697           | 14.63     | 0. |
| 6  | 20.34            | 0.616          | 8.68            | 7.44        | 377.959           | 11.25     | 0. |
| 7  | 23.11            | 0.551          | 28.23           | 24.20       | 1025.166          | 30.51     | 0. |
| 7  | 26.00            | 0.499          | 116.65          | 100.00      | 3360.304          | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 12/20/2016 12:37:54 PM  
Printed by: wet


Ident: ICV  
Analysis from: 12/19/2016 2:14:42 PM  
File: \_2016-12-19\_

Last save: 12/20/2016 12:22:16 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56235

Last save: 12/19/2016 8:39:58 AM

SAMPLE:  
: AP/AK  
Vial number: 12  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.77             | 0.289          | 5.36            | 9.58        | 112.303           | 7.08      | 0. |
| 2  | 8.44             | 0.274          | 6.50            | 11.60       | 120.055           | 7.57      | 0. |
| 3  | 9.96             | 0.306          | 13.58           | 24.25       | 279.010           | 17.59     | 0. |
| 4  | 12.65            | 0.822          | 3.82            | 6.83        | 173.500           | 10.94     | 0. |
| 5  | 14.09            | 0.383          | 10.06           | 17.97       | 252.788           | 15.94     | 0. |
| 6  | 20.44            | 0.654          | 3.83            | 6.85        | 177.344           | 11.18     | 0. |
| 7  | 23.17            | 0.557          | 12.82           | 22.91       | 470.978           | 29.70     | 0. |
| 7  | 26.00            | 0.469          | 55.98           | 99.99       | 1585.979          | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 12/20/2016 12:38:07 PM  
Printed by: wet


Ident: ICB  
Analysis from: 12/19/2016 2:43:39 PM  
File: \_2016-12-19\_

Last save: 12/20/2016 12:36:08 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56236

Last save: 12/19/2016 8:39:58 AM

SAMPLE:  
: AP/AK  
Vial number: 13  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:21:16 PM  
Printed by: wet


Ident: CCV  
Analysis from: 1/4/2017 5:08:01 PM  
File: \_2017-01-04\_

Last save: 1/4/2017 5:33:52 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56438

Last save: 1/4/2017 10:00:52 AM

SAMPLE:  
: AP/AK  
Vial number: 2  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.77             | 0.278          | 6.33            | 10.89       | 125.144           | 7.84      | 0. |
| 2  | 8.42             | 0.261          | 6.85            | 11.79       | 119.730           | 7.50      | 0. |
| 3  | 9.92             | 0.291          | 14.55           | 25.03       | 283.687           | 17.77     | 0. |
| 4  | 12.66            | 0.822          | 3.90            | 6.70        | 175.379           | 10.98     | 0. |
| 5  | 14.03            | 0.373          | 9.33            | 16.05       | 227.087           | 14.22     | 0. |
| 6  | 20.48            | 0.633          | 4.27            | 7.34        | 193.372           | 12.11     | 0. |
| 7  | 23.37            | 0.547          | 12.90           | 22.19       | 472.256           | 29.58     | 0. |
| 7  | 26.00            | 0.458          | 58.14           | 99.99       | 1596.656          | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:21:28 PM  
Printed by: wet


Ident: CCB  
Analysis from: 1/4/2017 5:36:58 PM  
File: \_2017-01-04\_

Last save: 1/5/2017 4:16:22 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56439

Last save: 1/4/2017 10:00:52 AM

SAMPLE:  
: AP/AK  
Vial number: 3  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:22:22 PM  
Printed by: wet


Ident: LB85206BLW  
Analysis from: 1/4/2017 6:05:54 PM  
File: \_2017-01-04\_

Last save: 1/5/2017 4:22:14 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56440

Last save: 1/4/2017 10:00:52 AM

SAMPLE:  
: AP/AK  
Vial number: 4  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:22:32 PM  
 Printed by: wet  
 Ident: LB85206BSW  
 Analysis from: 1/4/2017 6:34:51 PM  
 File: 2017-01-04


241  
 1/21/18

Last save: 1/5/2017 4:22:14 PM

Method: ANIONS 12-19-16.mtw  
 Run operator: wet  
 Analysis number: 56441

Last save: 1/4/2017 10:00:52 AM

SAMPLE:  
 : AP/AK  
 Vial number: 5  
 Volume: 20.0 µL  
 Dilution: 1.00  
 Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.85             | 0.280          | 6.16            | 10.48       | 121.162           | 7.40      | 0. |
| 2  | 8.57             | 0.267          | 6.87            | 11.69       | 122.748           | 7.50      | 0. |
| 3  | 10.12            | 0.299          | 14.44           | 24.57       | 288.832           | 17.64     | 0. |
| 4  | 12.89            | 0.848          | 3.77            | 6.41        | 175.300           | 10.71     | 0. |
| 5  | 14.27            | 0.372          | 10.45           | 17.78       | 253.816           | 15.50     | 0. |
| 6  | 20.78            | 0.644          | 4.26            | 7.25        | 196.167           | 11.98     | 0. |
| 7  | 23.73            | 0.558          | 12.82           | 21.81       | 479.011           | 29.26     | 0. |
| 7  | 26.00            | 0.467          | 58.78           | 100.00      | 1637.036          | 100.00    | 0. |

This report has been created by IC Net  
 METROHM LTD

Report date: 1/5/2017 4:22:48 PM  
Printed by: wet


Ident: I1002-08  
Analysis from: 1/4/2017 7:03:48 PM  
File: \_2017-01-04\_

Last save: 1/5/2017 10:24:48 AM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56442

Last save: 12/20/2016 12:26:32 PM

SAMPLE:  
: AP/AK  
Vial number: 10  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.77             | 0.282          | 0.09            | 8.29        | 1.860             | 5.34      | 0. |
| 2  | 8.48             | 0.274          | 0.14            | 12.66       | 2.610             | 7.49      | 0. |
| 3  | 10.04            | 0.303          | 0.23            | 20.25       | 4.681             | 13.43     | 0. |
| 4  | 12.71            | 0.577          | 0.20            | 17.47       | 7.773             | 22.30     | 0. |
| 5  | 14.23            | 0.371          | 0.17            | 14.47       | 3.942             | 11.31     | 0. |
| 6  | 20.89            | 0.860          | 0.08            | 7.35        | 5.086             | 14.59     | 0. |
| 7  | 23.74            | 0.607          | 0.22            | 19.06       | 8.900             | 25.54     | 0. |
| 7  | 26.00            | 0.468          | 1.14            | 99.55       | 34.852            | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:23:00 PM  
Printed by: wet


Ident: I1002-09  
Analysis from: 1/4/2017 7:32:44 PM  
File: \_2017-01-04\_

Last save: 1/5/2017 10:24:48 AM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56443

Last save: 12/20/2016 12:26:32 PM

SAMPLE:  
: AP/AK  
Vial number: 11  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.77             | 0.282          | 0.09            | 8.44        | 1.843             | 5.36      | 0. |
| 2  | 8.47             | 0.273          | 0.12            | 10.67       | 2.154             | 6.26      | 0. |
| 3  | 10.04            | 0.303          | 0.25            | 22.14       | 5.063             | 14.72     | 0. |
| 4  | 12.68            | 0.581          | 0.20            | 17.75       | 7.864             | 22.87     | 0. |
| 5  | 14.24            | 0.377          | 0.16            | 13.84       | 3.795             | 11.04     | 0. |
| 6  | 20.93            | 0.819          | 0.08            | 7.29        | 4.798             | 13.95     | 0. |
| 7  | 23.71            | 0.590          | 0.22            | 19.51       | 8.872             | 25.80     | 0. |
| 7  | 26.00            | 0.461          | 1.12            | 99.63       | 34.389            | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:23:10 PM  
Printed by: wet


Ident: I1002-10  
Analysis from: 1/4/2017 8:01:41 PM  
File: \_2017-01-04\_

Last save: 1/5/2017 10:24:48 AM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56444

Last save: 12/20/2016 12:26:32 PM

SAMPLE:  
: AP/AK  
Vial number: 12  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.76             | 0.290          | 0.09            | 7.97        | 1.736             | 5.12      | 0. |
| 2  | 8.43             | 0.279          | 0.12            | 10.93       | 2.236             | 6.60      | 0. |
| 3  | 9.94             | 0.294          | 0.25            | 22.05       | 4.796             | 14.15     | 0. |
| 4  | 12.57            | 0.586          | 0.20            | 17.56       | 7.739             | 22.83     | 0. |
| 5  | 14.12            | 0.382          | 0.16            | 14.00       | 3.857             | 11.38     | 0. |
| 6  | 20.83            | 0.839          | 0.08            | 6.75        | 4.447             | 13.12     | 0. |
| 7  | 23.60            | 0.588          | 0.23            | 20.36       | 9.081             | 26.79     | 0. |
| 7  | 26.00            | 0.465          | 1.11            | 99.63       | 33.891            | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:23:21 PM  
Printed by: wet


Ident: I1002-11  
Analysis from: 1/4/2017 8:30:38 PM  
File: \_2017-01-04\_

Last save: 1/5/2017 10:24:48 AM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56445

Last save: 12/20/2016 12:26:32 PM

SAMPLE:  
: AP/AK  
Vial number: 13  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.76             | 0.287          | 0.09            | 8.35        | 1.899             | 5.49      | 0. |
| 2  | 8.46             | 0.274          | 0.12            | 10.81       | 2.178             | 6.30      | 0. |
| 3  | 10.02            | 0.304          | 0.25            | 22.19       | 5.180             | 14.97     | 0. |
| 4  | 12.77            | 0.584          | 0.20            | 17.66       | 7.800             | 22.54     | 0. |
| 5  | 14.30            | 0.371          | 0.15            | 13.82       | 3.954             | 11.43     | 0. |
| 6  | 20.84            | 0.825          | 0.08            | 7.29        | 4.747             | 13.72     | 0. |
| 7  | 23.66            | 0.593          | 0.22            | 19.46       | 8.841             | 25.55     | 0. |
| 7  | 26.00            | 0.463          | 1.11            | 99.58       | 34.599            | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:23:32 PM  
Printed by: wet


Ident: I1002-12  
Analysis from: 1/4/2017 8:59:35 PM  
File: \_2017-01-04\_

Last save: 1/5/2017 10:24:50 AM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56446

Last save: 12/20/2016 12:26:32 PM

SAMPLE:  
: AP/AK  
Vial number: 10  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.88             | 0.307          | 0.10            | 8.40        | 2.205             | 6.42      | 0. |
| 2  | 8.58             | 0.266          | 0.15            | 12.65       | 2.514             | 7.32      | 0. |
| 3  | 10.11            | 0.295          | 0.23            | 20.25       | 4.573             | 13.32     | 0. |
| 4  | 12.73            | 0.570          | 0.20            | 17.33       | 7.732             | 22.52     | 0. |
| 5  | 14.25            | 0.374          | 0.17            | 14.47       | 4.053             | 11.81     | 0. |
| 6  | 20.94            | 0.813          | 0.08            | 7.14        | 4.416             | 12.86     | 0. |
| 7  | 23.81            | 0.590          | 0.22            | 19.27       | 8.836             | 25.74     | 0. |
| 7  | 26.00            | 0.459          | 1.15            | 99.52       | 34.330            | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:23:42 PM  
Printed by: wet


Ident: I1002-13  
Analysis from: 1/4/2017 9:28:32 PM  
File: \_2017-01-04\_

Last save: 1/5/2017 10:24:50 AM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56447

Last save: 12/20/2016 12:26:32 PM

SAMPLE:  
: AP/AK  
Vial number: 11  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.77             | 0.320          | 0.08            | 7.59        | 1.965             | 5.93      | 0. |
| 2  | 8.42             | 0.268          | 0.12            | 10.86       | 2.130             | 6.43      | 0. |
| 3  | 9.95             | 0.296          | 0.24            | 22.21       | 4.717             | 14.24     | 0. |
| 4  | 12.56            | 0.572          | 0.19            | 17.94       | 7.535             | 22.76     | 0. |
| 5  | 14.09            | 0.370          | 0.15            | 14.08       | 3.650             | 11.02     | 0. |
| 6  | 20.64            | 0.848          | 0.08            | 7.12        | 4.534             | 13.69     | 0. |
| 7  | 23.46            | 0.592          | 0.21            | 19.79       | 8.582             | 25.92     | 0. |
| 7  | 26.00            | 0.467          | 1.08            | 99.58       | 33.114            | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:23:52 PM  
Printed by: wet


Ident: I1002-14  
Analysis from: 1/4/2017 9:57:28 PM  
File: \_2017-01-04\_

Last save: 1/5/2017 10:24:50 AM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56448

Last save: 12/20/2016 12:26:32 PM

SAMPLE:  
: AP/AK  
Vial number: 12  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.79             | 0.288          | 0.09            | 8.01        | 1.796             | 5.13      | 0. |
| 2  | 8.43             | 0.267          | 0.13            | 10.95       | 2.214             | 6.33      | 0. |
| 3  | 9.95             | 0.294          | 0.25            | 21.80       | 4.935             | 14.10     | 0. |
| 4  | 12.56            | 0.571          | 0.20            | 17.37       | 7.730             | 22.09     | 0. |
| 5  | 14.09            | 0.372          | 0.16            | 13.89       | 3.855             | 11.01     | 0. |
| 6  | 20.65            | 0.837          | 0.09            | 7.59        | 5.110             | 14.60     | 0. |
| 7  | 23.53            | 0.600          | 0.23            | 20.03       | 9.360             | 26.74     | 0. |
| 7  | 26.00            | 0.461          | 1.14            | 99.63       | 35.000            | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:24:02 PM  
Printed by: wet


Ident: I1002-14  
Analysis from: 1/4/2017 10:26:25 PM  
File: \_2017-01-04\_

Last save: 1/5/2017 10:24:50 AM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56449

Last save: 12/20/2016 12:26:32 PM

SAMPLE:  
: AP/AK  
Vial number: 13  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.76             | 0.288          | 0.10            | 8.53        | 1.933             | 5.54      | 0. |
| 2  | 8.44             | 0.277          | 0.12            | 10.64       | 2.189             | 6.28      | 0. |
| 3  | 10.00            | 0.301          | 0.25            | 21.93       | 5.001             | 14.34     | 0. |
| 4  | 12.68            | 0.583          | 0.20            | 17.43       | 7.799             | 22.37     | 0. |
| 5  | 14.23            | 0.373          | 0.16            | 13.84       | 3.777             | 10.83     | 0. |
| 6  | 20.82            | 0.847          | 0.09            | 7.70        | 5.087             | 14.59     | 0. |
| 7  | 23.69            | 0.603          | 0.22            | 19.46       | 9.083             | 26.05     | 0. |
| 7  | 26.00            | 0.467          | 1.13            | 99.52       | 34.868            | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:24:13 PM  
Printed by: wet


Ident: CCV  
Analysis from: 1/4/2017 10:55:22 PM  
File: \_2017-01-04\_

Last save: 1/4/2017 11:21:12 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56450

Last save: 1/4/2017 10:00:52 AM

SAMPLE:  
: AP/AK  
Vial number: 2  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>uS/cm | Height<br>% | Area<br>uS/cm*sec | Area<br>% |    |
|----|------------------|----------------|-----------------|-------------|-------------------|-----------|----|
| 1  | 5.83             | 0.286          | 6.23            | 10.60       | 127.019           | 7.62      | 0. |
| 2  | 8.53             | 0.268          | 6.92            | 11.76       | 123.795           | 7.43      | 0. |
| 3  | 10.07            | 0.299          | 14.69           | 24.98       | 294.015           | 17.64     | 0. |
| 4  | 12.85            | 0.831          | 3.91            | 6.65        | 176.549           | 10.59     | 0. |
| 5  | 14.33            | 0.437          | 9.45            | 16.07       | 267.765           | 16.06     | 0. |
| 6  | 20.93            | 0.632          | 4.46            | 7.58        | 201.735           | 12.10     | 0. |
| 7  | 23.83            | 0.543          | 13.14           | 22.35       | 476.300           | 28.57     | 0. |
| 7  | 26.00            | 0.471          | 58.79           | 100.00      | 1667.178          | 100.00    | 0. |

This report has been created by IC Net  
METROHM LTD

Report date: 1/5/2017 4:24:26 PM  
Printed by: wet


Ident: CCB  
Analysis from: 1/4/2017 11:24:19 PM  
File: \_2017-01-04\_

Last save: 1/5/2017 4:16:22 PM

Method: ANIONS 12-19-16.mtw  
Run operator: wet  
Analysis number: 56451

Last save: 1/4/2017 10:00:52 AM

SAMPLE:  
: AP/AK  
Vial number: 3  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD