

Run Number: LB86564

BalanceID: WC SC-3

DATE & TIME IN: 03/27/2017 11:50

DATE & TIME OUT: 03/28/2017 16:15

OVEN TEMP: 104 °C

OVEN TEMP: 105 °C

Dish #	Lab ID	Client ID	Empty Dish Weight (g)	Final Empty Dish Weight (g)	Sample Volume (ml)	1st Empty Dish+Sample weight after 1.5hr drying @105°C (g)	2nd Empty Dish+Sample weight after 1.5hr drying @105°C (g)	Final Empty Dish+Sample weight after 1.5hr drying @105°C (g)	Weight (g)	Result mg/L
1	LB86564BL	LB86564BL	3.1403	3.1403	100	3.1403	3.1403	3.1403	0	0
2	LB86564BS	LB86564BS	1.4305	1.4305	100	1.4851	1.4851	1.4851	0.0546	546
3	I2387-01	032317-PRRT-E-D-M	1.4157	1.4157	400	1.4347	1.4347	1.4347	0.019	47.5
4	I2387-02	032317-PRRT-E-U-M	1.4112	1.4112	400	1.4657	1.4657	1.4657	0.0545	136.3
5	I2387-03	032317-TURB-E-D-M	1.4245	1.4245	400	1.4722	1.4722	1.4722	0.0477	119.3
6	I2388-01	032217-PRRT-E-D-M	1.4145	1.4145	400	1.4342	1.4342	1.4342	0.0197	49.3
7	I2388-02	032217-PRRT-E-U-M	1.4123	1.4123	400	1.5977	1.5977	1.5977	0.1854	463.5
8	I2396-01	032317-TURB-E-U-M	1.4152	1.4152	400	1.4693	1.4693	1.4693	0.0541	135.3
9	I2421-01	HHD-GW1	1.4108	1.4108	400	1.4536	1.4536	1.4536	0.0428	107
10	I2421-01DUP	HHD-GW1DUP	1.4203	1.4203	400	1.465	1.465	1.465	0.0447	111.7
11	I2421-03	HHD-GW2	1.4176	1.4176	200	1.486	1.486	1.486	0.0684	342

A = Sample Volume (ml)

B = Empty Dish Weight (g)

C = Final Empty Dish + Sample weight after 1.5 hr drying @105°C(g)

D = Weight (g)

$$\text{Weight (g)} = C - B$$

$$\text{Result mg/L} = \frac{D}{A} * 1000 * 1000$$

WORKLIST(Hardcopy Internal Chain)

1386564

WorkList Name : TSS I2387

WorkList ID : 96856

Date : 3/27/2017 10:28:34 AM

Due Date	Matrix	Sample	Test	Preservative	Customer	Storage Location	Customer Sample	Collect Date	Method
	Water	I2387-01 -B	TSS	Cool 4 deg C	HDR103	M31	032317-PRRT-E-D-M	03/23/2017	SM2540 D
	Water	I2387-02 -B	TSS	Cool 4 deg C	HDR103	M31	032317-PRRT-E-U-M	03/23/2017	SM2540 D
	Water	I2387-03	TSS	Cool 4 deg C	HDR103	M31	032317-TURB-E-D-M	03/23/2017	SM2540 D
	Water	I2388-01 -B	TSS	Cool 4 deg C	HDR103	M33	032217-PRRT-E-D-M	03/22/2017	SM2540 D
	Water	I2388-02 -B	TSS	Cool 4 deg C	HDR103	M33	032217-PRRT-E-U-M	03/22/2017	SM2540 D
	Water	I2396-01	TSS	Cool 4 deg C	HDR103	M12	032317-TURB-E-U-M	03/23/2017	SM2540 D
	Water	I2421-01 -G	TSS	Cool 4 deg C	DEWB01	M41	HHD-GW1	03/24/2017	SM2540 D
	Water	I2421-03 -G	TSS	Cool 4 deg C	DEWB01	M41	HHD-GW2	03/24/2017	SM2540 D

Date/Time 03/27/17 10:28 AM
 Received by: JP
 Relinquished by: JP

Date/Time 03/27/17 1:00 PM
 Received by: CD
 Relinquished by: JP