

DATE & TIME IN: 05/09/2017 17:30

OVEN TEMP: 98 °C

Run Number: LB87361

DATE & TIME OUT: 05/10/2017 12:50

OVEN TEMP: 105 °C

BalanceID: WC SC-3

OvenID: WC Oven-2

Dish #	Lab ID	Client ID	Empty Dish Weight (g)	Final Empty Dish Weight (g)	Sample Vol (ml)	1st Dish+Sample weight after Drying @103-@105°C (g)	2nd Dish+Sample weight after Drying @103-@105°C (g)	Final Dish+Sample weight after Drying @103-@105°C (g)	Weight (g)	Result (mg/L)
1	LB87361BL	LB87361BL	68.0149	68.0149	100	68.0149	68.0149	68.0149	0	0
2	I3046-01	TWP-13	75.261	75.261	100	75.2997	75.2997	75.2997	0.0387	387
3	I3046-02	TWP-11	79.5217	79.5217	100	79.5986	79.5986	79.5986	0.0769	769
4	I3046-02DUP	TWP-11DUP	85.7632	85.7632	100	85.8404	85.8404	85.8404	0.0772	772
5	I3062-15	TWP-02	84.6074	84.6074	100	85.6816	85.6816	85.6816	1.0742	10742

A = Sample Volume (ml)
 B = Empty Dish Weight (g)
 C = Final Empty Dish + Sample weight after drying @105°C(g)
 D = Weight (g)

Weight (g) = C - B

Result mg/L = $\frac{D}{A} \times 1000 \times 1000$

WORKLIST(Hardcopy Internal Chain)

WorkList Name : TS I3046

WorkList ID : 98260

Date : 5/9/2017 10:34:00 AM

Due Date	Matrix	Sample	Test	Preservative	Customer	Storage Location	Customer Sample	Collect Date	Method
	Water	I3046-01 -F	TS	Cool 4 deg C	LIRO01	F62	TWP-13	05/08/2017	SM2540 B
	Water	I3046-02 -F	TS	Cool 4 deg C	LIRO01	F62	TWP-11	05/08/2017	SM2540 B
	Water	I3062-15 -F	TS	Cool 4 deg C	LIRO01	G52	TWP-02	05/09/2017	SM2540 B

Date/Time 05/09/17 11:00 AM
 Received by: TS
 Relinquished by: CD

Date/Time 05/09/17 6:35 PM
 Received by: CD
 Relinquished by: TS