



PERCENT SOLID

Analyst Name: sweta
Date: 7/27/2018

OVENTEMP IN Celsius(°C): 108
Time IN: 15:45
In Date: 12/26/2017
Weight Check 1.0g: 1.00 g
Weight Check 10g: 10.00 g
OvenID: M Oven-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:10
Out Date: 12/27/2017
Weight Check 1.0g: 1.00 g
Weight Check 10g: 10.00 g
BalanceID: M SC-1

QC: LB92380

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
I7073-01	BE4R5	1	1.14	9.96	8.08	78.7
I7073-02	BE4R5MS	2	1.14	9.96	8.08	78.7
I7073-03	BE4R5MSD	3	1.14	9.96	8.08	78.7
I7073-04	BE4R6	4	1.14	9.83	9.24	93.2
I7073-05	BE4R7	5	1.15	9.98	7.89	76.3
I7073-06	BE4R8	6	1.13	9.77	8.85	89.4
I7073-07	BE4R9	7	1.16	9.88	9.30	93.3
I7073-08	BE4S0	8	1.15	9.96	9.31	92.6
I7073-09	BE4S2	9	1.14	9.87	8.55	84.9
I7073-10	BE4S3	10	1.14	9.94	8.20	80.2
I7073-11	BE4S4	11	1.13	9.93	8.58	84.7
I7073-12	BE4S5	12	1.15	9.68	8.40	85
I7073-13	BE4S6	13	1.15	9.81	8.44	84.2
I7073-14	BE4S7	14	1.14	9.90	8.34	82.2
I7073-15	BE4S8	15	1.14	9.90	8.47	83.7
I7073-16	BE4S9	16	1.14	9.96	8.36	81.9
I7073-17	BE4T0	17	1.14	9.90	7.69	74.8
I7073-18	BE4T1	18	1.16	9.89	7.83	76.4
I7073-19	BE4T2	19	1.15	9.70	7.61	75.6

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$