

Clear table


ident	concentrati on F-	concentrati on CL-	concentrati on NO2	concentrati on BR-	concentrati on NO3	concentrati on HPO4	concentrati on SO4	file name	date time
STD1	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	2018-04-04	4/4/18 11:13
STD2	0.0910	0.2170	0.2450	0.7070	0.1820	-0.3600	1.3260	2018-04-04	4/4/18 11:44
STD3	0.4090	0.6350	0.6260	2.0800	0.5220	1.1120	3.1540	2018-04-04	4/4/18 12:15
STD4	0.7630	1.1880	1.1820	3.9530	0.9880	2.4190	5.8100	2018-04-04	4/4/18 12:45
STD5	1.8900	2.8590	2.8540	9.5410	2.3830	5.2270	14.0460	2018-04-04	4/4/18 13:16
STD6	3.9980	5.8710	5.8410	19.5730	4.8890	9.7730	29.2050	2018-04-04	4/4/18 13:47
STD7	5.0330	7.6390	7.6770	25.4480	6.3640	12.4030	38.3350	2018-04-04	4/4/18 14:18
ICV	1.9450	2.8980	2.9180	9.7820	2.4460	5.4650	15.1150	2018-04-04	4/4/18 15:51
ICB	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	2018-04-04	4/4/18 16:22
CCV	2.0770	2.9870	2.9640	10.0130	2.5020	5.6200	14.9420	2018-04-05	4/5/18 15:47
CCB	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	2018-04-05	4/5/18 16:18
LB94432BLW	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	2018-04-05	4/5/18 16:48
LB94432BSW	2.0090	2.8640	2.8620	9.6100	2.4080	5.5020	14.9000	2018-04-05	4/5/18 17:19
J2133-01	0.1120	0.2000	0.2090	0.5840	0.1580	-0.4860	1.1940	2018-04-05	4/5/18 17:50
J2133-03	0.1060	0.1860	0.2110	0.5610	0.1510	-0.4910	1.2060	2018-04-05	4/5/18 18:21
J2133-02	0.2500	0.3930	0.4010	1.2880	0.3250	0.3750	2.1460	2018-04-05	4/5/18 18:52
CCV	2.1320	2.9970	2.9880	10.0490	2.5120	5.6790	14.9460	2018-04-05	4/5/18 20:25
CCB	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	2018-04-05	4/5/18 20:56

IC-2 9056
4-6-18 ALS

CALIBRATION OF COMPONENT F-

Method: AnionsIC2-040418.mtw
 Equation: $Q = 0.0173137 \cdot A + 1.31036$
 RSD: 2.839 %
 Correlation coefficient: 0.999671

IC-2 cal
 4/4/18 AK




K3 = 0 K2 = 0 K1 = 0.0173137 K0 = 1.31036
 Base: Area
 Ref.channel: ch1
 ISTD:
 Formula: Linear
 Weight: 1

Level	Height	Area	Conc.	Vol/Dil	Retention	Used	File
1	0	0	0	0	0		
2	4.414	93.57	0.1	20	5.099		
3	21.15	417.4	0.4	20	5.099		
4	42.68	835.1	0.8	20	5.099		
5	108.8	2123	2	20	5.099		
6	236.9	4547	4	20	5.099		
7	299.8	5739	5	20	5.099		

CALIBRATION OF COMPONENT CL-

Method: AnionsIC2-040418.mtw
 Equation: $Q = 0.0274608 \cdot A + 2.15555$
 RSD: 3.949 %
 Correlation coefficient: 0.999363




K3 = 0 K2 = 0 K1 = 0.0274608 K0 = 2.15555
 Base: Area
 Ref.channel: chl
 ISTD:
 Formula: Linear
 Weight: 1

Level	Height	Area	Conc.	Vol/Dil	Retention	Used	File
1	0	0	0	0	0		0
2	5.3	85.83	0.15	20	8.753		
3	21.46	389.6	0.6	20	8.753		
4	44	792	1.2	20	8.753		
5	113.2	2014	3	20	8.753		
6	238.6	4200	6	20	8.753		
7	313.6	5485	7.5	20	8.753		

CALIBRATION OF COMPONENT NO2

Method: AnionsIC2-040418.mtw
 Equation: $Q = 0.0118752 \cdot A + 3.14318$
 RSD: 4.965 %
 Correlation coefficient: 0.998993




K3 = 0 K2 = 0 K1 = 0.0118752 K0 = 3.14318
 Base: Area
 Ref.channel: ch1
 ISTD:
 Formula: Linear
 Weight: 1

Level	Height	Area	Conc.	Vol/Dil	Retention	Used	File
1	0	0	0	0	0		0
2	9.067	172.6	0.15	20	11.21		
3	40.62	804.9	0.6	20	11.21		
4	86.96	1735	1.2	20	11.21		
5	228.7	4556	3	20	11.21		
6	492.2	9551	6	20	11.21		
7	651.2	1.266e+04	7.5	20	11.21		

CALIBRATION OF COMPONENT BR-

Method: AnionsIC2-040418.mtw
 Equation: $Q = 0.062906 \cdot A + 6.52841$
 RSD: 3.798 %
 Correlation coefficient: 0.999411




K3 = 0 K2 = 0 K1 = 0.062906 K0 = 6.52841
 Base: Area
 Ref.channel: ch1
 ISTD:
 Formula: Linear
 Weight: 1

Level	Height	Area	Conc.	Vol/Dil	Retention	Used	File
1	0	0	0	0	0		0
2	4.778	134.5	0.5	20	14.51		
3	19.35	572.1	2	20	14.51		
4	39.52	1163	4	20	14.51		
5	99.75	2947	10	20	14.51		
6	205.6	6126	20	20	14.51		
7	267.5	7987	25	20	14.51		

CALIBRATION OF COMPONENT NO3

Method: AnionsIC2-040418.mtw
 Equation: $Q = 0.0112644 \cdot A + 1.68428$
 RSD: 3.858 %
 Correlation coefficient: 0.999392




K3 = 0 K2 = 0 K1 = 0.0112644 K0 = 1.68428
 Base: Area
 Ref.channel: ch1
 ISTD:
 Formula: Linear
 Weight: 1

Level	Height	Area	Conc.	Vol/Dil	Retention	Used	File
1	0	0	0	0	0		0
2	5.55	188.2	0.125	20	17.18		
3	22.4	792.4	0.5	20	17.18		
4	45.33	1618	1	20	17.18		
5	111.7	4108	2.5	20	17.18		
6	222.7	8546	5	20	17.18		
7	285.2	1.115e+04	6.25	20	17.18		

CALIBRATION OF COMPONENT HPO4

Method: AnionsIC2-040418.mtw
 Equation: $Q = 0.026679 \cdot A - 20.7096$
 RSD: 8.060 %
 Correlation coefficient: 0.997346




K3 = 0 K2 = 0 K1 = 0.026679 K0 = -20.7096
 Base: Area
 Ref.channel: ch1
 ISTD:
 Formula: Linear
 Weight: 1

Level	Height	Area	Conc.	Vol/Dil	Retention	Used	File
1	0	0	0	0	0		0
2	11.05	522.3	0.25	20	23.13		
3	35.53	1627	1	20	23.13		
4	58.6	2606	2	20	23.13		
5	110.8	4740	5	20	23.13		
6	194.8	8140	10	20	23.13		
7	242.9	1.007e+04	12.5	20	23.13		

CALIBRATION OF COMPONENT SO4

Method: AnionsIC2-040418.mtw
 Equation: $Q = 0.0393001 \cdot A + 15.3906$
 RSD: 5.023 %
 Correlation coefficient: 0.998970



K3 = 0 K2 = 0 K1 = 0.0393001 K0 = 15.3906
 Base: Area
 Ref.channel: chl
 ISTD:
 Formula: Linear
 Weight: 1

Level	Height	Area	Conc.	Vol/Dil	Retention	Used	File
1	0	0	0	0	0		0
2	8.569	314.2	0.75	20	25.36		
3	33.33	1247	3	20	25.36		
4	68.27	2592	6	20	25.36		
5	176.5	6814	15	20	25.36		
6	373.3	1.451e+04	30	20	25.36		
7	489.1	1.912e+04	37.5	20	25.36		

Report date: 4/5/2018 1:11:03 PM
Printed by: wet


Ident: STD1
Analysis from: 4/4/2018 11:13:10 AM
File: _2018-04-04_

Last save: 4/4/2018 11:41:11 AM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17014

Last save: 4/4/2018 9:46:4

SAMPLE:
: AK/AP
Vial number: 2
Volume: 20.0 μ L
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net
METROHM LTD

Report date: 4/5/2018 1:14:52 PM
Printed by: wet


Ident: STD2
Analysis from: 4/4/2018 11:44:06 AM
File: _2018-04-04_

Last save: 4/4/2018 12:12:06 PM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17015

Last save: 4/4/2018 11:41:

SAMPLE:
: AK/AP
Vial number: 3
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.02	0.268	4.41	9.06	93.574	6.19
2	8.65	0.226	5.30	10.87	85.828	5.68
3	11.08	0.258	9.07	18.61	172.646	11.42
4	14.55	0.371	4.78	9.80	134.506	8.90
5	17.31	0.438	5.55	11.39	188.238	12.46
6	22.42	0.648	11.05	22.68	522.285	34.56
7	24.85	0.526	8.57	17.58	314.164	20.79
7	28.00	0.391	48.72	99.99	1511.240	100.00

Report date: 4/5/2018 1:15:03 PM
Printed by: wet


Ident: STD3
Analysis from: 4/4/2018 12:15:02 PM
File: _2018-04-04_

Last save: 4/4/2018 12:43:02 PM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17016

Last save: 4/4/2018 12:12:

SAMPLE:
: AK/AP
Vial number: 4
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.03	0.251	21.15	10.91	417.353	7.13
2	8.65	0.231	21.46	11.07	389.636	6.66
3	11.08	0.259	40.61	20.95	804.908	13.76
4	14.52	0.373	19.35	9.98	572.131	9.78
5	17.25	0.446	22.40	11.56	792.363	13.54
6	22.41	0.625	35.53	18.33	1627.429	27.82
7	24.87	0.527	33.33	17.20	1246.673	21.31
7	28.00	0.387	193.83	100.00	5850.494	100.00

Report date: 4/5/2018 1:15:13 PM
Printed by: wet


Ident: STD4
Analysis from: 4/4/2018 12:45:57 PM
File: _2018-04-04_

Last save: 4/4/2018 1:13:58 PM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17017

Last save: 4/4/2018 12:43:

SAMPLE:
: AK/AP
Vial number: 5
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.03	0.245	42.67	11.07	835.098	7.36
2	8.65	0.229	43.99	11.42	792.048	6.98
3	11.09	0.256	86.96	22.57	1735.400	15.30
4	14.49	0.372	39.51	10.25	1163.216	10.26
5	17.21	0.451	45.33	11.76	1617.543	14.26
6	22.39	0.606	58.59	15.21	2605.818	22.98
7	24.87	0.528	68.27	17.72	2592.349	22.86
7	28.00	0.384	385.33	100.00	11341.473	100.00

This report has been created by IC Net
METROHM LTD

Report date: 4/5/2018 1:15:23 PM
Printed by: wet


Ident: STD5
Analysis from: 4/4/2018 1:16:53 PM
File: _2018-04-04_

Last save: 4/4/2018 1:44:53 PM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17018

Last save: 4/4/2018 1:13:5

SAMPLE:
: AK/AP
Vial number: 6
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.03	0.240	108.80	11.46	2122.570	7.77
2	8.64	0.226	113.15	11.92	2013.510	7.38
3	11.09	0.254	228.66	24.09	4556.349	16.69
4	14.44	0.375	99.75	10.51	2946.942	10.79
5	17.13	0.470	111.68	11.76	4108.205	15.05
6	22.33	0.585	110.75	11.67	4739.868	17.36
7	24.84	0.529	176.47	18.59	6813.646	24.96
7	28.00	0.383	949.26	100.00	27301.089	100.00

This report has been created by IC Net
METROHM LTD

Report date: 4/5/2018 1:15:32 PM
Printed by: wet


Ident: STD6
Analysis from: 4/4/2018 1:47:48 PM
File: _2018-04-04_

Last save: 4/4/2018 2:15:48 PM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17019

Last save: 4/4/2018 1:44:5

SAMPLE:
: AK/AP
Vial number: 7
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.04	0.245	236.89	12.06	4546.982	8.17
2	8.64	0.224	238.63	12.15	4200.028	7.55
3	11.11	0.247	492.22	25.06	9550.815	17.17
4	14.40	0.379	205.58	10.47	6126.396	11.01
5	17.04	0.495	222.70	11.34	8545.599	15.36
6	22.28	0.576	194.81	9.92	8140.026	14.63
7	24.81	0.531	373.26	19.00	14512.769	26.09
7	28.00	0.385	1964.09	100.00	55622.615	100.00

This report has been created by IC Net
METROHM LTD

Report date: 4/5/2018 1:15:44 PM
Printed by: wet


Ident: STD7
Analysis from: 4/4/2018 2:18:43 PM
File: _2018-04-04_

Last save: 4/4/2018 2:46:43 PM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17020

Last save: 4/4/2018 2:15:4

SAMPLE:
: AK/AP
Vial number: 8
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.05	0.244	299.83	11.76	5738.732	7.95
2	8.64	0.222	313.57	12.30	5485.294	7.60
3	11.13	0.249	651.17	25.54	12664.999	17.54
4	14.39	0.381	267.53	10.49	7987.104	11.06
5	17.02	0.507	285.20	11.19	11149.350	15.44
6	22.25	0.574	242.91	9.53	10073.915	13.95
7	24.80	0.534	489.11	19.19	19117.198	26.47
7	28.00	0.387	2549.32	100.00	72216.591	100.00

This report has been created by IC Net
METROHM LTD

Report date: 4/5/2018 1:16:28 PM
Printed by: wet


Ident: ICV
Analysis from: 4/4/2018 3:51:30 PM
File: _2018-04-04_

Last save: 4/4/2018 4:19:30 PM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17023

Last save: 4/4/2018 2:46:4

SAMPLE:
: AK/AP
Vial number: 11
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.03	0.243	111.87	11.43	2171.098	7.69
2	8.64	0.226	113.94	11.64	2031.808	7.20
3	11.09	0.250	236.24	24.14	4649.853	16.47
4	14.43	0.374	101.76	10.40	3006.153	10.65
5	17.11	0.468	114.13	11.66	4193.650	14.86
6	22.34	0.591	113.04	11.55	4872.771	17.26
7	24.85	0.534	187.58	19.17	7300.652	25.87
7	28.00	0.384	978.55	100.00	28225.986	100.00

This report has been created by IC Net
METROHM LTD

Report date: 4/5/2018 1:16:49 PM
Printed by: wet


Ident: ICB
Analysis from: 4/4/2018 4:22:25 PM
File: _2018-04-04_

Last save: 4/5/2018 1:16:45 PM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17024

Last save: 4/4/2018 2:46:4

SAMPLE:
: AK/AP
Vial number: 12
Volume: 20.0 μ L
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net
METROHM LTD

Report date: 4/6/2018 11:37:57 AM
Printed by: wet


Ident: CCV
Analysis from: 4/5/2018 3:47:07 PM
File: _2018-04-05_

Last save: 4/5/2018 4:15:08 PM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17038

Last save: 4/4/2018 2:46:4

SAMPLE:
: AK/AP
Vial number: 2
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.03	0.245	119.17	11.89	2323.831	8.09
2	8.63	0.224	118.33	11.81	2097.231	7.30
3	11.08	0.247	241.83	24.13	4727.177	16.46
4	14.42	0.369	104.96	10.47	3079.617	10.72
5	17.11	0.463	117.35	11.71	4291.931	14.94
6	22.28	0.590	115.03	11.48	4989.561	17.37
7	24.86	0.533	185.45	18.51	7212.237	25.11
7	28.00	0.382	1002.11	100.00	28721.585	100.00

This report has been created by IC Net
METROHM LTD

Report date: 4/6/2018 11:38:10 AM
Printed by: wet


Ident: CCB
Analysis from: 4/5/2018 4:18:03 PM
File: _2018-04-05_

Last save: 4/6/2018 11:34:22 AM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17039

Last save: 4/4/2018 2:46:4

SAMPLE:
: AK/AP
Vial number: 3
Volume: 20.0 μ L
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net
METROHM LTD

Report date: 4/6/2018 11:38:21 AM
Printed by: wet

Ident: LB94832BLW
Analysis from: 4/5/2018 4:48:59 PM
File: _2018-04-05_


AK
4-6-18

Last save: 4/6/2018 11:31:54 AM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17040

Last save: 4/4/2018 2:46:4

SAMPLE:
: AK/AP
Vial number: 4
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net
METROHM LTD

Report date: 4/6/2018 11:38:30 AM
Printed by: wet


Ident: LB94732BSW
Analysis from: 4/5/2018 5:19:55 PM
File: 2018-04-05

Last save: 4/6/2018 11:31:54 AM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17041

Last save: 4/4/2018 2:46:4

SAMPLE:
: AK/AP
Vial number: 5
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.03	0.246	114.26	11.81	2245.232	8.03
2	8.64	0.225	112.84	11.66	2007.381	7.17
3	11.09	0.248	231.54	23.93	4555.439	16.28
4	14.44	0.369	100.56	10.39	2951.556	10.55
5	17.12	0.461	113.09	11.69	4126.071	14.75
6	22.29	0.597	111.97	11.57	4900.969	17.52
7	24.88	0.539	183.37	18.95	7191.099	25.70
7	28.00	0.384	967.63	100.00	27977.747	100.00

Report date: 4/6/2018 11:38:40 AM
Printed by: wet


Ident: J2133-01
Analysis from: 4/5/2018 5:50:50 PM
File: _2018-04-05_

Last save: 4/6/2018 11:36:08 AM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17042

Last save: 4/4/2018 2:46:4

SAMPLE:
: AK/AP
Vial number: 6
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.03	0.268	2.46	7.80	53.400	5.10
2	8.65	0.227	3.62	11.46	67.263	6.42
3	11.08	0.246	5.10	16.15	86.882	8.29
4	14.55	0.368	2.81	8.91	81.928	7.82
5	17.32	0.436	3.47	10.98	130.657	12.47
6	22.41	0.681	8.27	26.18	411.925	39.31
7	24.92	0.541	5.84	18.50	215.774	20.59
7	28.00	0.395	31.57	99.97	1047.830	100.00

This report has been created by IC Net
METROHM LTD

Report date: 4/6/2018 11:38:50 AM
Printed by: wet


Ident: J2133-03
Analysis from: 4/5/2018 6:21:47 PM
File: _2018-04-05_

Last save: 4/6/2018 11:36:08 AM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17043

Last save: 4/4/2018 2:46:4

SAMPLE:
: AK/AP
Vial number: 7
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.02	0.262	2.39	7.59	47.181	4.64
2	8.64	0.218	3.65	11.60	57.087	5.61
3	11.07	0.249	5.05	16.03	89.974	8.84
4	14.54	0.360	2.78	8.81	74.688	7.34
5	17.30	0.428	3.48	11.05	118.206	11.62
6	22.40	0.675	8.27	26.26	408.189	40.13
7	24.91	0.544	5.87	18.64	221.965	21.82
7	28.00	0.391	31.49	99.98	1017.289	100.00

Report date: 4/6/2018 11:39:00 AM
Printed by: wet


Ident: J2133-02
Analysis from: 4/5/2018 6:52:43 PM
File: _2018-04-05_

Last save: 4/5/2018 7:20:43 PM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17044

Last save: 4/4/2018 2:46:4

SAMPLE:
: AK/AP
Vial number: 8
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.02	0.251	10.75	10.07	213.339	6.42
2	8.65	0.226	11.61	10.87	207.436	6.24
3	11.07	0.254	21.33	19.97	410.487	12.35
4	14.52	0.366	10.45	9.79	305.564	9.20
5	17.27	0.433	12.20	11.43	427.954	12.88
6	22.38	0.652	22.07	20.67	1057.499	31.82
7	24.91	0.545	18.36	17.20	700.712	21.09
7	28.00	0.390	106.76	100.00	3322.991	100.00

Report date: 4/6/2018 11:39:23 AM
Printed by: wet


Ident: CCV
Analysis from: 4/5/2018 8:25:36 PM
File: _2018-04-05_

Last save: 4/5/2018 8:53:37 PM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17047

Last save: 4/4/2018 2:46:4

SAMPLE:
: AK/AP
Vial number: 2
Volume: 20.0 µL
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No	Retention min	Width/2 min	Height mV	Height %	Area mV*sec	Area %
1	5.04	0.246	120.04	11.97	2387.018	8.26
2	8.64	0.224	119.14	11.88	2104.283	7.28
3	11.09	0.248	243.72	24.31	4767.883	16.49
4	14.43	0.368	105.91	10.56	3091.051	10.69
5	17.11	0.460	118.48	11.82	4309.719	14.91
6	22.29	0.606	113.68	11.34	5033.784	17.41
7	24.88	0.548	181.68	18.12	7214.449	24.96
7	28.00	0.386	1002.65	100.00	28908.186	100.00

Report date: 4/6/2018 11:39:42 AM
Printed by: wet


Ident: CCB
Analysis from: 4/5/2018 8:56:32 PM
File: _2018-04-05_

Last save: 4/6/2018 11:39:39 AM

Method: AnionsIC2-040418.mtw
Run operator: wet
Analysis number: 17048

Last save: 4/4/2018 2:46:4

SAMPLE:
: AK/AP
Vial number: 3
Volume: 20.0 μ L
Dilution: 1.00
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net
METROHM LTD