

SOP ID : M1312-SPLP--8

SDG No : N/A

Weigh By : JP

Balance ID : WC SC-4

pH Meter ID : WC PH METER-1

Extraction By : JP

Filter By : JP

Pipette ID : WC

Tumbler ID : T-1

TCLP Filter ID : 113075

Start Prep Date : 04/01/2021 Time : 15:00

End Prep Date : 04/02/2021 Time : 08:25

Combination Ratio : 20

ZHE Cleaning Batch : N/A

Initial Room Temperature: 22 °C

Final Room Temperature: 22 °C

TCLP Technician Signature : JP

Supervisor By : AD

Standardized Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE US	Lot Number
SPLP FLUID	N/A	WP90237
N/A	N/A	N/A
N/A	N/A	N/A
pH Strips	N/A	W1931,W1934,W2755,W2350
pH Strips	N/A	W1937,W1938,W1940,W1941,W1942
N/A	N/A	N/A
DI WATER	PH 5.03	W2606
N/A	N/A	N/A

Extraction Conformance/Non-Conformance Comments:

Marix spikes are added after filtration and before preservation. m1794-04 is used for MS-MSD. Tumbler T-1 CHECKED,30 RPM.

Date / Time	Received By	Relinquished By	Location
04/02/21 9:20	KS	JP	Met: -
	Analysis Group	Preparation Group	Ext: -