

SOP ID : M1311-TCLP-13
 SDG No : N/A
 Weigh By : JP
 Balance ID : WC SC-4
 pH Meter ID : WC PH METER-1
 Extraction By : JP
 Filter By : RS
 Pipette ID : WC
 Tumbler ID : T-1
 TCLP Filter ID : 113075

Start Prep Date : 08/19/2021 Time : 17:30
 End Prep Date : 08/20/2021 Time : 09:35
 Combination Ratio : 20
 ZHE Cleaning Batch : N/A
 Initial Room Temperature : 24 °C
 Final Room Temperature : 23 °C
 TCLP Technician Signature : JP
 Supervisor By : AS

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE US	Lot Number
TCLP-FLUID-1	N/A	WP92122
HCL-TCLP,1N	N/A	WP92125
HNO3-TCLP,1N	N/A	WP92126
pH Strips	N/A	W1931,W1934,W2755,W2350
pH Strips	N/A	W1937,W1938,W1939,W1940,W1941,W1942
1L Amber bottle	N/A	19948
120ml Plastic bottle	N/A	210210
1:1 HNO3	N/A	N/A

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. m3489-19 is used for MS-MSD. Tumbler T-1 checked, 30 rpm.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
08/20/21 11:00	JP / 1124 Dom	RJ / Ext
	Preparation Group	Analysis Group SA / met D'y