

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : AS

Instrument ID : Konelab

3/26/2019 7:13

Test: ~~Total~~ CN

Reactive
myh

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	96.289	0.0	0.088	
ICB1	0.399	0.0	0.004	
CCV1	241.030	0.0	0.214	
CCB1	0.051	0.0	0.004	
PB118187BL	0.306	0.0	0.004	
K2040-02	0.736	0.0	0.004	
K2040-02DUP	0.294	0.0	0.004	
K2040-04	0.505	0.0	0.004	
K2040-06	0.407	0.0	0.004	
K2050-07	0.575	0.0	0.004	
K2050-08	0.038	0.0	0.004	
K2050-09	0.639	0.0	0.004	
K2050-10	0.369	0.0	0.004	
K2050-11	0.686	0.0	0.004	
CCV2	240.482	0.0	0.213	
CCB2	0.189	0.0	0.004	
PB118248BL	0.005	0.0	0.004	
K2058-02	-0.232	0.0	0.004	
K2058-02DUP	-0.316	0.0	0.003	
K2058-04	0.185	0.0	0.004	
K2061-02	0.027	0.0	0.004	
K2071-04	0.281	0.0	0.004	
K2071-06	0.326	0.0	0.004	
K2071-03	-0.153	0.0	0.004	
CCV3	255.460	0.0	0.226	
CCB3	0.201	0.0	0.004	

N	26
Mean	32.261
SD	80.8461
CV%	250.60

Aquakem v. 7.2AQ1

Results from time period:

Mon Mar 25 15:38:41 2019

Mon Mar 25 16:06:20 2019

AS

Sample Id	Sam/Ctr/c/	Test short nar	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive CN	P	0.0096	µg/l	3/25/2019 8:17:11	
5.OPPBCN	A	Reactive CN	P	4.2328	µg/l	3/25/2019 8:17:12	
10PPBCN	A	Reactive CN	P	10.0929	µg/l	3/25/2019 8:17:13	
50PPBCN	A	Reactive CN	P	50.1604	µg/l	3/25/2019 8:17:14	
100PPBCN	A	Reactive CN	P	101.5942	µg/l	3/25/2019 8:17:15	
250PPBCN	A	Reactive CN	P	248.4784	µg/l	3/25/2019 8:17:16	
500PPBCN	A	Reactive CN	P	500.4318	µg/l	3/25/2019 8:17:17	
ICV1	S	Reactive CN	P	96.289	µg/l	3/25/2019 15:38:41	
ICB1	S	Reactive CN	P	0.3988	µg/l	3/25/2019 15:38:42	
CCV1	S	Reactive CN	P	241.0302	µg/l	3/25/2019 15:38:43	
CCB1	S	Reactive CN	P	0.0509	µg/l	3/25/2019 15:38:44	
PB118187BL	S	Reactive CN	P	0.3058	µg/l	3/25/2019 15:38:45	
K2040-02	S	Reactive CN	P	0.7361	µg/l	3/25/2019 15:38:46	
K2040-02DUP	S	Reactive CN	P	0.2945	µg/l	3/25/2019 15:38:47	
K2040-04	S	Reactive CN	P	0.5046	µg/l	3/25/2019 15:38:48	
K2040-06	S	Reactive CN	P	0.4065	µg/l	3/25/2019 15:38:49	
K2050-07	S	Reactive CN	P	0.5751	µg/l	3/25/2019 15:38:50	
K2050-08	S	Reactive CN	P	0.0382	µg/l	3/25/2019 15:38:51	
K2050-09	S	Reactive CN	P	0.6391	µg/l	3/25/2019 15:46:14	
K2050-10	S	Reactive CN	P	0.3688	µg/l	3/25/2019 15:46:15	
K2050-11	S	Reactive CN	P	0.6864	µg/l	3/25/2019 15:46:16	
CCV2	S	Reactive CN	P	240.4825	µg/l	3/25/2019 15:46:17	
CCB2	S	Reactive CN	P	0.1888	µg/l	3/25/2019 15:46:18	
PB118248BL	S	Reactive CN	P	0.0054	µg/l	3/25/2019 15:46:19	
K2058-02	S	Reactive CN	P	-0.2316	µg/l	3/25/2019 15:46:20	
K2058-02DUP	S	Reactive CN	P	-0.3163	µg/l	3/25/2019 15:46:21	
K2058-04	S	Reactive CN	P	0.185	µg/l	3/25/2019 15:46:22	
K2061-02	S	Reactive CN	P	0.0271	µg/l	3/25/2019 15:46:23	
K2071-04	S	Reactive CN	P	0.2815	µg/l	3/25/2019 15:50:33	
K2071-06	S	Reactive CN	P	0.3261	µg/l	3/25/2019 15:50:34	
K2071-03	S	Reactive CN	P	-0.1531	µg/l	3/25/2019 16:06:18	
CCV3	S	Reactive CN	P	255.4602	µg/l	3/25/2019 16:06:19	
CCB3	S	Reactive CN	P	0.2007	µg/l	3/25/2019 16:06:20	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

3/25/2019 15:31


Test Total CN

Accepted 3/25/2019 8:19

Factor 1149
Bias 0.004

Coeff. of det. 0.999972

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.004	0.0096	0.0000	
2	5.0PPBCN	0.007	4.2328	5.0000	
3	10PPBCN	0.012	10.0929	10.0000	
4	50PPBCN	0.047	50.1604	50.0000	
5	100PPBCN	0.092	101.5942	100.0000	
6	250PPBCN	0.220	248.4784	250.0000	
7	500PPBCN	0.439	500.4318	500.0000	