

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

4/3/2019 13:58

Test: ~~Total~~ ^{Reactive} CN

mg/d

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	98.460	0.0	0.093	
ICB1	1.568	0.0	0.005	
CCV1	246.482	0.0	0.229	
CCB1	1.361	0.0	0.005	
PB118582BL	1.436	0.0	0.005	
K2197-02	1.155	0.0	0.004	
K2197-02DUP	0.923	0.0	0.004	
K2197-04	1.352	0.0	0.004	
K2199-02	1.324	0.0	0.004	
K2199-04	1.419	0.0	0.005	
K2203-02	0.972	0.0	0.004	
K2203-04	1.142	0.0	0.004	
K2214-04	1.064	0.0	0.004	
K2214-08	1.211	0.0	0.004	
CCV2	246.389	0.0	0.229	
CCB2	1.461	0.0	0.005	
K2221-02	1.085	0.0	0.004	
K2221-04	0.842	0.0	0.004	
CCV3	246.255	0.0	0.229	
CCB3	1.465	0.0	0.005	

N	20
Mean	42.868
SD	90.3409
CV%	210.74

Aquakem v. 7.2AQ1

Results from time period:

Wed Apr 03 13:42:26 2019

Wed Apr 03 13:49:42 2019

AS

Sample Id	Sam/Ctr/c	Test short nan	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive CN	P	0.7302	µg/l	4/3/2019 8:03:09	
5.OPPBCN	A	Reactive CN	P	5.4303	µg/l	4/3/2019 8:03:10	
10PPBCN	A	Reactive CN	P	9.2007	µg/l	4/3/2019 8:03:11	
50PPBCN	A	Reactive CN	P	49.2089	µg/l	4/3/2019 8:03:12	
100PPBCN	A	Reactive CN	P	99.9098	µg/l	4/3/2019 8:03:13	
250PPBCN	A	Reactive CN	P	250.8227	µg/l	4/3/2019 8:03:14	
500PPBCN	A	Reactive CN	P	499.6975	µg/l	4/3/2019 8:03:15	
ICV1	S	Reactive CN	P	98.4599	µg/l	4/3/2019 13:42:26	
ICB1	S	Reactive CN	P	1.5679	µg/l	4/3/2019 13:42:27	
CCV1	S	Reactive CN	P	246.4821	µg/l	4/3/2019 13:42:28	
CCB1	S	Reactive CN	P	1.3612	µg/l	4/3/2019 13:42:29	
PB118582BL	S	Reactive CN	P	1.4359	µg/l	4/3/2019 13:42:30	
K2197-02	S	Reactive CN	P	1.1548	µg/l	4/3/2019 13:42:31	
K2197-02DUP	S	Reactive CN	P	0.9235	µg/l	4/3/2019 13:42:32	
K2197-04	S	Reactive CN	P	1.3519	µg/l	4/3/2019 13:42:33	
K2199-02	S	Reactive CN	P	1.3244	µg/l	4/3/2019 13:42:34	
K2199-04	S	Reactive CN	P	1.4188	µg/l	4/3/2019 13:42:35	
K2203-02	S	Reactive CN	P	0.9717	µg/l	4/3/2019 13:42:36	
K2203-04	S	Reactive CN	P	1.1417	µg/l	4/3/2019 13:49:33	
K2214-04	S	Reactive CN	P	1.0639	µg/l	4/3/2019 13:49:34	
K2214-08	S	Reactive CN	P	1.211	µg/l	4/3/2019 13:49:35	
CCV2	S	Reactive CN	P	246.3885	µg/l	4/3/2019 13:49:36	
CCB2	S	Reactive CN	P	1.4613	µg/l	4/3/2019 13:49:37	
K2221-02	S	Reactive CN	P	1.0853	µg/l	4/3/2019 13:49:38	
K2221-04	S	Reactive CN	P	0.8419	µg/l	4/3/2019 13:49:39	
CCV3	S	Reactive CN	P	246.255	µg/l	4/3/2019 13:49:41	
CCB3	S	Reactive CN	P	1.4654	µg/l	4/3/2019 13:49:42	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

4/3/2019 9:24


Test Total CN

Accepted 4/3/2019 9:21

Factor 1092
Bias 0.003

Coeff. of det. 0.999987

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.004	0.7302	0.0000	
2	5.0PPBCN	0.008	5.4303	5.0000	
3	10PPBCN	0.012	9.2007	10.0000	
4	50PPBCN	0.048	49.2089	50.0000	
5	100PPBCN	0.095	99.9098	100.0000	
6	250PPBCN	0.233	250.8227	250.0000	
7	500PPBCN	0.461	499.6975	500.0000	