

LB101966

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

4/9/2019 13:42

Test: *Reactive*
~~Total~~ CN

Sample Id	Result <i>Ng/L</i>	Dil. 1 +	Response	Errors
ICV1	99.380	0.0	0.094	
ICB1	0.614	0.0	0.004	
CCV1	249.852	0.0	0.230	
CCB1	0.201	0.0	0.004	
PB118730BL	0.271	0.0	0.004	
K2198-02	0.152	0.0	0.004	
K2198-02DUP	-0.099	0.0	0.004	
K2199-06	0.065	0.0	0.004	
K2199-08	0.123	0.0	0.004	
K2334-03	0.044	0.0	0.004	
K2337-04	-0.108	0.0	0.004	
K2338-02	0.139	0.0	0.004	
K2338-04	0.159	0.0	0.004	
K2338-09	0.445	0.0	0.004	
CCV2	248.504	0.0	0.229	
CCB2	0.328	0.0	0.004	
K2338-11	0.295	0.0	0.004	
K2338-13	0.214	0.0	0.004	
PB118732BL	0.075	0.0	0.004	
K2337-02	0.245	0.0	0.004	
K2337-02DUP	0.108	0.0	0.004	
K2338-06	-0.031	0.0	0.004	
CCV3	263.625	0.0	0.242	
CCB3	0.485	0.0	0.004	

N	24
Mean	36.045
SD	86.5709
CV%	240.17

Aquakem v. 7.2AQ1

Results from time period:

Tue Apr 09 12:59:20 2019

Tue Apr 09 13:41:23 2019

✓AS

Sample Id	Sam/Ctr/c/	Test short narr	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive CN	P	0.2064	µg/l	4/9/2019 12:59:20	
5.OPPBCN	A	Reactive CN	P	4.7473	µg/l	4/9/2019 12:59:21	
10PPBCN	A	Reactive CN	P	10.2037	µg/l	4/9/2019 12:59:22	
50PPBCN	A	Reactive CN	P	50.5317	µg/l	4/9/2019 12:59:23	
100PPBCN	A	Reactive CN	P	101.0259	µg/l	4/9/2019 12:59:24	
250PPBCN	A	Reactive CN	P	247.0896	µg/l	4/9/2019 12:59:25	
500PPBCN	A	Reactive CN	P	501.1953	µg/l	4/9/2019 12:59:26	
ICV1	S	Reactive CN	P	99.3797	µg/l	4/9/2019 13:29:55	
ICB1	S	Reactive CN	P	0.6141	µg/l	4/9/2019 13:29:56	
CCV1	S	Reactive CN	P	249.8521	µg/l	4/9/2019 13:29:57	
CCB1	S	Reactive CN	P	0.2014	µg/l	4/9/2019 13:29:58	
PB118730BL	S	Reactive CN	P	0.2712	µg/l	4/9/2019 13:29:59	
K2198-02	S	Reactive CN	P	0.1518	µg/l	4/9/2019 13:30:00	
K2198-02DUP	S	Reactive CN	P	-0.0994	µg/l	4/9/2019 13:30:01	
K2199-06	S	Reactive CN	P	0.0649	µg/l	4/9/2019 13:30:02	
K2199-08	S	Reactive CN	P	0.1229	µg/l	4/9/2019 13:30:03	
K2334-03	S	Reactive CN	P	0.0437	µg/l	4/9/2019 13:30:04	
K2337-04	S	Reactive CN	P	-0.1084	µg/l	4/9/2019 13:30:05	
K2338-02	S	Reactive CN	P	0.1385	µg/l	4/9/2019 13:37:30	
K2338-04	S	Reactive CN	P	0.159	µg/l	4/9/2019 13:37:31	
K2338-09	S	Reactive CN	P	0.4455	µg/l	4/9/2019 13:37:32	
CCV2	S	Reactive CN	P	248.5036	µg/l	4/9/2019 13:37:33	
CCB2	S	Reactive CN	P	0.3282	µg/l	4/9/2019 13:37:34	
K2338-11	S	Reactive CN	P	0.295	µg/l	4/9/2019 13:37:35	
K2338-13	S	Reactive CN	P	0.2136	µg/l	4/9/2019 13:37:36	
PB118732BL	S	Reactive CN	P	0.0753	µg/l	4/9/2019 13:37:37	
K2337-02	S	Reactive CN	P	0.2447	µg/l	4/9/2019 13:37:38	
K2337-02DUP	S	Reactive CN	P	0.1083	µg/l	4/9/2019 13:37:39	
K2338-06	S	Reactive CN	P	-0.0309	µg/l	4/9/2019 13:37:40	
CCV3	S	Reactive CN	P	263.6252	µg/l	4/9/2019 13:41:22	
CCB3	S	Reactive CN	P	0.4853	µg/l	4/9/2019 13:41:23	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

4/9/2019 13:01


Test Total CN

Accepted 4/9/2019 13:00

Factor 1105
Bias 0.004

Coeff. of det. 0.999945

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.004	0.2064	0.0000	
2	5.0PPBCN	0.008	4.7473	5.0000	
3	10PPBCN	0.013	10.2037	10.0000	
4	50PPBCN	0.050	50.5317	50.0000	
5	100PPBCN	0.095	101.0259	100.0000	
6	250PPBCN	0.227	247.0896	250.0000	
7	500PPBCN	0.457	501.1953	500.0000	