

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

6/3/2019 14:31

Test: Total CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	100.683	0.0	0.102	
ICB1	0.109	0.0	0.004	
CCV1	238.910	0.0	0.237	
CCB1	0.081	0.0	0.004	
PB120282BL	-0.116	0.0	0.004	
K2641-14	-0.002	0.0	0.004	
K2641-14DUP	-0.246	0.0	0.003	
K2641-16	-0.308	0.0	0.003	
K2641-18	-0.286	0.0	0.003	
K2641-20	-0.256	0.0	0.003	
K2641-22	-0.421	0.0	0.003	
K3135-02	-0.037	0.0	0.004	
K3138-11	-0.416	0.0	0.003	
K3138-12	-0.477	0.0	0.003	
CCV2	239.430	0.0	0.237	
CCB2	-0.092	0.0	0.004	
K3138-13	-0.187	0.0	0.003	
K3138-14	0.312	0.0	0.004	
K3138-15	0.053	0.0	0.004	
K3139-04	0.489	0.0	0.004	
K3139-08	-0.257	0.0	0.003	
CCV3	241.588	0.0	0.240	
CCB3	0.216	0.0	0.004	

N	23
Mean	35.599
SD	83.6010
CV%	234.84

Aquakem v. 7.2AQ1

Results from time period:

Mon Jun 03 14:08:29 2019

Mon Jun 03 14:31:02 2019

✓ AS

Sample Id	Sam/Ctr/c/	Test short na	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive CN	P	-0.2019	µg/l	6/3/2019 11:48:35	
5.OPPBCN	A	Reactive CN	P	4.3799	µg/l	6/3/2019 11:48:36	
10PPBCN	A	Reactive CN	P	9.6099	µg/l	6/3/2019 11:48:37	
50PPBCN	A	Reactive CN	P	49.9935	µg/l	6/3/2019 11:48:38	
100PPBCN	A	Reactive CN	P	101.2708	µg/l	6/3/2019 11:48:39	
250PPBCN	A	Reactive CN	P	250.3746	µg/l	6/3/2019 11:48:40	
500PPBCN	A	Reactive CN	P	499.5732	µg/l	6/3/2019 11:48:41	
ICV1	S	Reactive CN	P	100.6834	µg/l	6/3/2019 14:08:29	
ICB1	S	Reactive CN	P	0.1091	µg/l	6/3/2019 14:08:30	
CCV1	S	Reactive CN	P	238.9102	µg/l	6/3/2019 14:08:31	
CCB1	S	Reactive CN	P	0.0814	µg/l	6/3/2019 14:08:32	
PB120282BL	S	Reactive CN	P	-0.1157	µg/l	6/3/2019 14:08:33	
K2641-14	S	Reactive CN	P	-0.0022	µg/l	6/3/2019 14:08:34	
K2641-14DUP	S	Reactive CN	P	-0.2456	µg/l	6/3/2019 14:08:35	
K2641-16	S	Reactive CN	P	-0.308	µg/l	6/3/2019 14:08:36	
K2641-18	S	Reactive CN	P	-0.286	µg/l	6/3/2019 14:08:37	
K2641-20	S	Reactive CN	P	-0.2564	µg/l	6/3/2019 14:08:38	
K2641-22	S	Reactive CN	P	-0.421	µg/l	6/3/2019 14:08:39	
K3135-02	S	Reactive CN	P	-0.0366	µg/l	6/3/2019 14:12:20	
K3138-11	S	Reactive CN	P	-0.4158	µg/l	6/3/2019 14:12:21	
K3138-12	S	Reactive CN	P	-0.4771	µg/l	6/3/2019 14:12:22	
CCV2	S	Reactive CN	P	239.4296	µg/l	6/3/2019 14:30:53	
CCB2	S	Reactive CN	P	-0.0921	µg/l	6/3/2019 14:30:54	
K3138-13	S	Reactive CN	P	-0.1866	µg/l	6/3/2019 14:30:55	
K3138-14	S	Reactive CN	P	0.3118	µg/l	6/3/2019 14:30:56	
K3138-15	S	Reactive CN	P	0.0534	µg/l	6/3/2019 14:30:57	
K3139-04	S	Reactive CN	P	0.4887	µg/l	6/3/2019 14:30:58	
K3139-08	S	Reactive CN	P	-0.2572	µg/l	6/3/2019 14:30:59	
CCV3	S	Reactive CN	P	241.5885	µg/l	6/3/2019 14:31:01	
CCB3	S	Reactive CN	P	0.2155	µg/l	6/3/2019 14:31:02	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

6/3/2019 12:31


Test Total CN

Accepted 6/3/2019 12:31

Factor 1024
Bias 0.004

Coeff. of det. 0.999988

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.003	-0.2019	0.0000	
2	5.0PPBCN	0.008	4.3799	5.0000	
3	10PPBCN	0.013	9.6099	10.0000	
4	50PPBCN	0.052	49.9935	50.0000	
5	100PPBCN	0.103	101.2708	100.0000	
6	250PPBCN	0.248	250.3746	250.0000	
7	500PPBCN	0.491	499.5732	500.0000	