

LB106203

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : AK

Instrument ID : Konelab

11/11/2019 16:41

Test: *Reactive*
~~Total~~ CN *WU*

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	94.916	0.0	0.086	
ICB1	1.331	0.0	0.003	
CCV1	243.703	0.0	0.218	
CCB1	0.295	0.0	0.003	
PB124697BL	0.017	0.0	0.002	
K5673-08	0.069	0.0	0.002	
K5673-08DUP	0.224	0.0	0.002	
K5676-15	-0.096	0.0	0.002	
K5746-07	-0.008	0.0	0.002	
K5758-02 c3	0.055	0.0	0.002	
K5760-02	0.086	0.0	0.002	
K5760-08	-0.032	0.0	0.002	
K5760-13	0.024	0.0	0.002	
K5763-04	0.431	0.0	0.003	
CCV2	243.846	0.0	0.218	
CCB2	0.178	0.0	0.002	
K5763-08	-0.209	0.0	0.002	
K5763-12	-0.012	0.0	0.002	
K5778-04	-0.062	0.0	0.002	
K5778-08	0.075	0.0	0.002	
CCV3	244.672	0.0	0.219	
CCB3	0.324	0.0	0.003	

N 22
 Mean 37.719
 SD 86.3070
 CV% 228.81

AK
11/12/19

Aquakem v. 7.2AQ1

Results from time period:

Mon Nov 11 15:33:07 2019

Mon Nov 11 16:35:42 2019

Sample Id	Sam/Ctr/c/	Test short Test type	Result	Result unit	Result date and time	Stat
0.0PPBCN	A	Reactive Cl P	-0.2441	µg/l	11/11/2019 13:01:05	
5.0PPBCN	A	Reactive Cl P	3.589	µg/l	11/11/2019 13:01:06	
10PPBCN	A	Reactive Cl P	9.8897	µg/l	11/11/2019 13:01:07	
50PPBCN	A	Reactive Cl P	49.2409	µg/l	11/11/2019 13:01:08	
100PPBCN	A	Reactive Cl P	100.2388	µg/l	11/11/2019 13:01:09	
250PPBCN	A	Reactive Cl P	254.4825	µg/l	11/11/2019 13:01:10	
500PPBCN	A	Reactive Cl P	497.8032	µg/l	11/11/2019 13:01:11	
ICV1	S	Reactive Cl P	94.916	µg/l	11/11/2019 15:33:07	
ICB1	S	Reactive Cl P	1.3309	µg/l	11/11/2019 15:33:08	
CCV1	S	Reactive Cl P	243.7028	µg/l	11/11/2019 16:27:59	
CCB1	S	Reactive Cl P	0.2955	µg/l	11/11/2019 16:28:00	
PB124697BL	S	Reactive Cl P	0.0167	µg/l	11/11/2019 16:28:01	
K5673-08	S	Reactive Cl P	0.0687	µg/l	11/11/2019 16:28:02	
K5673-08DUP	S	Reactive Cl P	0.2239	µg/l	11/11/2019 16:28:03	
K5676-15	S	Reactive Cl P	-0.0957	µg/l	11/11/2019 16:28:04	
K5746-07	S	Reactive Cl P	-0.0076	µg/l	11/11/2019 16:28:05	
K5758-02	S	Reactive Cl P	0.0547	µg/l	11/11/2019 16:28:06	
K5760-03	S	Reactive Cl P	0.0862	µg/l	11/11/2019 16:28:07	
K5760-08	S	Reactive Cl P	-0.0318	µg/l	11/11/2019 16:35:32	
K5760-13	S	Reactive Cl P	0.0245	µg/l	11/11/2019 16:35:33	
K5763-04	S	Reactive Cl P	0.4311	µg/l	11/11/2019 16:35:34	
CCV2	S	Reactive Cl P	243.8457	µg/l	11/11/2019 16:35:35	
CCB2	S	Reactive Cl P	0.1781	µg/l	11/11/2019 16:35:36	
K5763-08	S	Reactive Cl P	-0.2094	µg/l	11/11/2019 16:35:37	
K5763-12	S	Reactive Cl P	-0.0119	µg/l	11/11/2019 16:35:38	
K5778-04	S	Reactive Cl P	-0.0619	µg/l	11/11/2019 16:35:39	
K5778-08	S	Reactive Cl P	0.0755	µg/l	11/11/2019 16:35:40	
CCV3	S	Reactive Cl P	244.672	µg/l	11/11/2019 16:35:41	
CCB3	S	Reactive Cl P	0.3243	µg/l	11/11/2019 16:35:42	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by :

Instrument ID : Konelab

11/11/2019 13:21


Test Total CN

Accepted 11/11/2019 13:21

Factor 1130
Bias 0.002

Coeff. of det. 0.999866

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.002	-0.2441	0.0000	
2	5.0PPBCN	0.005	3.5890	5.0000	
3	10PPBCN	0.011	9.8897	10.0000	
4	50PPBCN	0.046	49.2409	50.0000	
5	100PPBCN	0.091	100.2388	100.0000	
6	250PPBCN	0.227	254.4825	250.0000	
7	500PPBCN	0.443	497.8032	500.0000	