

UB106493

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AC

Instrument ID : Konelab

11/26/2019 16:02

Test: Reactive
~~Total~~ CN ug/L

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	92.514	0.0	0.099	
ICB1	0.239	0.0	0.002	
CCV1	248.479	0.0	0.263	
CCB1	0.163	0.0	0.002	
PB125120BL	-0.025	0.0	0.002	
K5667-04	-0.134	0.0	0.002	
K5667-04DUP	-0.175	0.0	0.002	
K6010-01	-0.148	0.0	0.002	
K6013-03	-0.233	0.0	0.002	
K6022-04	-0.149	0.0	0.002	
K6022-08	-0.229	0.0	0.002	
K6022-12	-0.286	0.0	0.002	
K6022-16	-0.215	0.0	0.002	
K6025-01	-0.187	0.0	0.002	
CCV2	249.584	0.0	0.264	
CCB2	-0.089	0.0	0.002	
PB125135BL	-0.370	0.0	0.002	
K6008-01	-0.429	0.0	0.002	
K6008-01DUP	-0.496	0.0	0.002	
K6025-02	-0.282	0.0	0.002	
CCV3	246.704	0.0	0.261	
CCB3	-0.006	0.0	0.002	

N 22
Mean 37.920
SD 87.7831
CV% 231.50

Aquakem v. 7.2AQ1

Results from time period:

Tue Nov 26 15:42:50 2019

Tue Nov 26 15:50:35 2019

Sample Id	Sam/Ctr/c	Test short	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive	Cl P	-0.3465	µg/l	11/26/2019 13:28:03	
5.OPPBCN	A	Reactive	Cl P	5.1707	µg/l	11/26/2019 13:28:04	
10PPBCN	A	Reactive	Cl P	10.4045	µg/l	11/26/2019 13:28:05	
50PPBCN	A	Reactive	Cl P	49.7294	µg/l	11/26/2019 13:28:06	
100PPBCN	A	Reactive	Cl P	97.0983	µg/l	11/26/2019 13:28:07	
250PPBCN	A	Reactive	Cl P	254.6921	µg/l	11/26/2019 13:28:08	
500PPBCN	A	Reactive	Cl P	498.2516	µg/l	11/26/2019 13:28:09	
ICV1	S	Reactive	Cl P	92.5144	µg/l	11/26/2019 15:42:50	
ICB1	S	Reactive	Cl P	0.239	µg/l	11/26/2019 15:42:51	
CCV1	S	Reactive	Cl P	248.4792	µg/l	11/26/2019 15:42:52	
CCB1	S	Reactive	Cl P	0.1625	µg/l	11/26/2019 15:42:53	
PB125120BL	S	Reactive	Cl P	-0.0249	µg/l	11/26/2019 15:42:54	
K5667-04	S	Reactive	Cl P	-0.1338	µg/l	11/26/2019 15:42:55	
K5667-04DUP	S	Reactive	Cl P	-0.1754	µg/l	11/26/2019 15:42:56	
K6010-01	S	Reactive	Cl P	-0.1481	µg/l	11/26/2019 15:42:57	
K6013-03	S	Reactive	Cl P	-0.2329	µg/l	11/26/2019 15:42:58	
K6022-04	S	Reactive	Cl P	-0.1494	µg/l	11/26/2019 15:42:59	
K6022-08	S	Reactive	Cl P	-0.2293	µg/l	11/26/2019 15:43:00	
K6022-12	S	Reactive	Cl P	-0.2859	µg/l	11/26/2019 15:50:25	
K6022-16	S	Reactive	Cl P	-0.2149	µg/l	11/26/2019 15:50:26	
K6025-01	S	Reactive	Cl P	-0.1869	µg/l	11/26/2019 15:50:27	
CCV2	S	Reactive	Cl P	249.5839	µg/l	11/26/2019 15:50:28	
CCB2	S	Reactive	Cl P	-0.0892	µg/l	11/26/2019 15:50:29	
PB125135BL	S	Reactive	Cl P	-0.37	µg/l	11/26/2019 15:50:30	
K6008-01	S	Reactive	Cl P	-0.429	µg/l	11/26/2019 15:50:31	
K6008-01DUP	S	Reactive	Cl P	-0.4963	µg/l	11/26/2019 15:50:32	
K6025-02	S	Reactive	Cl P	-0.2819	µg/l	11/26/2019 15:50:33	
CCV3	S	Reactive	Cl P	246.7036	µg/l	11/26/2019 15:50:34	
CCB3	S	Reactive	Cl P	-0.0057	µg/l	11/26/2019 15:50:35	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AR

Instrument ID : Konelab

11/26/2019 16:16

Test Total CN


Accepted 11/26/2019 13:42

Factor 954.1

Bias 0.002

Coeff. of det. 0.999835

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.002	-0.3465	0.0000	
2	5.0PPBCN	0.008	5.1707	5.0000	
3	10PPBCN	0.013	10.4045	10.0000	
4	50PPBCN	0.054	49.7294	50.0000	
5	100PPBCN	0.104	97.0983	100.0000	
6	250PPBCN	0.269	254.6921	250.0000	
7	500PPBCN	0.524	498.2516	500.0000	