

LB106904

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : AK

Instrument ID : Konelab

12/23/2019 15:10

Test: *Reactive*
~~Total~~ CN

Sample Id	<i>ng/L</i> Result	Dil. 1 +	Response	Errors
ICV1	94.672	0.0	0.109	
ICB1	1.141	0.0	0.002	
CCV1	242.405	0.0	0.277	
CCB1	0.856	0.0	0.002	
PB125675BL	0.756	0.0	0.002	
K6117-09	0.764	0.0	0.002	
K6117-09DUP	0.760	0.0	0.002	
K6386-02	0.745	0.0	0.002	
K6395-04	0.843	0.0	0.002	
CCV2	239.856	0.0	0.274	
CCB2	1.199	0.0	0.002	

N	11
Mean	53.091
SD	97.0854
CV%	182.87

Aquakem v. 7.2AQ1

Results from time period:

Mon Dec 23 15:03:34 2019

Mon Dec 23 15:03:44 2019

AR

Sample Id	Sam/Ctr/c/	Test short	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive	Cl P	0.6273	µg/l	12/23/2019 13:25:13	
5.OPPBCN	A	Reactive	Cl P	5.5772	µg/l	12/23/2019 13:25:14	
10PPBCN	A	Reactive	Cl P	10.5037	µg/l	12/23/2019 13:25:15	
50PPBCN	A	Reactive	Cl P	48.3803	µg/l	12/23/2019 13:25:16	
100PPBCN	A	Reactive	Cl P	99.4048	µg/l	12/23/2019 13:25:17	
250PPBCN	A	Reactive	Cl P	250.4829	µg/l	12/23/2019 13:25:18	
500PPBCN	A	Reactive	Cl P	500.0237	µg/l	12/23/2019 13:25:19	
ICV1	S	Reactive	Cl P	94.6721	µg/l	12/23/2019 15:03:34	
ICB1	S	Reactive	Cl P	1.1409	µg/l	12/23/2019 15:03:35	
CCV1	S	Reactive	Cl P	242.4052	µg/l	12/23/2019 15:03:36	
CCB1	S	Reactive	Cl P	0.856	µg/l	12/23/2019 15:03:37	
PB125675BL	S	Reactive	Cl P	0.7564	µg/l	12/23/2019 15:03:38	
K6117-09	S	Reactive	Cl P	0.7639	µg/l	12/23/2019 15:03:39	
K6117-09DUP	S	Reactive	Cl P	0.7598	µg/l	12/23/2019 15:03:40	
K6386-02	S	Reactive	Cl P	0.7453	µg/l	12/23/2019 15:03:41	
K6395-04	S	Reactive	Cl P	0.8434	µg/l	12/23/2019 15:03:42	
CCV2	S	Reactive	Cl P	239.8564	µg/l	12/23/2019 15:03:43	
CCB2	S	Reactive	Cl P	1.1994	µg/l	12/23/2019 15:03:44	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AK

Instrument ID : Konelab

12/23/2019 13:33


Test Total CN

Accepted 12/23/2019 13:33

Factor 878.5
Bias 0.001

Coeff. of det. 0.999980

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.002	0.6273	0.0000	
2	5.0PPBCN	0.007	5.5772	5.0000	
3	10PPBCN	0.013	10.5037	10.0000	
4	50PPBCN	0.056	48.3803	50.0000	
5	100PPBCN	0.114	99.4048	100.0000	
6	250PPBCN	0.286	250.4829	250.0000	
7	500PPBCN	0.570	500.0237	500.0000	