

LB106933

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AR

Instrument ID : Konelab

12/26/2019 12:28

Test: ^{Reactive} ~~Total~~ CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	95.652	0.0	0.097	
ICB1	0.926	0.0	0.003	
CCV1	238.260	0.0	0.239	
CCB1	0.682	0.0	0.002	
PB125710BL	0.641	0.0	0.002	
K6118-13	0.592	0.0	0.002	
K6118-13DUP	0.680	0.0	0.002	
K6413-14	0.697	0.0	0.002	
K6413-15	0.660	0.0	0.002	
K6413-16	0.767	0.0	0.002	
K6414-09	0.607	0.0	0.002	
K6414-10	0.636	0.0	0.002	
K6436-08	0.764	0.0	0.002	
K6436-09	0.895	0.0	0.003	
CCV2	245.202	0.0	0.246	
CCB2	0.665	0.0	0.002	
K6437-12	0.619	0.0	0.002	
K6437-13	0.578	0.0	0.002	
K6437-15	0.581	0.0	0.002	
CCV3	241.566	0.0	0.243	
CCB3	0.654	0.0	0.002	

N 21
Mean 39.634
SD 87.0100
CV% 219.53

K6413
K6413
K6413

AK
1430/19

Aquakem v. 7.2AQ1

Results from time period:

Thu Dec 26 12:18:52 2019

Thu Dec 26 12:26:08 2019

Sample Id	Sam/Ctr/c/	Test short	Test type	Result	Result unit	Result date and time	Stat
0.0PPBCN	A	Reactive	Cl P	0.1662	µg/l	12/26/2019 10:07:30	
5.0PPBCN	A	Reactive	Cl P	5.1777	µg/l	12/26/2019 10:07:31	
10PPBCN	A	Reactive	Cl P	10.1515	µg/l	12/26/2019 10:07:32	
50PPBCN	A	Reactive	Cl P	50.2619	µg/l	12/26/2019 10:07:33	
100PPBCN	A	Reactive	Cl P	101.0775	µg/l	12/26/2019 10:07:34	
250PPBCN	A	Reactive	Cl P	246.8235	µg/l	12/26/2019 10:07:35	
500PPBCN	A	Reactive	Cl P	501.3417	µg/l	12/26/2019 10:07:36	
ICV1	S	Reactive	Cl P	95.6516	µg/l	12/26/2019 12:18:52	
ICB1	S	Reactive	Cl P	0.9264	µg/l	12/26/2019 12:18:53	
CCV1	S	Reactive	Cl P	238.2601	µg/l	12/26/2019 12:18:54	
CCB1	S	Reactive	Cl P	0.682	µg/l	12/26/2019 12:18:55	
PB125710BL	S	Reactive	Cl P	0.6412	µg/l	12/26/2019 12:18:56	
K6118-13	S	Reactive	Cl P	0.5919	µg/l	12/26/2019 12:18:57	
K6118-13DUP	S	Reactive	Cl P	0.6803	µg/l	12/26/2019 12:18:58	
K6413-14	S	Reactive	Cl P	0.6974	µg/l	12/26/2019 12:18:59	
K6413-15	S	Reactive	Cl P	0.6598	µg/l	12/26/2019 12:19:00	
K6413-16	S	Reactive	Cl P	0.7672	µg/l	12/26/2019 12:19:01	
K6414-09	S	Reactive	Cl P	0.607	µg/l	12/26/2019 12:19:02	
K6414-10	S	Reactive	Cl P	0.6357	µg/l	12/26/2019 12:25:59	
K6436-08	S	Reactive	Cl P	0.7637	µg/l	12/26/2019 12:26:00	
K6436-09	S	Reactive	Cl P	0.8948	µg/l	12/26/2019 12:26:01	
CCV2	S	Reactive	Cl P	245.2021	µg/l	12/26/2019 12:26:02	
CCB2	S	Reactive	Cl P	0.665	µg/l	12/26/2019 12:26:03	
K6437-12	S	Reactive	Cl P	0.619	µg/l	12/26/2019 12:26:04	
K6437-13	S	Reactive	Cl P	0.5782	µg/l	12/26/2019 12:26:05	
K6437-15	S	Reactive	Cl P	0.5809	µg/l	12/26/2019 12:26:06	
CCV3	S	Reactive	Cl P	241.5658	µg/l	12/26/2019 12:26:07	
CCB3	S	Reactive	Cl P	0.6536	µg/l	12/26/2019 12:26:08	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AK

Instrument ID : Konelab

12/26/2019 10:10


Test Total CN

Accepted 12/26/2019 10:10

Factor 1003
Bias 0.002

Coeff. of det. 0.999936

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.002	0.1662	0.0000	
2	5.0PPBCN	0.007	5.1777	5.0000	
3	10PPBCN	0.012	10.1515	10.0000	
4	50PPBCN	0.052	50.2619	50.0000	
5	100PPBCN	0.103	101.0775	100.0000	
6	250PPBCN	0.248	246.8235	250.0000	
7	500PPBCN	0.502	501.3417	500.0000	