

UB106965

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : AK

Instrument ID : Konelab

12/30/2019 12:56

Test: ~~Total~~ ^{Reactive} CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	95.446	0.0	0.092	
ICB1	1.320	0.0	0.002	
CCV1	239.787	0.0	0.229	
CCB1	1.165	0.0	0.002	
PB125766BL	1.043	0.0	0.002	
K6113-15	0.908	0.0	0.002	
K6113-15DUP	1.024	0.0	0.002	
K6117-12	0.985	0.0	0.002	
K6441-08	1.033	0.0	0.002	
K6441-09	1.133	0.0	0.002	
K6442-07	0.995	0.0	0.002	
PB125790BL	0.863	0.0	0.002	
K6114-02	0.958	0.0	0.002	
K6114-02DUP	1.093	0.0	0.002	
CCV2	243.881	0.0	0.233	
CCB2	1.034	0.0	0.002	
K6455-07	0.909	0.0	0.002	
K6456-02	0.946	0.0	0.002	
K6456-04	1.041	0.0	0.002	
CCV3	248.002	0.0	0.237	
CCB3	1.300	0.0	0.002	

N	21
Mean	40.232
SD	87.6411
CV%	217.84

Aquakem v. 7.2AQ1

AK

Results from time period:

Mon Dec 30 12:16:07 2019

Mon Dec 30 12:23:23 2019

Sample Id	Sam/Ctr/c/	Test short	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive	Cl P	0.7527	µg/l	12/30/2019 8:57:03	
5.OPPBCN	A	Reactive	Cl P	5.5156	µg/l	12/30/2019 8:57:04	
10PPBCN	A	Reactive	Cl P	10.3731	µg/l	12/30/2019 8:57:05	
50PPBCN	A	Reactive	Cl P	49.7973	µg/l	12/30/2019 8:57:06	
100PPBCN	A	Reactive	Cl P	98.1742	µg/l	12/30/2019 8:57:07	
250PPBCN	A	Reactive	Cl P	250.0288	µg/l	12/30/2019 8:57:08	
500PPBCN	A	Reactive	Cl P	500.3585	µg/l	12/30/2019 8:57:09	
ICV1	S	Reactive	Cl P	95.4458	µg/l	12/30/2019 12:16:07	
ICB1	S	Reactive	Cl P	1.3197	µg/l	12/30/2019 12:16:08	
CCV1	S	Reactive	Cl P	239.787	µg/l	12/30/2019 12:16:09	
CCB1	S	Reactive	Cl P	1.1646	µg/l	12/30/2019 12:16:10	
PB125766BL	S	Reactive	Cl P	1.0429	µg/l	12/30/2019 12:16:11	
K6113-15	S	Reactive	Cl P	0.9085	µg/l	12/30/2019 12:16:12	
K6113-15DUP	S	Reactive	Cl P	1.0235	µg/l	12/30/2019 12:16:13	
K6117-12	S	Reactive	Cl P	0.9846	µg/l	12/30/2019 12:16:14	
K6441-08	S	Reactive	Cl P	1.0327	µg/l	12/30/2019 12:16:15	
K6441-09	S	Reactive	Cl P	1.1334	µg/l	12/30/2019 12:16:16	
K6442-07	S	Reactive	Cl P	0.9954	µg/l	12/30/2019 12:16:17	
PB125790BL	S	Reactive	Cl P	0.8633	µg/l	12/30/2019 12:23:14	
K6114-02	S	Reactive	Cl P	0.9582	µg/l	12/30/2019 12:23:15	
K6114-02DUP	S	Reactive	Cl P	1.0932	µg/l	12/30/2019 12:23:16	
CCV2	S	Reactive	Cl P	243.8814	µg/l	12/30/2019 12:23:17	
CCB2	S	Reactive	Cl P	1.0335	µg/l	12/30/2019 12:23:18	
K6455-07	S	Reactive	Cl P	0.9089	µg/l	12/30/2019 12:23:19	
K6456-02	S	Reactive	Cl P	0.9456	µg/l	12/30/2019 12:23:20	
K6456-04	S	Reactive	Cl P	1.0414	µg/l	12/30/2019 12:23:21	
CCV3	S	Reactive	Cl P	248.0023	µg/l	12/30/2019 12:23:22	
CCB3	S	Reactive	Cl P	1.3005	µg/l	12/30/2019 12:23:23	

CHEMTECH
 284 Sheffield Street,
 Mountainside, NJ 07092
 Reviewed by : AK

Instrument ID : Konelab

12/30/2019 9:08


Test Total CN

Accepted 12/30/2019 9:08

Factor 1051
 Bias 0.001

Coeff. of det. 0.999978

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.002	0.7527	0.0000	
2	5.0PPBCN	0.006	5.5156	5.0000	
3	10PPBCN	0.011	10.3731	10.0000	
4	50PPBCN	0.048	49.7973	50.0000	
5	100PPBCN	0.094	98.1742	100.0000	
6	250PPBCN	0.239	250.0288	250.0000	
7	500PPBCN	0.477	500.3585	500.0000	