

LB107281

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : AK

Instrument ID : Konelab

1/21/2020 13:58

Test: Total CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	95.116	0.0	0.096	
ICB1	0.322	0.0	0.003	
CCV1	240.979	0.0	0.239	
CCB1	0.296	0.0	0.003	
PB126248BL	0.201	0.0	0.002	
L1004-15	0.121	0.0	0.002	
L1004-15DUP	0.247	0.0	0.002	
L1210-07	0.168	0.0	0.002	
L1211-02	0.135	0.0	0.002	
L1211-04	0.028	0.0	0.002	
CCV2	238.702	0.0	0.237	
CCB2	0.059	0.0	0.002	

N	12
Mean	48.031
SD	93.6208
CV%	194.92

AK

Aquakem v. 7.2AQ1

Results from time period:

Tue Jan 21 13:53:26 2020

Tue Jan 21 13:56:14 2020

Sample Id	Sam/Ctr/c/	Test short Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive Cl P	-0.1993	µg/l	1/21/2020 12:03:11	
5.OPPBCN	A	Reactive Cl P	4.8217	µg/l	1/21/2020 12:03:12	
10PPBCN	A	Reactive Cl P	10.0842	µg/l	1/21/2020 12:03:13	
50PPBCN	A	Reactive Cl P	50.5439	µg/l	1/21/2020 12:03:14	
100PPBCN	A	Reactive Cl P	100.0913	µg/l	1/21/2020 12:03:15	
250PPBCN	A	Reactive Cl P	249.4614	µg/l	1/21/2020 12:03:16	
500PPBCN	A	Reactive Cl P	500.1967	µg/l	1/21/2020 12:03:17	
ICV1	S	Reactive Cl P	95.1161	µg/l	1/21/2020 13:53:26	
ICB1	S	Reactive Cl P	0.3216	µg/l	1/21/2020 13:53:27	
CCV1	S	Reactive Cl P	240.9792	µg/l	1/21/2020 13:53:28	
CCB1	S	Reactive Cl P	0.2955	µg/l	1/21/2020 13:53:29	
PB126248BL	S	Reactive Cl P	0.201	µg/l	1/21/2020 13:53:30	
L1004-15	S	Reactive Cl P	0.1209	µg/l	1/21/2020 13:53:31	
L1004-15DUP	S	Reactive Cl P	0.2467	µg/l	1/21/2020 13:53:32	
L1210-07	S	Reactive Cl P	0.1682	µg/l	1/21/2020 13:53:33	
L1211-02	S	Reactive Cl P	0.1352	µg/l	1/21/2020 13:53:34	
L1211-04	S	Reactive Cl P	0.0278	µg/l	1/21/2020 13:53:35	
CCV2	S	Reactive Cl P	238.7025	µg/l	1/21/2020 13:53:36	
CCB2	S	Reactive Cl P	0.0593	µg/l	1/21/2020 13:56:14	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092

Reviewed by : AK Instrument ID : Konelab

1/21/2020 12:29


Test Total CN

Accepted 1/21/2020 12:29

Factor 1018
Bias 0.002

Coeff. of det. 0.999997

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.002	-0.1993	0.0000	
2	5.0PPBCN	0.007	4.8217	5.0000	
3	10PPBCN	0.012	10.0842	10.0000	
4	50PPBCN	0.052	50.5439	50.0000	
5	100PPBCN	0.101	100.0913	100.0000	
6	250PPBCN	0.247	249.4614	250.0000	
7	500PPBCN	0.493	500.1967	500.0000	