

LB107537

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : AK

Instrument ID : Konelab

2/4/2020 15:19

Test: ~~Total~~ CN

Reactive
WJL

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	96.757	0.0	0.106	
ICB1	-1.374	0.0	0.003	
CCV1	250.183	0.0	0.267	
CCB1	-1.626	0.0	0.003	
PB126544BL	-1.741	0.0	0.003	
L1016-19	-1.807	0.0	0.003	
L1016-19DUP	-1.906	0.0	0.002	
L1343-02	-1.485	0.0	0.003	
L1363-05	-1.931	0.0	0.002	
L1368-06	-1.736	0.0	0.003	
L1383-08	-0.785	0.0	0.004	
L1383-12	-1.600	0.0	0.003	
PB126590BL	-1.565	0.0	0.003	
CCV2	251.437	0.0	0.268	
CCB2	-1.819	0.0	0.003	
L1377-01	-1.816	0.0	0.003	
L1377-01DUP	-1.737	0.0	0.003	
L1387-02	-1.728	0.0	0.003	
L1388-02	-1.618	0.0	0.003	
L1383-04	-1.601	0.0	0.003	
CCV3	268.868	0.0	0.286	
CCB3	-1.835	0.0	0.002	

N	22
Mean	38.070
SD	91.4511
CV%	240.22

ALC

Aquakem v. 7.2AQ1

Results from time period:

Tue Feb 04 14:50:59 2020

Tue Feb 04 15:18:12 2020

Sample Id	Sam/Ctr/c/	Test short	Test type	Result	Result unit	Result date and time	Stat
0.0PPBCN	A	Reactive	Cl P	-1.8799	µg/l	2/4/2020 9:39:00	
5.0PPBCN	A	Reactive	Cl P	3.3933	µg/l	2/4/2020 9:39:01	
10PPBCN	A	Reactive	Cl P	8.399	µg/l	2/4/2020 9:39:02	
50PPBCN	A	Reactive	Cl P	46.9428	µg/l	2/4/2020 9:39:03	
100PPBCN	A	Reactive	Cl P	102.0544	µg/l	2/4/2020 9:39:04	
250PPBCN	A	Reactive	Cl P	262.295	µg/l	2/4/2020 9:39:05	
500PPBCN	A	Reactive	Cl P	493.7954	µg/l	2/4/2020 9:39:06	
ICV1	S	Reactive	Cl P	96.7568	µg/l	2/4/2020 14:50:59	
ICB1	S	Reactive	Cl P	-1.3738	µg/l	2/4/2020 14:51:00	
CCV1	S	Reactive	Cl P	250.1833	µg/l	2/4/2020 14:51:01	
CCB1	S	Reactive	Cl P	-1.6264	µg/l	2/4/2020 14:51:02	
PB126544BL	S	Reactive	Cl P	-1.7411	µg/l	2/4/2020 14:51:03	
L1016-19	S	Reactive	Cl P	-1.8068	µg/l	2/4/2020 14:51:04	
L1016-19DUP	S	Reactive	Cl P	-1.9062	µg/l	2/4/2020 14:51:05	
L1343-02	S	Reactive	Cl P	-1.4846	µg/l	2/4/2020 14:51:06	
L1363-05	S	Reactive	Cl P	-1.9311	µg/l	2/4/2020 14:51:07	
L1368-06	S	Reactive	Cl P	-1.7361	µg/l	2/4/2020 14:51:08	
L1383-08	S	Reactive	Cl P	-0.7854	µg/l	2/4/2020 14:58:34	
L1383-12	S	Reactive	Cl P	-1.5997	µg/l	2/4/2020 14:58:35	
PB126590BL	S	Reactive	Cl P	-1.5651	µg/l	2/4/2020 14:58:36	
CCV2	S	Reactive	Cl P	251.4369	µg/l	2/4/2020 14:58:37	
CCB2	S	Reactive	Cl P	-1.8195	µg/l	2/4/2020 14:58:38	
L1377-01	S	Reactive	Cl P	-1.8157	µg/l	2/4/2020 14:58:39	
L1377-01DUP	S	Reactive	Cl P	-1.7369	µg/l	2/4/2020 14:58:40	
L1387-02	S	Reactive	Cl P	-1.7284	µg/l	2/4/2020 14:58:41	
L1388-02	S	Reactive	Cl P	-1.6183	µg/l	2/4/2020 14:58:42	
L1383-04	S	Reactive	Cl P	-1.6007	µg/l	2/4/2020 15:18:10	
CCV3	S	Reactive	Cl P	268.8677	µg/l	2/4/2020 15:18:11	
CCB3	S	Reactive	Cl P	-1.8349	µg/l	2/4/2020 15:18:12	

CHEMTECH
 284 Sheffield Street,
 Mountainside, NJ 07092
 Reviewed by : AL

Instrument ID : Konelab

2/4/2020 15:53


Test Total CN

Accepted 2/4/2020 9:46

Factor 953.8
 Bias 0.004

Coeff. of det. 0.998970

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors

1	0.0PPBCN	0.002	-1.8799	0.0000	
2	5.0PPBCN	0.008	3.3933	5.0000	
3	10PPBCN	0.013	8.3990	10.0000	
4	50PPBCN	0.054	46.9428	50.0000	
5	100PPBCN	0.111	102.0544	100.0000	
6	250PPBCN	0.279	262.2950	250.0000	
7	500PPBCN	0.522	493.7954	500.0000	