

CB108067

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : 

Instrument ID : Konelab

3/3/2020 13:47

Test: Total CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	93.995	0.0	0.078	
ICB1	0.832	0.0	0.003	
CCV1	241.857	0.0	0.197	
CCB1	0.652	0.0	0.003	
PB127242BL	0.344	0.0	0.003	
L1723-03	0.407	0.0	0.003	
L1723-03DUP	0.429	0.0	0.003	
L1743-19	0.324	0.0	0.003	
L1743-20	0.299	0.0	0.003	
L1743-21	0.471	0.0	0.003	
L1743-22	0.728	0.0	0.003	
L1743-24	0.996	0.0	0.003	
L1743-23	0.951	0.0	0.003	
L1743-25	1.154	0.0	0.003	
CCV2	244.156	0.0	0.199	
CCB2	0.009	0.0	0.002	
L1743-26	0.592	0.0	0.003	
L1743-27	0.489	0.0	0.003	
L1757-06	0.749	0.0	0.003	
L1762-03	0.579	0.0	0.003	
L1762-06	0.608	0.0	0.003	
L1762-09	0.498	0.0	0.003	
CCV3	243.278	0.0	0.198	
CCB3	0.461	0.0	0.003	

N	24
Mean	34.786
SD	82.6447
CV%	237.58

Aquakem v. 7.2AQ1

Results from time period:

Tue Mar 03 13:06:10 2020

Tue Mar 03 13:45:08 2020

AL

Sample Id	Sam/Ctr/c	Test short	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A		Reactive Cl P	0.893	µg/l	3/3/2020 11:42:16	
5.OPPBCN	A		Reactive Cl P	5.4813	µg/l	3/3/2020 11:42:17	
10PPBCN	A		Reactive Cl P	10.5926	µg/l	3/3/2020 11:42:18	
50PPBCN	A		Reactive Cl P	48.9544	µg/l	3/3/2020 11:42:19	
100PPBCN	A		Reactive Cl P	97.5863	µg/l	3/3/2020 11:42:20	
250PPBCN	A		Reactive Cl P	251.8436	µg/l	3/3/2020 11:42:21	
500PPBCN	A		Reactive Cl P	499.6489	µg/l	3/3/2020 11:42:22	
ICV1	S		Reactive Cl P	93.9954	µg/l	3/3/2020 13:06:10	
ICB1	S		Reactive Cl P	0.8321	µg/l	3/3/2020 13:06:11	
CCV1	S		Reactive Cl P	241.8568	µg/l	3/3/2020 13:06:12	
CCB1	S		Reactive Cl P	0.6517	µg/l	3/3/2020 13:06:13	
PB127242BL	S		Reactive Cl P	0.3444	µg/l	3/3/2020 13:06:14	
L1723-03	S		Reactive Cl P	0.4071	µg/l	3/3/2020 13:06:15	
L1723-03DUP	S		Reactive Cl P	0.4285	µg/l	3/3/2020 13:06:16	
L1743-19	S		Reactive Cl P	0.3241	µg/l	3/3/2020 13:06:17	
L1743-20	S		Reactive Cl P	0.2986	µg/l	3/3/2020 13:06:18	
L1743-21	S		Reactive Cl P	0.4715	µg/l	3/3/2020 13:06:19	
L1743-22	S		Reactive Cl P	0.7279	µg/l	3/3/2020 13:06:20	
L1743-24	S		Reactive Cl P	0.9958	µg/l	3/3/2020 13:13:46	
L1743-23	S		Reactive Cl P	0.9512	µg/l	3/3/2020 13:42:20	
L1743-25	S		Reactive Cl P	1.1542	µg/l	3/3/2020 13:42:21	
CCV2	S		Reactive Cl P	244.1561	µg/l	3/3/2020 13:42:22	
CCB2	S		Reactive Cl P	0.0094	µg/l	3/3/2020 13:42:23	
L1743-26	S		Reactive Cl P	0.5921	µg/l	3/3/2020 13:42:24	
L1743-27	S		Reactive Cl P	0.4893	µg/l	3/3/2020 13:42:25	
L1757-06	S		Reactive Cl P	0.7489	µg/l	3/3/2020 13:42:26	
L1762-03	S		Reactive Cl P	0.5785	µg/l	3/3/2020 13:42:27	
L1762-06	S		Reactive Cl P	0.6078	µg/l	3/3/2020 13:42:28	
L1762-09	S		Reactive Cl P	0.498	µg/l	3/3/2020 13:42:29	
CCV3	S		Reactive Cl P	243.2783	µg/l	3/3/2020 13:42:30	
CCB3	S		Reactive Cl P	0.4614	µg/l	3/3/2020 13:45:08	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AM

Instrument ID : Konelab

3/3/2020 11:47


Test Total CN

Accepted 3/3/2020 11:47

Factor 1244
Bias 0.002

Coeff. of det. 0.999942

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.003	0.8930	0.0000	
2	5.0PPBCN	0.007	5.4813	5.0000	
3	10PPBCN	0.011	10.5926	10.0000	
4	50PPBCN	0.042	48.9544	50.0000	
5	100PPBCN	0.081	97.5863	100.0000	
6	250PPBCN	0.205	251.8436	250.0000	
7	500PPBCN	0.404	499.6489	500.0000	