

LB108544

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092

Reviewed by : Dr

Instrument ID : Konelab

4/1/2020 18:20

Test: ^{Plumbe}~~Total~~ CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	97.440	0.0	0.084	
ICB1	1.166	0.0	0.003	
CCV1	235.391	0.0	0.200	
CCB1	1.141	0.0	0.003	
PB127954BL	0.912	0.0	0.003	
L2101-01	1.034	0.0	0.003	
L2101-01DUP	1.103	0.0	0.003	
L2101-02	0.271	0.0	0.002	
L2101-03	1.025	0.0	0.003	
L2103-06	0.769	0.0	0.003	
L2104-01	-0.068	0.0	0.002	
L2104-02	0.960	0.0	0.003	
L2104-03	0.951	0.0	0.003	
L2108-06	0.971	0.0	0.003	
CCV2	234.654	0.0	0.199	
CCB2	1.030	0.0	0.003	
L2109-01	0.944	0.0	0.003	
L2109-02	0.941	0.0	0.003	
PB127955BL	0.539	0.0	0.003	
L2101-04	0.846	0.0	0.003	
L2101-04DUP	0.832	0.0	0.003	
L1762-25	0.978	0.0	0.003	
CCV3	234.532	0.0	0.199	
CCB3	1.316	0.0	0.003	

N	24
Mean	34.153
SD	79.9440
CV%	234.08

Aquakem v. 7.2AQ1

Results from time period:

Wed Apr 01 17:07:02 2020

Wed Apr 01 18:18:09 2020

DK

Sample Id	Sam/Ctr/c/	Test short name	Test type	Result	Result unit	Result date and time
0.0PPBCN	A	REACTIVE CN	P	0.7645	µg/l	4/1/2020 17:30:01
5.0PPBCN	A	REACTIVE CN	P	5.9005	µg/l	4/1/2020 17:30:02
10PPBCN	A	REACTIVE CN	P	11.1097	µg/l	4/1/2020 17:30:03
50PPBCN	A	REACTIVE CN	P	50.3197	µg/l	4/1/2020 17:30:04
100PPBCN	A	REACTIVE CN	P	101.2571	µg/l	4/1/2020 17:30:05
250PPBCN	A	REACTIVE CN	P	241.9265	µg/l	4/1/2020 17:30:06
500PPBCN	A	REACTIVE CN	P	503.7222	µg/l	4/1/2020 17:30:07
ICV1	S	REACTIVE CN	P	97.4397	µg/l	4/1/2020 18:02:53
ICB1	S	REACTIVE CN	P	1.1661	µg/l	4/1/2020 18:02:54
CCV1	S	REACTIVE CN	P	235.3906	µg/l	4/1/2020 18:02:57
CCB1	S	REACTIVE CN	P	1.1407	µg/l	4/1/2020 18:02:58
PB127954BL	S	REACTIVE CN	P	0.9123	µg/l	4/1/2020 18:03:00
L2101-01	S	REACTIVE CN	P	1.0339	µg/l	4/1/2020 18:03:02
L2101-01DUP	S	REACTIVE CN	P	1.1029	µg/l	4/1/2020 18:10:24
L2101-02	S	REACTIVE CN	P	0.2705	µg/l	4/1/2020 18:10:25
L2101-03	S	REACTIVE CN	P	1.0245	µg/l	4/1/2020 18:10:26
L2103-06	S	REACTIVE CN	P	0.7694	µg/l	4/1/2020 18:10:27
L2104-01	S	REACTIVE CN	P	-0.0685	µg/l	4/1/2020 18:10:28
L2104-02	S	REACTIVE CN	P	0.96	µg/l	4/1/2020 18:10:29
L2104-03	S	REACTIVE CN	P	0.9505	µg/l	4/1/2020 18:10:30
L2108-06	S	REACTIVE CN	P	0.9711	µg/l	4/1/2020 18:10:31
CCV2	S	REACTIVE CN	P	234.654	µg/l	4/1/2020 18:10:33
CCB2	S	REACTIVE CN	P	1.0304	µg/l	4/1/2020 18:17:59
L2109-01	S	REACTIVE CN	P	0.9437	µg/l	4/1/2020 18:18:00
L2109-02	S	REACTIVE CN	P	0.9407	µg/l	4/1/2020 18:18:01
PB127955BL	S	REACTIVE CN	P	0.5387	µg/l	4/1/2020 18:18:02
L2101-04	S	REACTIVE CN	P	0.8457	µg/l	4/1/2020 18:18:03
L2101-04DUP	S	REACTIVE CN	P	0.8318	µg/l	4/1/2020 18:18:04
L1762-25	S	REACTIVE CN	P	0.9785	µg/l	4/1/2020 18:18:05
CCV3	S	REACTIVE CN	P	234.5317	µg/l	4/1/2020 18:18:07
CCB3	S	REACTIVE CN	P	1.3159	µg/l	4/1/2020 18:18:08

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : DK

Instrument ID : Konelab

4/1/2020 17:32


Test ^{Residue} ~~Total~~ CN

Accepted 4/1/2020 17:31

Factor 1193
Bias 0.002

Coeff. of det. 0.999595

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.003	0.7645	0.0000	
2	5.0PPBCN	0.007	5.9005	5.0000	
3	10PPBCN	0.012	11.1097	10.0000	
4	50PPBCN	0.044	50.3197	50.0000	
5	100PPBCN	0.087	101.2571	100.0000	
6	250PPBCN	0.205	241.9265	250.0000	
7	500PPBCN	0.425	503.7222	500.0000	