

LB108767

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : DL

Instrument ID : Konelab

4/17/2020 10:27

Test: ^{peache}~~Total~~ CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	97.106	0.0	0.085	
ICB1	1.130	0.0	0.004	
CCV1	244.695	0.0	0.210	
CCB1	0.816	0.0	0.004	
PB128285BL	0.906	0.0	0.004	
L2279-08	0.796	0.0	0.004	
L2279-08DUP	0.932	0.0	0.004	
L2298-02	0.673	0.0	0.004	
PB128271BL	0.836	0.0	0.004	
L2272-01	0.846	0.0	0.004	
L2272-01DUP	0.909	0.0	0.004	
L2272-02	1.055	0.0	0.004	
L2284-04	0.953	0.0	0.004	
L2284-08	1.082	0.0	0.004	
CCV2	240.743	0.0	0.206	
CCB2	0.640	0.0	0.004	
L2284-12	0.581	0.0	0.003	
L2284-16	0.529	0.0	0.003	
L2284-20	0.410	0.0	0.003	
L2287-05	0.553	0.0	0.003	
L2298-01	0.628	0.0	0.004	
L2300-01	0.564	0.0	0.003	
L2300-02	1.046	0.0	0.004	
L2303-01	0.958	0.0	0.004	
L2303-03	0.916	0.0	0.004	
PB128272BL	0.949	0.0	0.004	
CCV3	243.384	0.0	0.209	
CCB3	0.968	0.0	0.004	
L2304-02	0.766	0.0	0.004	
L2304-02DUP	0.656	0.0	0.004	
L2304-06	0.795	0.0	0.004	
L2304-10	0.530	0.0	0.003	
L2304-14	0.863	0.0	0.004	
L2306-01	1.153	0.0	0.004	
L2307-03	0.497	0.0	0.003	
CCV4	255.658	0.0	0.219	
CCB4	0.880	0.0	0.004	

N	37
Mean	29.930
SD	77.9488
CV%	260.44

Aquakem v. 7.2AQ1

Results from time period:

Fri Apr 17 08:57:55 2020

Fri Apr 17 10:24:40 2020

DL

Sample Id	Sam/Ctr/c/	Test short name	Test type	Result	Result unit	Result date and time
0.OPPBCN	A	REACTIVE CN	P	0.3663	µg/l	4/17/2020 8:57:55
5.OPPBCN	A	REACTIVE CN	P	5.1141	µg/l	4/17/2020 8:57:56
10PPBCN	A	REACTIVE CN	P	9.9569	µg/l	4/17/2020 8:57:57
50PPBCN	A	REACTIVE CN	P	48.5273	µg/l	4/17/2020 8:57:58
100PPBCN	A	REACTIVE CN	P	100.23	µg/l	4/17/2020 8:57:59
250PPBCN	A	REACTIVE CN	P	251.4087	µg/l	4/17/2020 8:58:00
500PPBCN	A	REACTIVE CN	P	499.3966	µg/l	4/17/2020 8:58:01
ICV1	S	REACTIVE CN	P	97.1056	µg/l	4/17/2020 10:05:11
ICB1	S	REACTIVE CN	P	1.1302	µg/l	4/17/2020 10:05:12
CCV1	S	REACTIVE CN	P	244.6952	µg/l	4/17/2020 10:05:13
CCB1	S	REACTIVE CN	P	0.8161	µg/l	4/17/2020 10:05:14
PB128285BL	S	REACTIVE CN	P	0.9057	µg/l	4/17/2020 10:05:15
L2279-08	S	REACTIVE CN	P	0.7964	µg/l	4/17/2020 10:05:16
L2279-08DUP	S	REACTIVE CN	P	0.9315	µg/l	4/17/2020 10:05:17
L2298-02	S	REACTIVE CN	P	0.6728	µg/l	4/17/2020 10:05:18
PB128271BL	S	REACTIVE CN	P	0.8356	µg/l	4/17/2020 10:05:19
L2272-01	S	REACTIVE CN	P	0.8457	µg/l	4/17/2020 10:05:20
L2272-01DUP	S	REACTIVE CN	P	0.9092	µg/l	4/17/2020 10:05:21
L2272-02	S	REACTIVE CN	P	1.0546	µg/l	4/17/2020 10:12:45
L2284-04	S	REACTIVE CN	P	0.9529	µg/l	4/17/2020 10:12:46
L2284-08	S	REACTIVE CN	P	1.0815	µg/l	4/17/2020 10:12:47
CCV2	S	REACTIVE CN	P	240.7428	µg/l	4/17/2020 10:12:48
CCB2	S	REACTIVE CN	P	0.6395	µg/l	4/17/2020 10:12:49
L2284-12	S	REACTIVE CN	P	0.5812	µg/l	4/17/2020 10:12:50
L2284-16	S	REACTIVE CN	P	0.5293	µg/l	4/17/2020 10:12:51
L2284-20	S	REACTIVE CN	P	0.41	µg/l	4/17/2020 10:12:52
L2287-05	S	REACTIVE CN	P	0.553	µg/l	4/17/2020 10:12:53
L2298-01	S	REACTIVE CN	P	0.6281	µg/l	4/17/2020 10:12:54
L2300-01	S	REACTIVE CN	P	0.5638	µg/l	4/17/2020 10:12:55
L2300-02	S	REACTIVE CN	P	1.0457	µg/l	4/17/2020 10:20:18
L2303-01	S	REACTIVE CN	P	0.958	µg/l	4/17/2020 10:20:19
L2303-03	S	REACTIVE CN	P	0.9155	µg/l	4/17/2020 10:20:20
PB128272BL	S	REACTIVE CN	P	0.9486	µg/l	4/17/2020 10:20:21
CCV3	S	REACTIVE CN	P	243.3836	µg/l	4/17/2020 10:20:22
CCB3	S	REACTIVE CN	P	0.9682	µg/l	4/17/2020 10:20:23
L2304-02	S	REACTIVE CN	P	0.7661	µg/l	4/17/2020 10:20:24
L2304-02DUP	S	REACTIVE CN	P	0.6556	µg/l	4/17/2020 10:20:25
L2304-06	S	REACTIVE CN	P	0.795	µg/l	4/17/2020 10:20:26
L2304-10	S	REACTIVE CN	P	0.5301	µg/l	4/17/2020 10:20:27
L2304-14	S	REACTIVE CN	P	0.863	µg/l	4/17/2020 10:20:28
L2306-01	S	REACTIVE CN	P	1.1535	µg/l	4/17/2020 10:24:37
L2307-03	S	REACTIVE CN	P	0.4966	µg/l	4/17/2020 10:24:38

CCV4	S	REACTIVE CN	P	255.6576 µg/l	4/17/2020 10:24:39
CCB4	S	REACTIVE CN	P	0.8795 µg/l	4/17/2020 10:24:40

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092

Reviewed by : JK

Instrument ID : Konelab

4/17/2020 8:58


Test ~~Total~~ ^{Reactive} CN

Accepted 4/17/2020 8:58

Factor 1184
Bias 0.003

Coeff. of det. 0.999977

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.003	0.3663	0.0000	
2	5.0PPBCN	0.007	5.1141	5.0000	
3	10PPBCN	0.011	9.9569	10.0000	
4	50PPBCN	0.044	48.5273	50.0000	
5	100PPBCN	0.088	100.2300	100.0000	
6	250PPBCN	0.215	251.4087	250.0000	
7	500PPBCN	0.425	499.3966	500.0000	