

LB109759

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : Dh

Instrument ID : Konelab

6/26/2020 13:40

Test: ^{Reache}
~~Total~~ CN

My/c

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	93.940	0.0	0.071	
ICB1	1.452	0.0	0.002	
CCV1	243.145	0.0	0.183	
CCB1	1.489	0.0	0.002	
PB129768BL	1.285	0.0	0.002	
L2810-16	1.237	0.0	0.002	
L2810-16DUP	1.599	0.0	0.002	
L2815-15	1.253	0.0	0.002	
L2817-14	1.173	0.0	0.002	
L3111-04	1.307	0.0	0.002	
L3111-09	1.231	0.0	0.002	
L3111-10	0.987	0.0	0.001	
L3113-02	1.346	0.0	0.002	
L3114-02	1.384	0.0	0.002	
CCV2	242.325	0.0	0.183	
CCB2	1.423	0.0	0.002	
L3115-02	1.264	0.0	0.002	
L3126-04	1.185	0.0	0.002	
L3126-08	1.282	0.0	0.002	
L3126-12	1.229	0.0	0.002	
CCV3	243.004	0.0	0.183	
CCB3	1.407	0.0	0.002	

N	22
Mean	38.452
SD	85.4175
CV%	222.14

Aquakem v. 7.2AQ1

Results from time period:

Fri Jun 26 13:32:24 2020

Fri Jun 26 13:40:09 2020

h

Sample Id	Sam/Ctr/c/	Test short name	Test type	Result	Result unit	Result date and time
0.OPPBCN	A	REACTIVE CN	P	0.7006	µg/l	6/26/2020 11:06:19
5.OPPBCN	A	REACTIVE CN	P	5.7347	µg/l	6/26/2020 11:06:20
10PPBCN	A	REACTIVE CN	P	10.8302	µg/l	6/26/2020 11:06:21
50PPBCN	A	REACTIVE CN	P	49.7182	µg/l	6/26/2020 11:06:22
100PPBCN	A	REACTIVE CN	P	97.0969	µg/l	6/26/2020 11:06:23
250PPBCN	A	REACTIVE CN	P	250.6691	µg/l	6/26/2020 11:06:24
500PPBCN	A	REACTIVE CN	P	500.2503	µg/l	6/26/2020 11:06:25
ICV1	S	REACTIVE CN	P	93.9404	µg/l	6/26/2020 13:32:24
ICB1	S	REACTIVE CN	P	1.4522	µg/l	6/26/2020 13:32:25
CCV1	S	REACTIVE CN	P	243.1451	µg/l	6/26/2020 13:32:26
CCB1	S	REACTIVE CN	P	1.4889	µg/l	6/26/2020 13:32:27
PB129768BL	S	REACTIVE CN	P	1.2847	µg/l	6/26/2020 13:32:28
L2810-16	S	REACTIVE CN	P	1.2366	µg/l	6/26/2020 13:32:29
L2810-16DUP	S	REACTIVE CN	P	1.5989	µg/l	6/26/2020 13:32:30
L2815-15	S	REACTIVE CN	P	1.2532	µg/l	6/26/2020 13:32:31
L2817-14	S	REACTIVE CN	P	1.173	µg/l	6/26/2020 13:32:32
L3111-04	S	REACTIVE CN	P	1.3066	µg/l	6/26/2020 13:32:33
L3111-09	S	REACTIVE CN	P	1.2312	µg/l	6/26/2020 13:32:34
L3111-10	S	REACTIVE CN	P	0.9872	µg/l	6/26/2020 13:39:59
L3113-02	S	REACTIVE CN	P	1.3458	µg/l	6/26/2020 13:40:00
L3114-02	S	REACTIVE CN	P	1.3839	µg/l	6/26/2020 13:40:01
CCV2	S	REACTIVE CN	P	242.3246	µg/l	6/26/2020 13:40:02
CCB2	S	REACTIVE CN	P	1.4233	µg/l	6/26/2020 13:40:03
L3115-02	S	REACTIVE CN	P	1.2638	µg/l	6/26/2020 13:40:04
L3126-04	S	REACTIVE CN	P	1.1852	µg/l	6/26/2020 13:40:05
L3126-08	S	REACTIVE CN	P	1.2818	µg/l	6/26/2020 13:40:06
L3126-12	S	REACTIVE CN	P	1.2286	µg/l	6/26/2020 13:40:07
CCV3	S	REACTIVE CN	P	243.0037	µg/l	6/26/2020 13:40:08
CCB3	S	REACTIVE CN	P	1.4066	µg/l	6/26/2020 13:40:09

LB109759

Calibration results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : DK

Instrument ID : Konelab

6/26/2020 13:15


Test ~~Total~~ ^{Reuche} CN

Accepted 6/26/2020 11:06

Factor 1331
Bias 0.001

Coeff. of det. 0.999948

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.001	0.7006	0.0000	
2	5.0PPBCN	0.005	5.7347	5.0000	
3	10PPBCN	0.009	10.8302	10.0000	
4	50PPBCN	0.038	49.7182	50.0000	
5	100PPBCN	0.074	97.0969	100.0000	
6	250PPBCN	0.189	250.6691	250.0000	
7	500PPBCN	0.377	500.2503	500.0000	