

LB109986

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : *DM*

Instrument ID : Konelab

7/13/2020 14:16

Test: ^{*Reche*} ~~Total~~ CN *My/L*

Sample Id	Result	Dil. 1 +	Response	Errors

ICV1	95.384	0.0	0.078	
ICB1	0.718	0.0	0.003	
CCV1	242.135	0.0	0.195	
CCB1	0.795	0.0	0.003	
PB130104BL	0.718	0.0	0.003	
L3181-05	0.572	0.0	0.002	
L3184-06	0.629	0.0	0.002	
L3184-06DUP	0.629	0.0	0.002	
L3256-04	0.769	0.0	0.003	
L3259-02	0.512	0.0	0.002	
CCV2	243.409	0.0	0.196	
CCB2	0.282	0.0	0.002	

N	12
Mean	48.879
SD	94.5366
CV%	193.41

Aquakem v. 7.2AQ1

Results from time period:

Mon Jul 13 14:12:39 2020

Mon Jul 13 14:15:27 2020

OK

Sample Id	Sam/Ctr/c/	Test short name	Test type	Result	Result unit	Result date and time
0.OPPBCN	A	REACTIVE CN	P	-0.2766	µg/l	7/13/2020 12:29:15
5.OPPBCN	A	REACTIVE CN	P	4.8843	µg/l	7/13/2020 12:29:16
10PPBCN	A	REACTIVE CN	P	9.9455	µg/l	7/13/2020 12:29:17
50PPBCN	A	REACTIVE CN	P	49.3385	µg/l	7/13/2020 12:29:18
100PPBCN	A	REACTIVE CN	P	99.5998	µg/l	7/13/2020 12:29:19
250PPBCN	A	REACTIVE CN	P	252.7199	µg/l	7/13/2020 12:29:20
500PPBCN	A	REACTIVE CN	P	498.7885	µg/l	7/13/2020 12:29:21
ICV1	S	REACTIVE CN	P	95.3844	µg/l	7/13/2020 14:12:39
ICB1	S	REACTIVE CN	P	0.7183	µg/l	7/13/2020 14:12:40
CCV1	S	REACTIVE CN	P	242.1347	µg/l	7/13/2020 14:12:41
CCB1	S	REACTIVE CN	P	0.7953	µg/l	7/13/2020 14:12:42
PB130104BL	S	REACTIVE CN	P	0.7179	µg/l	7/13/2020 14:12:43
L3181-05	S	REACTIVE CN	P	0.5716	µg/l	7/13/2020 14:12:44
L3184-06	S	REACTIVE CN	P	0.6289	µg/l	7/13/2020 14:12:45
L3184-06DUP	S	REACTIVE CN	P	0.6293	µg/l	7/13/2020 14:12:46
L3256-04	S	REACTIVE CN	P	0.7689	µg/l	7/13/2020 14:12:47
L3259-02	S	REACTIVE CN	P	0.5122	µg/l	7/13/2020 14:12:48
CCV2	S	REACTIVE CN	P	243.4095	µg/l	7/13/2020 14:12:49
CCB2	S	REACTIVE CN	P	0.2821	µg/l	7/13/2020 14:15:27

LB69486

Calibration results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
 284 Sheffield Street,
 Mountainside, NJ 07092
 Reviewed by : Dh

Instrument ID : Konelab

7/13/2020 14:00


Test ^{Revised} ~~Total~~ CN

Accepted 7/13/2020 12:29

Factor 1254
 Bias 0.002

Coeff. of det. 0.999954

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.002	-0.2766	0.0000	
2	5.0PPBCN	0.006	4.8843	5.0000	
3	10PPBCN	0.010	9.9455	10.0000	
4	50PPBCN	0.041	49.3385	50.0000	
5	100PPBCN	0.081	99.5998	100.0000	
6	250PPBCN	0.203	252.7199	250.0000	
7	500PPBCN	0.400	498.7885	500.0000	