

LB110086

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : Dr

Instrument ID : Konelab

7/17/2020 14:17

Test: Total CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	97.971	0.0	0.079	
ICB1	-0.482	0.0	0.001	
CCV1	243.555	0.0	0.195	
CCB1	-0.200	0.0	0.002	
PB130283BL	-0.612	0.0	0.001	
L3188-16	-0.658	0.0	0.001	
L3188-16DUP	-0.288	0.0	0.001	
L3335-04	-0.499	0.0	0.001	
L3335-08	-0.328	0.0	0.001	
L3335-12	-0.354	0.0	0.001	
CCV2	246.164	0.0	0.197	
CCB2	-0.743	0.0	0.001	

N	12
Mean	48.627
SD	95.8886
CV%	197.19

Aquakem v. 7.2AQ1

Results from time period:

Fri Jul 17 14:14:12 2020

Fri Jul 17 14:17:00 2020

De

Sample Id	Sam/Ctr/c	Test short name	Test type	Result	Result unit	Result date and time
0.OPPBCN	A	REACTIVE CN	P	-0.8909	µg/l	7/17/2020 11:40:53
5.OPPBCN	A	REACTIVE CN	P	4.1543	µg/l	7/17/2020 11:40:54
10PPBCN	A	REACTIVE CN	P	11.309	µg/l	7/17/2020 11:40:55
50PPBCN	A	REACTIVE CN	P	47.4186	µg/l	7/17/2020 11:40:56
100PPBCN	A	REACTIVE CN	P	100.317	µg/l	7/17/2020 11:40:57
250PPBCN	A	REACTIVE CN	P	255.0299	µg/l	7/17/2020 11:40:58
500PPBCN	A	REACTIVE CN	P	497.6621	µg/l	7/17/2020 11:40:59
ICV1	S	REACTIVE CN	P	97.9708	µg/l	7/17/2020 14:14:12
ICB1	S	REACTIVE CN	P	-0.4821	µg/l	7/17/2020 14:14:13
CCV1	S	REACTIVE CN	P	243.5554	µg/l	7/17/2020 14:14:14
CCB1	S	REACTIVE CN	P	-0.1997	µg/l	7/17/2020 14:14:15
PB130283BL	S	REACTIVE CN	P	-0.6118	µg/l	7/17/2020 14:14:16
L3188-16	S	REACTIVE CN	P	-0.6576	µg/l	7/17/2020 14:14:17
L3188-16DUP	S	REACTIVE CN	P	-0.2875	µg/l	7/17/2020 14:14:18
L3335-04	S	REACTIVE CN	P	-0.4989	µg/l	7/17/2020 14:14:19
L3335-08	S	REACTIVE CN	P	-0.3277	µg/l	7/17/2020 14:14:20
L3335-12	S	REACTIVE CN	P	-0.3538	µg/l	7/17/2020 14:14:21
CCV2	S	REACTIVE CN	P	246.1638	µg/l	7/17/2020 14:14:22
CCB2	S	REACTIVE CN	P	-0.7431	µg/l	7/17/2020 14:17:00

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : DN

Instrument ID : Konelab

7/17/2020 15:02


Test Total CN

Accepted 7/17/2020 11:46

Factor 1262
Bias 0.002

Coeff. of det. 0.999802

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.001	-0.8909	0.0000	
2	5.0PPBCN	0.005	4.1543	5.0000	
3	10PPBCN	0.011	11.3090	10.0000	
4	50PPBCN	0.039	47.4186	50.0000	
5	100PPBCN	0.081	100.3170	100.0000	
6	250PPBCN	0.204	255.0299	250.0000	
7	500PPBCN	0.396	497.6621	500.0000	