

LB110683

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : AM

Instrument ID : Konelab

8/25/2020 11:57

Test: Total CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	98.513	0.0	0.073	
ICB1	0.263	0.0	0.001	
CCV1	242.931	0.0	0.177	
CCB1	0.238	0.0	0.001	
PB131162BL	0.122	0.0	0.001	
L3773-02	0.071	0.0	0.001	
L3773-02DUP	0.041	0.0	0.001	
L3773-04	0.034	0.0	0.001	
L3780-08	0.255	0.0	0.001	
L3780-16	-0.132	0.0	0.001	
L3780-23	0.099	0.0	0.001	
L3784-02	0.138	0.0	0.001	
L3784-05	0.027	0.0	0.001	
L3784-08	0.836	0.0	0.002	
CCV2	242.650	0.0	0.177	
CCB2	0.350	0.0	0.001	
L3784-11	0.248	0.0	0.001	
L3784-14	-0.133	0.0	0.001	
L3784-17	0.250	0.0	0.001	
L3785-02	-0.410	0.0	0.001	
L3785-05	0.095	0.0	0.001	
CCV3	240.501	0.0	0.176	
CCB3	-0.320	0.0	0.001	

N	23
Mean	35.942
SD	84.1345
CV%	234.08

Aquakem v. 7.2AQ1

Results from time period:

Tue Aug 25 09:47:22 2020

Tue Aug 25 11:54:05 2020

Sample Id	Sam/Ctr/c/	Test short name	Test type	Result	Result unit	Result date and time	Stat
0.0PPBCN	A	Reactive CN	P	-0.5408	µg/l	8/25/2020 10:39:10	
5.0PPBCN	A	Reactive CN	P	4.9101	µg/l	8/25/2020 10:39:11	
10PPBCN	A	Reactive CN	P	10.0209	µg/l	8/25/2020 10:39:12	
50PPBCN	A	Reactive CN	P	50.7358	µg/l	8/25/2020 10:39:13	
100PPBCN	A	Reactive CN	P	101.512	µg/l	8/25/2020 10:39:14	
250PPBCN	A	Reactive CN	P	247.475	µg/l	8/25/2020 10:39:15	
500PPBCN	A	Reactive CN	P	500.887	µg/l	8/25/2020 10:39:16	
ICV1	S	Reactive CN	P	98.5128	µg/l	8/25/2020 11:43:42	
ICB1	S	Reactive CN	P	0.2634	µg/l	8/25/2020 11:43:43	
CCV1	S	Reactive CN	P	242.9307	µg/l	8/25/2020 11:43:44	
CCB1	S	Reactive CN	P	0.2384	µg/l	8/25/2020 11:43:45	
PB131162BL	S	Reactive CN	P	0.1225	µg/l	8/25/2020 11:43:46	
L3773-02	S	Reactive CN	P	0.0711	µg/l	8/25/2020 11:43:47	
L3773-02DUP	S	Reactive CN	P	0.0409	µg/l	8/25/2020 11:43:48	
L3773-04	S	Reactive CN	P	0.0343	µg/l	8/25/2020 11:43:49	
L3780-08	S	Reactive CN	P	0.2551	µg/l	8/25/2020 11:43:50	
L3780-16	S	Reactive CN	P	-0.132	µg/l	8/25/2020 11:43:51	
L3780-23	S	Reactive CN	P	0.0988	µg/l	8/25/2020 11:43:52	
L3784-02	S	Reactive CN	P	0.1384	µg/l	8/25/2020 11:51:17	
L3784-05	S	Reactive CN	P	0.0273	µg/l	8/25/2020 11:51:18	
L3784-08	S	Reactive CN	P	0.836	µg/l	8/25/2020 11:51:19	
CCV2	S	Reactive CN	P	242.6505	µg/l	8/25/2020 11:51:20	
CCB2	S	Reactive CN	P	0.3495	µg/l	8/25/2020 11:51:21	
L3784-11	S	Reactive CN	P	0.2481	µg/l	8/25/2020 11:51:22	
L3784-14	S	Reactive CN	P	-0.1335	µg/l	8/25/2020 11:51:23	
L3784-17	S	Reactive CN	P	0.2503	µg/l	8/25/2020 11:51:24	
L3785-02	S	Reactive CN	P	-0.4096	µg/l	8/25/2020 11:51:25	
L3785-05	S	Reactive CN	P	0.0953	µg/l	8/25/2020 11:51:26	
CCV3	S	Reactive CN	P	240.501	µg/l	8/25/2020 11:51:27	
CCB3	S	Reactive CN	P	-0.32	µg/l	8/25/2020 11:54:05	

CHEMTECH
 284 Sheffield Street,
 Mountainside, NJ 07092
 Reviewed by : AM

Instrument ID : Konelab

8/25/2020 10:43


Test Total CN

Accepted 8/25/2020 10:43

Factor 1376
 Bias 0.001

Coeff. of det. 0.999950

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors

1	0.0PPBCN	0.001	-0.5408	0.0000	
2	5.0PPBCN	0.005	4.9101	5.0000	
3	10PPBCN	0.008	10.0209	10.0000	
4	50PPBCN	0.038	50.7358	50.0000	
5	100PPBCN	0.075	101.5120	100.0000	
6	250PPBCN	0.181	247.4750	250.0000	
7	500PPBCN	0.365	500.8870	500.0000	