

LB110878

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : AM

Instrument ID : Konelab

9/9/2020 12:36

Test: Total CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	95.110	0.0	0.080	
ICB1	0.633	0.0	0.000	
CCV1	240.396	0.0	0.203	
CCB1	0.416	0.0	0.000	
PB131475BL	0.509	0.0	0.000	
L3808-11	0.544	0.0	0.000	
L3808-11DUP	0.073	0.0	0.000	
L3808-13	0.178	0.0	0.000	
L3808-14	0.471	0.0	0.000	
L3870-05	0.293	0.0	0.000	
L3872-04	0.317	0.0	0.000	
L3915-02	0.442	0.0	0.000	
L3915-06	0.457	0.0	0.000	
L3915-10	0.488	0.0	0.000	
CCV2	238.737	0.0	0.202	
CCB2	0.404	0.0	0.000	
L3915-14	0.424	0.0	0.000	
L3931-02	0.423	0.0	0.000	
L3931-04	0.202	0.0	0.000	
L3932-02	0.292	0.0	0.000	
L3932-04	0.251	0.0	0.000	
CCV3	235.694	0.0	0.199	
CCB3	0.142	0.0	0.000	

N	23
Mean	35.517
SD	82.6747
CV%	232.77

L3939
AM
9/09/2020

Aquakem v. 7.2AQ1

Results from time period:

Wed Sep 09 11:44:19 2020

Wed Sep 09 12:34:44 2020

Sample Id	Sam/Ctr/c/	Test short nan	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive CN	P	0.0196	µg/l	9/9/2020 11:44:19	
5.OPPBCN	A	Reactive CN	P	5.086	µg/l	9/9/2020 11:44:20	
10PPBCN	A	Reactive CN	P	10.2138	µg/l	9/9/2020 11:44:21	
50PPBCN	A	Reactive CN	P	50.7701	µg/l	9/9/2020 11:44:22	
100PPBCN	A	Reactive CN	P	99.5094	µg/l	9/9/2020 11:44:23	
250PPBCN	A	Reactive CN	P	248.7704	µg/l	9/9/2020 11:44:24	
500PPBCN	A	Reactive CN	P	500.6308	µg/l	9/9/2020 11:44:25	
ICV1	S	Reactive CN	P	95.1098	µg/l	9/9/2020 12:19:34	
ICB1	S	Reactive CN	P	0.6327	µg/l	9/9/2020 12:19:35	
CCV1	S	Reactive CN	P	240.3961	µg/l	9/9/2020 12:19:36	
CCB1	S	Reactive CN	P	0.4159	µg/l	9/9/2020 12:19:37	
PB131475BL	S	Reactive CN	P	0.5095	µg/l	9/9/2020 12:19:38	
L3808-11	S	Reactive CN	P	0.5438	µg/l	9/9/2020 12:19:39	
L3808-11DUP	S	Reactive CN	P	0.0727	µg/l	9/9/2020 12:19:40	
L3808-13	S	Reactive CN	P	0.1784	µg/l	9/9/2020 12:19:41	
L3808-14	S	Reactive CN	P	0.4706	µg/l	9/9/2020 12:19:42	
L3870-05	S	Reactive CN	P	0.2933	µg/l	9/9/2020 12:19:43	
L3872-04	S	Reactive CN	P	0.3167	µg/l	9/9/2020 12:19:44	
L3915-02	S	Reactive CN	P	0.4421	µg/l	9/9/2020 12:27:09	
L3915-06	S	Reactive CN	P	0.4566	µg/l	9/9/2020 12:27:10	
L3915-10	S	Reactive CN	P	0.4881	µg/l	9/9/2020 12:27:11	
CCV2	S	Reactive CN	P	238.7373	µg/l	9/9/2020 12:27:12	
CCB2	S	Reactive CN	P	0.4041	µg/l	9/9/2020 12:27:13	
L3915-14	S	Reactive CN	P	0.4243	µg/l	9/9/2020 12:27:14	
L3931-02	S	Reactive CN	P	0.4226	µg/l	9/9/2020 12:27:15	
L3931-04	S	Reactive CN	P	0.2025	µg/l	9/9/2020 12:27:16	
L3932-02	S	Reactive CN	P	0.2924	µg/l	9/9/2020 12:27:17	
L3939-04	S	Reactive CN	P	0.2514	µg/l	9/9/2020 12:27:18	
CCV3	S	Reactive CN	P	235.694	µg/l	9/9/2020 12:27:19	
CCB3	S	Reactive CN	P	0.1416	µg/l	9/9/2020 12:34:44	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AM

Instrument ID : Konelab

9/9/2020 11:47


Test Total CN

Accepted 9/9/2020 11:47

Factor 1182
Bias 0

Coeff. of det. 0.999986

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	-0.000	0.0196	0.0000	
2	5.0PPBCN	0.004	5.0860	5.0000	
3	10PPBCN	0.009	10.2138	10.0000	
4	50PPBCN	0.043	50.7701	50.0000	
5	100PPBCN	0.084	99.5094	100.0000	
6	250PPBCN	0.210	248.7704	250.0000	
7	500PPBCN	0.423	500.6308	500.0000	