

LB111519

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AM

Instrument ID : Konelab

10/15/2020 12:45

Test: Total CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	97.110	0.0	0.085	
ICB1	0.066	0.0	0.000	
CCV1	246.661	0.0	0.214	
CCB1	0.970	0.0	0.001	
PB132362BL	-0.016	0.0	0.000	
L4223-08	0.139	0.0	0.000	
L4226-10	-0.094	0.0	0.000	
L4226-10DUP	-0.071	0.0	0.000	
L4226-16	-0.069	0.0	0.000	
L4227-04	-0.027	0.0	0.000	
L4393-01	0.046	0.0	0.000	
L4393-02	-0.099	0.0	0.000	
L4394-01	-0.260	0.0	0.000	
L4395-01	-0.148	0.0	0.000	
CCV2	245.021	0.0	0.213	
CCB2	0.151	0.0	0.000	
L4395-03	0.226	0.0	0.000	
L4404-02	0.137	0.0	0.000	
L4406-08	0.166	0.0	0.000	
L4406-13	0.173	0.0	0.000	
CCV3	248.920	0.0	0.216	
CCB3	0.752	0.0	0.001	

N	22
Mean	38.171
SD	87.3461
CV%	228.83

Aquakem v. 7.2AQ1

Results from time period:

Thu Oct 15 11:31:10 2020

Thu Oct 15 12:40:14 2020

Sample Id	Sam/Ctr/c/	Test short name	Test type	Result	Result unit	Result date and time
0.0PPBCN	A	Reactive CN	P	0.1385	µg/l	10/15/2020 11:31:10
5.0PPBCN	A	Reactive CN	P	4.8398	µg/l	10/15/2020 11:31:11
10PPBCN	A	Reactive CN	P	9.7556	µg/l	10/15/2020 11:31:12
50PPBCN	A	Reactive CN	P	48.7658	µg/l	10/15/2020 11:31:13
100PPBCN	A	Reactive CN	P	99.8547	µg/l	10/15/2020 11:31:14
250PPBCN	A	Reactive CN	P	252.9732	µg/l	10/15/2020 11:31:15
500PPBCN	A	Reactive CN	P	498.6724	µg/l	10/15/2020 11:31:16
ICV1	S	Reactive CN	P	97.1102	µg/l	10/15/2020 12:10:23
ICB1	S	Reactive CN	P	0.0658	µg/l	10/15/2020 12:10:24
CCV1	S	Reactive CN	P	246.6609	µg/l	10/15/2020 12:10:27
CCB1	S	Reactive CN	P	0.9701	µg/l	10/15/2020 12:10:28
PB132362BL	S	Reactive CN	P	-0.0163	µg/l	10/15/2020 12:10:30
L4223-08	S	Reactive CN	P	0.1394	µg/l	10/15/2020 12:10:31
L4226-10	S	Reactive CN	P	-0.0943	µg/l	10/15/2020 12:16:33
L4226-10DUP	S	Reactive CN	P	-0.0714	µg/l	10/15/2020 12:16:34
L4226-16	S	Reactive CN	P	-0.0685	µg/l	10/15/2020 12:16:35
L4227-04	S	Reactive CN	P	-0.0271	µg/l	10/15/2020 12:16:36
L4393-01	S	Reactive CN	P	0.0461	µg/l	10/15/2020 12:16:37
L4393-02	S	Reactive CN	P	-0.0995	µg/l	10/15/2020 12:16:38
L4394-01	S	Reactive CN	P	-0.26	µg/l	10/15/2020 12:16:39
L4395-01	S	Reactive CN	P	-0.1484	µg/l	10/15/2020 12:16:40
CCV2	S	Reactive CN	P	245.021	µg/l	10/15/2020 12:37:27
CCB2	S	Reactive CN	P	0.1513	µg/l	10/15/2020 12:37:28
L4395-03	S	Reactive CN	P	0.2261	µg/l	10/15/2020 12:37:30
L4404-02	S	Reactive CN	P	0.1369	µg/l	10/15/2020 12:37:31
L4406-08	S	Reactive CN	P	0.1663	µg/l	10/15/2020 12:37:32
L4406-13	S	Reactive CN	P	0.173	µg/l	10/15/2020 12:37:33
CCV3	S	Reactive CN	P	248.9198	µg/l	10/15/2020 12:37:35
CCB3	S	Reactive CN	P	0.7521	µg/l	10/15/2020 12:37:36

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AM

Instrument ID : Konelab

10/15/2020 11:32


Test Total CN

Accepted 10/15/2020 11:31

Factor 1152
Bias 0

Coeff. of det. 0.999940

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.000	0.1385	0.0000	
2	5.0PPBCN	0.004	4.8398	5.0000	
3	10PPBCN	0.009	9.7556	10.0000	
4	50PPBCN	0.043	48.7658	50.0000	
5	100PPBCN	0.087	99.8547	100.0000	
6	250PPBCN	0.220	252.9732	250.0000	
7	500PPBCN	0.433	498.6724	500.0000	