

LB112233

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : AM

Instrument ID : Konelab

12/4/2020 15:18

Test: Total CN

Sample Id	Result ^{ug/l}	Dil. 1 +	Response	Errors
ICV1	94.579	0.0	0.081	
ICB1	1.648	0.0	0.001	
CCV1	241.812	0.0	0.208	
CCB1	2.068	0.0	0.001	
PB133375BL	1.966	0.0	0.001	
L4910-04	1.968	0.0	0.001	
L4910-04DUP	1.247	0.0	0.000	
L4910-06	1.324	0.0	0.000	
L4936-02	1.573	0.0	0.001	
PB133414BL	1.413	0.0	0.000	
L4593-01	1.249	0.0	0.000	
L4593-01DUP	1.140	0.0	0.000	
CCV2	241.479	0.0	0.208	
CCB2	1.789	0.0	0.001	

N	14
Mean	42.518
SD	87.9028
CV%	206.74

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AM

Instrument ID : Konelab

12/4/2020 14:18


Test Total CN

Accepted 12/4/2020 14:18

Factor 1156
Bias -0.001

Coeff. of det. 0.999677

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.000	1.3366	0.0000	
2	5.0PPBCN	0.004	6.0790	5.0000	
3	10PPBCN	0.009	11.2604	10.0000	
4	50PPBCN	0.042	49.4764	50.0000	
5	100PPBCN	0.086	100.4359	100.0000	
6	250PPBCN	0.209	242.9651	250.0000	
7	500PPBCN	0.435	503.4467	500.0000	

Aquakem v. 7.2AQ1

Results from time period:

Fri Dec 04 13:02:53 2020

Fri Dec 04 15:14:59 2020

Sample Id	Sam/Ctr/c/	Test short name	Test type	Result	Result unit	Result date and time
0.OPPBCN	A	Reactive CN	P	1.3366	µg/l	12/4/2020 14:12:46
5.OPPBCN	A	Reactive CN	P	6.079	µg/l	12/4/2020 14:12:47
10PPBCN	A	Reactive CN	P	11.2604	µg/l	12/4/2020 14:12:48
50PPBCN	A	Reactive CN	P	49.4764	µg/l	12/4/2020 14:12:49
100PPBCN	A	Reactive CN	P	100.4359	µg/l	12/4/2020 14:12:50
250PPBCN	A	Reactive CN	P	242.9651	µg/l	12/4/2020 14:12:51
500PPBCN	A	Reactive CN	P	503.4467	µg/l	12/4/2020 14:12:52
ICV1	S	Reactive CN	P	94.5785	µg/l	12/4/2020 15:07:44
ICB1	S	Reactive CN	P	1.6476	µg/l	12/4/2020 15:07:45
CCV1	S	Reactive CN	P	241.8117	µg/l	12/4/2020 15:07:47
CCB1	S	Reactive CN	P	2.0683	µg/l	12/4/2020 15:07:49
PB133375BL	S	Reactive CN	P	1.9661	µg/l	12/4/2020 15:07:51
L4910-04	S	Reactive CN	P	1.9684	µg/l	12/4/2020 15:07:53
L4910-04DUP	S	Reactive CN	P	1.2471	µg/l	12/4/2020 15:14:50
L4910-06	S	Reactive CN	P	1.3236	µg/l	12/4/2020 15:14:51
L4936-02	S	Reactive CN	P	1.5733	µg/l	12/4/2020 15:14:52
PB133414BL	S	Reactive CN	P	1.413	µg/l	12/4/2020 15:14:53
L4593-01	S	Reactive CN	P	1.2494	µg/l	12/4/2020 15:14:54
L4593-01DUP	S	Reactive CN	P	1.1398	µg/l	12/4/2020 15:14:55
CCV2	S	Reactive CN	P	241.4794	µg/l	12/4/2020 15:14:57
CCB2	S	Reactive CN	P	1.7888	µg/l	12/4/2020 15:14:58