

LB88221

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092

Reviewed by : AL

6/20/2017 14:56

Test: ~~Total~~ CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	95.077	0.0	0.090	
ICB1	-0.212	0.0	0.004	
CCV1	237.735	0.0	0.219	
CCB1	-0.011	0.0	0.005	
PB99945BL	-0.223	0.0	0.004	
I3762-03	-0.863	0.0	0.004	
I3762-03DUP	-0.370	0.0	0.004	
I3777-01	-0.504	0.0	0.004	
I3781-07	-0.820	0.0	0.004	
I3781-08	-0.634	0.0	0.004	
CCV2	246.960	0.0	0.227	
CCB2	-0.331	0.0	0.004	

N	12
Mean	47.983
SD	94.8300
CV%	197.63

LB 88221

AK

Aquakem v. 7.2AQ1

Results from time period:

Tue Jun 20 14:37:42 2017

Tue Jun 20 14:40:30 2017

Sample Id	Sam/Ctr/c	Test short	Test type	Result	Result unit	Result date and time	Stat
0.0PPBCN	A	Reactive	CI P	-1.2161	µg/l	6/20/2017 9:39:16	
5.0PPBCN	A	Reactive	CI P	4.1682	µg/l	6/20/2017 9:39:17	
10PPBCN	A	Reactive	CI P	9.1535	µg/l	6/20/2017 9:39:18	
50PPBCN	A	Reactive	CI P	49.6063	µg/l	6/20/2017 9:39:19	
100PPBCN	A	Reactive	CI P	101.7352	µg/l	6/20/2017 9:39:20	
250PPBCN	A	Reactive	CI P	253.6708	µg/l	6/20/2017 9:39:21	
500PPBCN	A	Reactive	CI P	497.8822	µg/l	6/20/2017 9:39:22	
ICV1	S	Reactive	CI P	95.0765	µg/l	6/20/2017 14:37:42	
ICB1	S	Reactive	CI P	-0.212	µg/l	6/20/2017 14:37:43	
CCV1	S	Reactive	CI P	237.7347	µg/l	6/20/2017 14:37:44	
CCB1	S	Reactive	CI P	-0.011	µg/l	6/20/2017 14:37:45	
PB99945BL	S	Reactive	CI P	-0.2231	µg/l	6/20/2017 14:37:46	
I3762-03	S	Reactive	CI P	-0.8634	µg/l	6/20/2017 14:37:47	
I3762-03DUP	S	Reactive	CI P	-0.3704	µg/l	6/20/2017 14:37:48	
I3777-01	S	Reactive	CI P	-0.5041	µg/l	6/20/2017 14:37:49	
I3781-07	S	Reactive	CI P	-0.8204	µg/l	6/20/2017 14:37:50	
I3781-08	S	Reactive	CI P	-0.6341	µg/l	6/20/2017 14:37:51	
CCV2	S	Reactive	CI P	246.96	µg/l	6/20/2017 14:37:52	
CCB2	S	Reactive	CI P	-0.3309	µg/l	6/20/2017 14:40:30	

CHEMTECH

284 Sheffield Street,
Mountainside, NJ 07092Reviewed by : AK

6/20/2017 9:43


Test Total CN

Accepted 6/20/2017 9:43

Factor 1110
Bias 0.005

Coeff. of det. 0.999883

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.003	-1.2161	0.0000	
2	5.0PPBCN	0.008	4.1682	5.0000	
3	10PPBCN	0.013	9.1535	10.0000	
4	50PPBCN	0.049	49.6063	50.0000	
5	100PPBCN	0.096	101.7352	100.0000	
6	250PPBCN	0.233	253.6708	250.0000	
7	500PPBCN	0.453	497.8822	500.0000	