

SOP ID : M1311-TCLP-13
SDG No : N/A
Weigh By : JP
Balance ID : WC SC-4
pH Meter ID : WC PH METER-1
Extraction By : JP
Filter By : JP
Pipette ID : WC
Tumbler ID : T-1 /T-2
TCLP Filter ID : 113075

Start Prep Date : 06/11/2021 **Time :** 15:30
End Prep Date : 06/12/2021 **Time :** 07:25
Combination Ratio : 20
ZHE Cleaning Batch : N/A
Initial Room Temperature: 21 °C
Final Room Temperature: 20 °C
TCLP Technician Signature : *JP*
Supervisor By : *ASB*

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE US	Lot Number
TCLP-FLUID-1	N/A	WP90939
HCL-TCLP,1N	N/A	WP90941
HNO3-TCLP,1N	N/A	WP90942
pH Strips	N/A	W1931,W1934,W2755,W2350
pH Strips	N/A	W1937,W1938,W1939,W1940,W1941,W1942
1L Amber bottle	N/A	19948
120ml Plastic bottle	N/A	210210
N/A	N/A	N/A

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. M2672-21 is used for MS-MSD. Tumbler T-1/T-2 checked, 30 rpm.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
06-12-21 10:30	JP / TCLP Room	RJA / NOA 245 Ext
	Preparation Group	Analysis Group / Met Dig