

SOP ID : M1311-TCLP-13
SDG No : N/A
Weigh By : JP
Balance ID : WC SC-4
pH Meter ID : WC PH METER-1
Extraction By : JP
Filter By : JP
Pipette ID : WC
Tumbler ID : T-1 /T-2
TCLP Filter ID : 113075

Start Prep Date : 09/28/2021 **Time :** 16:35
End Prep Date : 09/29/2021 **Time :** 09:25
Combination Ratio : 20
ZHE Cleaning Batch : N/A
Initial Room Temperature: 24 °C
Final Room Temperature: 22 °C
TCLP Technician Signature : JP
Supervisor By : JS

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE US	Lot Number
TCLP-FLUID-1	N/A	WP92681
HCL-TCLP,1N	N/A	WP92683
HNO3-TCLP,1N	N/A	WP92684
pH Strips	N/A	W1931,W1934,W2755,W2350
pH Strips	N/A	W1937,W1938,W1939,W1940,W1941,W1942
1L Amber bottle	N/A	90521
120ml Plastic bottle	N/A	210210
1:1 HNO3	MP64937	N/A

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. M3971-11 IS USED FOR MS-MSD. Tumbler T-1/T-2 checked,30 rpm.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
09/29/21 10:45	JP /TCLP Room	Re Ext
	Preparation Group	Analysis Group <u>met dig</u>