

Method : H:\VOL1\HPLC\HPLC L\METHODS\HL010313.M (Chemstation Integrator)
 Title : SW8330A/B Primary (Extend C18 3.5u 4.6x10mm)
 Last Update : Fri Jan 04 13:26:48 2013

Calibration Files

CAL1 =HL003992.D CAL2 =HL003991.D CAL3 =HL003990.D
 CAL4 =HL003989.D CAL5 =HL003988.D CAL6 =HL003986.D

	Compound	CAL1	CAL2	CAL3	CAL4	CAL5	CAL6	Avg	%RSD
1) S	4-Nitroaniline	3.423	3.534	3.662	3.623	3.610	3.727	3.596 E3	2.95
2) T	HMX	5.475	4.218	4.017	3.982	3.980	4.111	4.297 E3	13.60
3) T	RDX	4.382	5.490	5.551	5.470	5.412	5.562	5.311 E3	8.63
4) T	1,3,5-Trinitrobenze	1.237	1.197	1.187	1.149	1.136	1.168	1.179 E4	3.10
5) T	Tetryl	8.387	8.517	8.093	8.236	8.107	8.357	8.283 E3	2.02
6) T	2,4,6-Trinitrotolue	1.131	1.150	1.093	1.106	1.088	1.121	1.115 E4	2.15
7) T	2-Amino-4,6-dinitro	1.003	1.057	1.055	1.072	1.053	1.081	1.054 E4	2.55
8) T	4-Amino-2,6-dinitro	7.333	7.557	7.135	7.223	7.057	7.349	7.276 E3	2.45
9) T	2,4-Dinitrotoluene	1.644	1.539	1.519	1.508	1.483	1.527	1.537 E4	3.63
10) T	2,6-Dinitrotoluene	8.710	8.309	8.593	8.410	8.215	8.546	8.464 E3	2.19
11) T	1,3-Dinitrobenzene	1.695	1.701	1.647	1.623	1.598	1.641	1.651 E4	2.45
12) T	Nitrobenzene	1.179	1.115	1.041	1.077	1.034	1.054	1.083 E4	5.11
13) T	2-Nitrotoluene	7.082	6.914	6.568	6.512	6.311	6.504	6.649 E3	4.35
14) T	3-Nitrotoluene	7.201	7.080	7.182	7.129	6.776	7.024	7.065 E3	2.21
15) T	4-Nitrotoluene	5.939	5.681	5.542	5.554	5.323	5.505	5.591 E3	3.69
16) T	3,5-Dinitroaniline	1.722	1.642	1.591	1.565	1.535	1.573	1.605 E4	4.21
17) T	Nitroqlvcerin		6.251	5.679	5.530	5.583	5.624	5.733 E1	5.14
18) T	PETN	9.419	9.017	9.112	8.888	8.809	9.076	9.053 E3	2.35