

Data Path : Z:\HPCHEM1\HPLC_L\Data\HL032013\
 Data File : HL005153.D
 Signal(s) : VWD1A.CH
 Acq On : 20-Mar-2013, 14:00
 Operator : QM
 Sample : HSTDICCC0500
 Misc :
 ALS Vial : 4 Sample Multiplier: 1

Instrument :
 HPLC_L
 LabSampleId :
 HSTDICCC0500

Integration File: autoint1.e
 Quant Time: Mar 22 22:40:34 2013
 Quant Method : Z:\HPCHEM1\HPLC_L\METHODS\HL032013.M
 Quant Title :
 QLast Update : Mon Mar 11 14:17:48 2013
 Response via : Initial Calibration
 Integrator: ChemStation

Volume Inj. : 50 ul
 Signal Phase : ZORBAX Extend-C18
 Signal Info : 4.6 X 100mm 3.5-Micron 80A

Compound	R.T.	Response	Conc Units

System Monitoring Compounds			
1) S 4-Nitroaniline	1.389	2365476	554.972 ug/L
Target Compounds			
2) T HMX	0.936	2388703	517.288 ug/L
3) T RDX	1.813	2815940	523.660 ug/L
4) T 1,3,5-Trinitrobenzene	2.818	6028470	525.139 ug/L
5) T 3,5-Dinitroaniline	3.747	8522411	501.884 ng/ml
6) T 1,3-Dinitrobenzene	4.115	7094739	479.483 ug/L
7) T Nitrobenzene	5.037	5297038	495.094 ug/L
8) T Tetryl	6.558	3967451	553.825 ug/L
9) T 2,4,6-Trinitrotoluene	6.863	5944302	488.783 ug/L
10) T 2-Amino-4,6-dinitroto...	7.385	5885080	485.841 ug/L
11) T 4-Amino-2,6-dinitroto...	7.760	3548238	552.330 ug/L
12) T 2,4-Dinitrotoluene	8.574	8122491	512.439 ug/L
13) T 2,6-Dinitrotoluene	9.001	4271921	538.694 ug/L
14) T 2-Nitrotoluene	12.388	3671564	510.797 ug/L
15) T 4-Nitrotoluene	13.176	3089805	529.604 ug/L
16) T 3-Nitrotoluene	14.516	3977559	524.069 ug/L
17) T PETN	20.331	4695128	500.691 ug/L m

(f)=RT Delta > 1/2 Window

(m)=manual int.

Data Path : Z:\HPCHEM1\HPLC_L\Data\HL032013\
Data File : HL005153.D
Signal(s) : VWD1A.CH
Acq On : 20-Mar-2013, 14:00
Operator : QM
Sample : HSTDICCC0500
Misc :
ALS Vial : 4 Sample Multiplier: 1

Instrument :
HPLC_L
LabSampleId :
HSTDICCC0500

Integration File: autoint1.e
Quant Time: Mar 22 22:40:34 2013
Quant Method : Z:\HPCHEM1\HPLC_L\METHODS\HL032013.M
Quant Title :
QLast Update : Mon Mar 11 14:17:48 2013
Response via : Initial Calibration
Integrator: ChemStation

Volume Inj. : 50 ul
Signal Phase : ZORBAX Extend-C18
Signal Info : 4.6 X 100mm 3.5-Micron 80A

