


Method : H:\VOL1\HPLC\HPLC L\METHODS\HL102412.M (Chemstation Integrator)
 Title : SW8330A/B Primary (Extend C18 3.5u 4.6x10mm)
 Last Update : Thu Oct 25 10:56:40 2012


Calibration Files

CAL1 =HL003621.D CAL2 =HL003620.D CAL3 =HL003619.D
 CAL4 =HL003618.D CAL5 =HL003617.D CAL6 =HL003615.D

	Compound	CAL1	CAL2	CAL3	CAL4	CAL5	CAL6	Avg	%RSD
1) S	4-Nitroaniline	2.926	5.364	4.420	4.074	3.893	3.878	4.092 E3	19.46
2) T	HMX	2.642	3.670	4.072	4.015	4.075	4.074	3.758 E3	15.14
3) T	RDX	4.728	5.317	5.581	5.568	5.551	5.549	5.382 E3	6.23
4) T	1,3,5-Trinitrobenze	1.262	1.153	1.185	1.173	1.172	1.172	1.186 E4	3.26
5) T	Tetryl	8.815	8.296	8.492	8.452	8.470	8.568	8.516 E3	2.02
6) T	2,4,6-Trinitrotolue	1.238	1.119	1.130	1.118	1.117	1.116	1.139 E4	4.25
7) T	2-Amino-4,6-dinitro	0.970	1.073	1.119	1.113	1.116	1.114	1.084 E4	5.39
8) T	4-Amino-2,6-dinitro	7.041	7.593	7.195	7.121	7.142	7.204	7.216 E3	2.69
9) T	2,4-Dinitrotoluene	1.780	1.625	1.548	1.543	1.543	1.542	1.597 E4	5.97
10) T	2,6-Dinitrotoluene	9.164	8.631	8.456	8.507	8.480	8.475	8.619 E3	3.18
11) T	1,3-Dinitrobenzene	1.694	1.599	1.678	1.664	1.647	1.646	1.655 E4	1.98
12) T	Nitrobenzene	1.094	1.029	1.089	1.079	1.065	1.065	1.070 E4	2.18
13) T	2-Nitrotoluene	6.837	7.158	6.700	6.680	6.619	6.645	6.773 E3	3.00
14) T	3-Nitrotoluene	6.873	7.601	7.238	7.199	7.167	7.158	7.206 E3	3.23
15) T	4-Nitrotoluene	6.213	6.141	5.687	5.633	5.595	5.624	5.816 E3	4.85
16) T	3,5-Dinitroaniline	1.651	1.544	1.578	1.564	1.544	1.540	1.570 E4	2.68
17) T	Nitroglycerin		5.891	5.776	5.708	5.717	5.715	5.762 E1	1.34
18) T	PETN	1.024	0.899	0.912	0.904	0.907	0.909	0.926 E4	5.23



Method Name: H:\VOL1\HPLC\HPLC_L\METHODS\HL102412.M
Calibration Table Last Updated: Thu Oct 25 10:56:40 2012



Method Name: H:\VOL1\HPLC\HPLC_L\METHODS\HL102412.M
Calibration Table Last Updated: Thu Oct 25 10:56:40 2012