

Data File : H:\VOL1\HPLC\HPLC\_L\121025\HL003665.D  
Acq On : 25-Oct-2012, 18:05  
Sample : 8330NG-CCV  
Misc :  
IntFile : events.e

Vial: 3  
Operator: HLM  
Inst : G1314A  
Multiplr: 1.00

Instrument :  
HPLC\_L  
LabSampleId :  
8330NG-CCV

Method : H:\VOL1\HPLC\HPLC\_L\METHODS\HL102412.M (Chemstation Integrator)  
Title : SW8330A/B Primary (Extend C18 3.5u 4.6x10mm)  
Last Update : Thu Oct 25 10:56:40 2012  
Response via : Multiple Level Calibration

Min. RRF : 0.000 Min. Rel. Area : 50% Max. R.T. Dev 0.50min  
Max. RRF Dev : 15% Max. Rel. Area : 150%

|      | Compound      | AvqRF  | CCRF   | %Dev | Area% | Dev(min) |
|------|---------------|--------|--------|------|-------|----------|
| 17 T | Nitroglycerin | 57.616 | 58.184 | -1.0 | 0#    | 0.03     |

Data File : H:\VOL1\HPLC\HPLC L\121025\HL003665.D  
 Acq On : 25-Oct-2012, 18:05  
 Sample : 8330NG-CCV  
 Misc :  
 IntFile : events.e

Vial: 3  
 Operator: HLM  
 Inst : G1314A  
 Multiplr: 1.00

Instrument :  
 HPLC\_L  
 LabSampleId :  
 8330NG-CCV

Method : H:\VOL1\HPLC\HPLC L\METHODS\HL102412.M (Chemstation Integrator)  
 Title : SW8330A/B Primary (Extend C18 3.5u 4.6x10mm)  
 Last Update : Thu Oct 25 10:56:40 2012  
 Response via : Multiple Level Calibration

Min. RRF : 0.000 Min. Rel. Area : 50% Max. R.T. Dev 0.50min  
 Max. RRF Dev : 15% Max. Rel. Area : 150%

|      | Compound                   | AvqRF  | CCRF     | %Dev   | Area% | Dev(min) |
|------|----------------------------|--------|----------|--------|-------|----------|
| 1 S  | 4-Nitroaniline             | 4.092  | 0.000 E3 | 100.0# | 0#    | -1.34#   |
| 2 T  | HMX                        | 3.758  | 0.000 E3 | 100.0# | 0#    | -0.91#   |
| 3 T  | RDX                        | 5.382  | 0.000 E3 | 100.0# | 0#    | -1.73#   |
| 4 T  | 1,3,5-Trinitrobenzene      | 11.863 | 0.000 E3 | 100.0# | 0#    | -2.70#   |
| 5 T  | Tetryl                     | 8.516  | 0.000 E3 | 100.0# | 0#    | -6.02#   |
| 6 T  | 2,4,6-Trinitrotoluene      | 11.395 | 0.000 E3 | 100.0# | 0#    | -6.43#   |
| 7 T  | 2-Amino-4,6-dinitrotoluene | 10.842 | 0.000 E3 | 100.0# | 0#    | -6.86#   |
| 8 T  | 4-Amino-2,6-dinitrotoluene | 7.216  | 0.000 E3 | 100.0# | 0#    | -7.15#   |
| 9 T  | 2,4-Dinitrotoluene         | 15.968 | 0.000 E3 | 100.0# | 0#    | -8.02#   |
| 10 T | 2,6-Dinitrotoluene         | 8.619  | 0.000 E3 | 100.0# | 0#    | -8.37#   |
| 11 T | 1,3-Dinitrobenzene         | 16.546 | 0.000 E3 | 100.0# | 0#    | -3.57#   |
| 12 T | Nitrobenzene               | 10.702 | 0.000 E3 | 100.0# | 0#    | -4.79#   |
| 13 T | 2-Nitrotoluene             | 6.773  | 0.000 E3 | 100.0# | 0#    | -11.54#  |
| 14 T | 3-Nitrotoluene             | 7.206  | 0.000 E3 | 100.0# | 0#    | -13.52#  |
| 15 T | 4-Nitrotoluene             | 5.816  | 0.000 E3 | 100.0# | 0#    | -12.28#  |
| 16 T | 3,5-Dinitroaniline         | 15.701 | 0.000 E3 | 100.0# | 0#    | -3.89#   |
| 18 T | PETN                       | 9.261  | 0.000 E3 | 100.0# | 0#    | -18.37#  |