


Method : H:\VOL1\HPLC\HPLC L\METHODS\HL123112.M (Chemstation Integrator)
 Title : SW8330A/B Primary (Extend C18 3.5u 4.6x10mm)
 Last Update : Wed Jan 02 11:23:27 2013

Calibration Files

CAL1 =HL003864.D CAL2 =HL003863.D CAL3 =HL003862.D
 CAL4 =HL003858.D CAL5 =HL003861.D CAL6 =HL003860.D

	Compound	CAL1	CAL2	CAL3	CAL4	CAL5	CAL6	Avg	%RSD
1) S	4-Nitroaniline	4.180	3.772	4.013	3.769	3.854	3.737	3.888 E3	4.49
2) T	HMX	1.581	0.725	0.458	0.435	0.430	0.400	0.671 E4	68.71
3) T	RDX	4.630	5.498	5.707	5.548	5.607	5.551	5.424 E3	7.29
4) T	1,3,5-Trinitrobenze	1.173	1.264	1.183	1.167	1.168	1.154	1.185 E4	3.37
5) T	Tetryl	1.011	0.902	0.877	0.860	0.851	0.849	0.892 E4	6.94
6) T	2,4,6-Trinitrotolue	1.368	1.179	1.161	1.114	1.111	1.098	1.172 E4	8.63
7) T	2-Amino-4,6-dinitro	1.213	1.136	1.140	1.099	1.094	1.085	1.128 E4	4.21
8) T	4-Amino-2,6-dinitro	6.832	7.604	7.721	7.272	7.319	7.175	7.320 E3	4.33
9) T	2,4-Dinitrotoluene	1.468	1.676	1.599	1.538	1.543	1.526	1.558 E4	4.56
10) T	2,6-Dinitrotoluene	8.175	9.417	9.085	8.413	8.492	8.362	8.657 E3	5.58
11) T	1,3-Dinitrobenzene	1.834	1.706	1.697	1.647	1.648	1.627	1.693 E4	4.46
12) T	Nitrobenzene	1.158	1.079	1.082	1.038	1.044	1.042	1.074 E4	4.25
13) T	2-Nitrotoluene	6.504	6.736	6.660	6.638	6.615	6.584	6.623 E3	1.17
14) T	3-Nitrotoluene	6.628	7.230	7.112	7.504	7.103	7.095	7.112 E3	3.99
15) T	4-Nitrotoluene	4.881	5.637	5.612	5.674	5.584	5.533	5.487 E3	5.48
16) T	3,5-Dinitroaniline	1.629	1.598	1.612	1.573	1.581	1.565	1.593 E4	1.54
17) T	Nitroqlvcerin		5.826	5.501	5.503	5.666	5.550	5.609 E1	2.47
18) T	PETN	9.526	9.213	9.070	9.110	9.086	9.046	9.175 E3	1.98



Method Name: H:\VOL1\HPLC\HPLC_L\METHODS\HL123112.M
Calibration Table Last Updated: Wed Jan 02 11:23:27 2013