

Data Path : Z:\HPCHEM1\HPLC_L\Data\HL053113\
 Data File : HL005354.D
 Signal(s) : VWD1A.CH
 Acq On : 01-Jun-2013, 00:19
 Operator : QM
 Sample : TEST-W-02
 Misc :
 ALS Vial : 14 Sample Multiplier: 1

Instrument :
 HPLC_L
 ClientSampleId :
 TEST-W-02

Integration File: autoint1.e
 Quant Time: Jun 04 06:14:35 2013
 Quant Method : Z:\HPCHEM1\HPLC_L\METHODS\HL053013.M
 Quant Title :
 QLast Update : Fri May 31 16:54:32 2013
 Response via : Initial Calibration
 Integrator: ChemStation

Volume Inj. : 50 ul
 Signal Phase : ZORBAX Extend-C18
 Signal Info : 4.6 X 100mm 3.5-Micron 80A

Compound	R.T.	Response	Conc Units

System Monitoring Compounds			
1) S 4-Nitroaniline	1.307	955288	229.488 ug/L
Target Compounds			
2) T HMX	0.879	823945	189.482 ug/L m
3) T RDX	1.705	1079993	209.146 ug/L
4) T 1,3,5-Trinitrobenzene	2.639	2509016	229.134 ug/L
5) T 1,3-Dinitrobenzene	3.506	3986167	248.392 ng/mL
6) T 3,5-Dinitroaniline	3.857	3400854	250.595 ug/L
7) T Nitrobenzene	4.718	2394621	238.636 ug/L
8) T Tetryl	6.154	1673645	229.956 ug/L
9) T 2,4,6-Trinitrotoluene	6.430	2827331	239.698 ug/L
10) T 2-Amino-4,6-dinitroto...	6.938	2745631	233.851 ug/L
11) T 4-Amino-2,6-dinitroto...	7.296	1837957	292.555 ug/L
12) T 2,4-Dinitrotoluene	8.030	3841492	243.538 ug/L
13) T 2,6-Dinitrotoluene	8.427	2184851	258.939 ug/L
14) T 2-Nitrotoluene	11.612	1593684	239.950 ug/L
15) T 4-Nitrotoluene	12.340	1325896	242.405 ug/L m
16) T 3-Nitrotoluene	13.608	1757059	250.733 ug/L m
17) T PETN	19.001	2251342	252.413 ug/L m

(f)=RT Delta > 1/2 Window

(m)=manual int.

Data Path : Z:\HPCHEM1\HPLC_L\Data\HL053113\
Data File : HL005354.D
Signal(s) : VWD1A.CH
Acq On : 01-Jun-2013, 00:19
Operator : QM
Sample : TEST-W-02
Misc :
ALS Vial : 14 Sample Multiplier: 1

Instrument :
HPLC_L
ClientSampleId :
TEST-W-02

Integration File: autoint1.e
Quant Time: Jun 04 06:14:35 2013
Quant Method : Z:\HPCHEM1\HPLC_L\METHODS\HL053013.M
Quant Title :
QLast Update : Fri May 31 16:54:32 2013
Response via : Initial Calibration
Integrator: ChemStation

Volume Inj. : 50 ul
Signal Phase : ZORBAX Extend-C18
Signal Info : 4.6 X 100mm 3.5-Micron 80A

