

Data Path : Z:\HPCHEM1\HPLC_L\Data\HL060513\
 Data File : HL005408.D
 Signal(s) : VWD1A.CH
 Acq On : 06-Jun-2013, 03:07
 Operator : QM
 Sample : HSTDCCC0500
 Misc :
 ALS Vial : 1 Sample Multiplier: 1

Instrument :
 HPLC_L
 LabSampleId :
 HSTDCCC0500

Integration File: autoint1.e
 Quant Time: Jun 07 04:18:32 2013
 Quant Method : Z:\HPCHEM1\HPLC_L\METHODS\HL060413.M
 Quant Title :
 QLast Update : Thu Jun 06 16:42:50 2013
 Response via : Initial Calibration
 Integrator: ChemStation

Volume Inj. : 50 ul
 Signal Phase : ZORBAX Extend-C18
 Signal Info : 4.6 X 100mm 3.5-Micron 80A

Compound	R.T.	Response	Conc Units

System Monitoring Compounds			
1) S 4-Nitroaniline	1.299	1998707	516.882 ug/L
Target Compounds			
2) T HMX	0.881	1951598	433.587 ug/L
3) T RDX	1.688	2414021	510.266 ug/L
4) T 1,3,5-Trinitrobenzene	2.610	5408217	479.229 ug/L
5) T 1,3-Dinitrobenzene	3.476	7499898	487.946 ng/mL
6) T 3,5-Dinitroaniline	3.822	6439765	476.929 ug/L
7) T Nitrobenzene	4.666	4710170	477.635 ug/L
8) T Tetryl	6.080	3625458	457.430 ug/L
9) T 2,4,6-Trinitrotoluene	6.345	5259362	462.657 ug/L
10) T 2-Amino-4,6-dinitroto...	6.857	5231633	467.710 ug/L
11) T 4-Amino-2,6-dinitroto...	7.210	3126411	513.667 ug/L
12) T 2,4-Dinitrotoluene	7.950	6874102	470.689 ug/L
13) T 2,6-Dinitrotoluene	8.346	3542173	455.488 ug/L
14) T 2-Nitrotoluene	11.439	2966739	446.550 ug/L
15) T 4-Nitrotoluene	12.155	2557501	470.978 ug/L
16) T 3-Nitrotoluene	13.411	3396781	471.595 ug/L
17) T PETN	18.846	5682897	597.793 ug/L

(f)=RT Delta > 1/2 Window

(m)=manual int.

Data Path : Z:\HPCHEM1\HPLC_L\Data\HL060513\
Data File : HL005408.D
Signal(s) : VWD1A.CH
Acq On : 06-Jun-2013, 03:07
Operator : QM
Sample : HSTDCCC0500
Misc :
ALS Vial : 1 Sample Multiplier: 1

Instrument :
HPLC_L
LabSampleId :
HSTDCCC0500

Integration File: autoint1.e
Quant Time: Jun 07 04:18:32 2013
Quant Method : Z:\HPCHEM1\HPLC_L\METHODS\HL060413.M
Quant Title :
QLast Update : Thu Jun 06 16:42:50 2013
Response via : Initial Calibration
Integrator: ChemStation

Volume Inj. : 50 ul
Signal Phase : ZORBAX Extend-C18
Signal Info : 4.6 X 100mm 3.5-Micron 80A

