

Data Path : Z:\HPCHEM1\HPLC_L\Data\HL060613\
 Data File : HL005416.D
 Signal(s) : VWD1A.CH
 Acq On : 06-Jun-2013, 18:57
 Operator : QM
 Sample : HSTDCCC0500
 Misc :
 ALS Vial : 1 Sample Multiplier: 1

Instrument :
 HPLC_L
 LabSampleId :
 HSTDCCC0500

Integration File: autoint1.e
 Quant Time: Jun 07 05:37:19 2013
 Quant Method : Z:\HPCHEM1\HPLC_L\METHODS\HL060413.M
 Quant Title :
 QLast Update : Thu Jun 06 16:42:50 2013
 Response via : Initial Calibration
 Integrator: ChemStation

Volume Inj. : 50 ul
 Signal Phase : ZORBAX Extend-C18
 Signal Info : 4.6 X 100mm 3.5-Micron 80A

Compound	R.T.	Response	Conc Units

System Monitoring Compounds			
1) S 4-Nitroaniline	1.306	1821097	470.951 ug/L
Target Compounds			
2) T HMX	0.885	1883232	418.398 ug/L
3) T RDX	1.698	2404491	508.252 ug/L
4) T 1,3,5-Trinitrobenzene	2.629	5368695	475.727 ug/L
5) T 1,3-Dinitrobenzene	3.503	7503917	488.207 ng/ml
6) T 3,5-Dinitroaniline	3.856	6485216	480.295 ug/L
7) T Nitrobenzene	4.698	4675365	474.105 ug/L
8) T Tetryl	6.134	3503768	442.076 ug/L
9) T 2,4,6-Trinitrotoluene	6.400	5445393	479.022 ug/L
10) T 2-Amino-4,6-dinitroto...	6.921	5106453	456.519 ug/L
11) T 4-Amino-2,6-dinitroto...	7.276	3008726	494.332 ug/L
12) T 2,4-Dinitrotoluene	8.009	6912433	473.314 ug/L
13) T 2,6-Dinitrotoluene	8.405	3649495	469.289 ug/L
14) T 2-Nitrotoluene	11.539	3182895	479.085 ug/L
15) T 4-Nitrotoluene	12.275	2688816	495.161 ug/L
16) T 3-Nitrotoluene	13.547	3456431	479.876 ug/L
17) T PETN	18.979	5590523	588.076 ug/L

(f)=RT Delta > 1/2 Window

(m)=manual int.

Data Path : Z:\HPCHEM1\HPLC_L\Data\HL060613\
Data File : HL005416.D
Signal(s) : VWD1A.CH
Acq On : 06-Jun-2013, 18:57
Operator : QM
Sample : HSTDCCC0500
Misc :
ALS Vial : 1 Sample Multiplier: 1

Instrument :
HPLC_L
LabSampleId :
HSTDCCC0500

Integration File: autoint1.e
Quant Time: Jun 07 05:37:19 2013
Quant Method : Z:\HPCHEM1\HPLC_L\METHODS\HL060413.M
Quant Title :
QLast Update : Thu Jun 06 16:42:50 2013
Response via : Initial Calibration
Integrator: ChemStation

Volume Inj. : 50 ul
Signal Phase : ZORBAX Extend-C18
Signal Info : 4.6 X 100mm 3.5-Micron 80A

