

Data File : H:\VOL1\HPLC\HPLC L\120615\HL003311.D  
Acq On : 16-Jun-2012, 00:12  
Sample : 8330NG-CCV  
Misc :  
IntFile : events.e

Vial: 3  
Operator: HLM  
Inst : G1314A  
Multiplr: 1.00

Instrument :  
HPLC\_L  
ClientSampleId :  
CCAL06

Method : H:\VOL1\HPLC\HPLC L\METHODS\HL060712.M (Chemstation Integrator)  
Title : SW8330A/B Primary (Extend C18 3.5u 4.6x10mm)  
Last Update : Wed Jun 13 15:12:30 2012  
Response via : Multiple Level Calibration

Min. RRF : 0.000 Min. Rel. Area : 50% Max. R.T. Dev 0.50min  
Max. RRF Dev : 20% Max. Rel. Area : 150%

|      | Compound      | AvqRF  | CCRF   | %Dev   | Area% | Dev(min) |
|------|---------------|--------|--------|--------|-------|----------|
| 17 T | Nitroglycerin | 59.494 | 71.804 | -20.7# | 0#    | 0.00     |

Data File : H:\VOL1\HPLC\HPLC L\120615\HL003311.D  
 Acq On : 16-Jun-2012, 00:12  
 Sample : 8330NG-CCV  
 Misc :  
 IntFile : events.e

Vial: 3  
 Operator: HLM  
 Inst : G1314A  
 Multiplr: 1.00

Instrument :  
 HPLC\_L  
 ClientSampleId :  
 CCAL06

Method : H:\VOL1\HPLC\HPLC L\METHODS\HL060712.M (Chemstation Integrator)  
 Title : SW8330A/B Primary (Extend C18 3.5u 4.6x10mm)  
 Last Update : Wed Jun 13 15:12:30 2012  
 Response via : Multiple Level Calibration

Min. RRF : 0.000 Min. Rel. Area : 50% Max. R.T. Dev 0.50min  
 Max. RRF Dev : 20% Max. Rel. Area : 150%

|      | Compound                   | AvqRF  | CCRF     | %Dev   | Area% | Dev(min) |
|------|----------------------------|--------|----------|--------|-------|----------|
| 1 S  | 4-Nitroaniline             | 3.637  | 0.000 E3 | 100.0# | 0#    | -1.38#   |
| 2 T  | HMX                        | 4.633  | 0.000 E3 | 100.0# | 0#    | -0.93#   |
| 3 T  | RDX                        | 5.679  | 0.000 E3 | 100.0# | 0#    | -1.79#   |
| 4 T  | 1,3,5-Trinitrobenzene      | 12.207 | 0.000 E3 | 100.0# | 0#    | -2.80#   |
| 5 T  | Tetryl                     | 8.916  | 0.000 E3 | 100.0# | 0#    | -6.42#   |
| 6 T  | 2,4,6-Trinitrotoluene      | 11.779 | 0.000 E3 | 100.0# | 0#    | -6.77#   |
| 7 T  | 2-Amino-4,6-dinitrotoluene | 11.168 | 0.000 E3 | 100.0# | 0#    | -7.26#   |
| 8 T  | 4-Amino-2,6-dinitrotoluene | 7.547  | 0.000 E3 | 100.0# | 0#    | -7.62#   |
| 9 T  | 2,4-Dinitrotoluene         | 15.838 | 0.000 E3 | 100.0# | 0#    | -8.46#   |
| 10 T | 2,6-Dinitrotoluene         | 8.802  | 0.000 E3 | 100.0# | 0#    | -8.86#   |
| 11 T | 1,3-Dinitrobenzene         | 16.606 | 0.000 E3 | 100.0# | 0#    | -3.71#   |
| 12 T | Nitrobenzene               | 11.201 | 0.000 E3 | 100.0# | 0#    | -4.99#   |
| 13 T | 2-Nitrotoluene             | 6.976  | 0.000 E3 | 100.0# | 0#    | -12.22#  |
| 14 T | 3-Nitrotoluene             | 7.385  | 0.000 E3 | 100.0# | 0#    | -14.32#  |
| 15 T | 4-Nitrotoluene             | 5.946  | 0.000 E3 | 100.0# | 0#    | -13.01#  |
| 16 T | 3,5-Dinitroaniline         | 19.018 | 0.000 E3 | 100.0# | 0#    | -4.07#   |
| 18 T | PETN                       | 9.632  | 0.000 E3 | 100.0# | 0#    | -19.59#  |