

Data Path : Z:\HPCHEM1\HPLC_L\Data\HL061813\
 Data File : HL005427.D
 Signal(s) : VWD1A.CH
 Acq On : 18-Jun-2013, 20:15
 Operator : QM
 Sample : HSTDICC0020
 Misc :
 ALS Vial : 6 Sample Multiplier: 1

Instrument :
 HPLC_L
 LabSampleId :
 HSTDICC0020

Integration File: autoint1.e
 Quant Time: Jun 20 07:28:50 2013
 Quant Method : Z:\HPCHEM1\HPLC_L\METHODS\HL061813.M
 Quant Title :
 QLast Update : Wed Jun 19 14:54:43 2013
 Response via : Initial Calibration
 Integrator: ChemStation

Volume Inj. : 50 ul
 Signal Phase : ZORBAX Extend-C18
 Signal Info : 4.6 X 100mm 3.5-Micron 80A

Compound	R.T.	Response	Conc Units

System Monitoring Compounds			
1) S 4-Nitroaniline	5.642	177638	44.743 ug/L
Target Compounds			
2) T HMX	4.128	159909	35.254 ug/L
3) T RDX	6.784	223479	42.179 ug/L
4) T 1,3,5-Trinitrobenzene	9.368	475830	43.985 ug/L m
5) T 1,3-Dinitrobenzene	12.387	699558	45.314 ng/mlm
6) T 3,5-Dinitroaniline	12.779	617943	45.674 ug/L m
7) T Teteryl	13.410	390811	38.151 ug/L m
8) T Nitrobenzene	14.611	418105	62.128 ug/L
9) T 2,4,6-Trinitrotoluene	15.877	573401	50.931 ug/L
10) T 4-Amino-2,6-dinitroto...	16.506	329048	30.711 ug/L m
11) T 2-Amino-4,6-dinitroto...	17.347	586670	105.207 ug/L m
12) T 2,6-Dinitrotoluene	19.783	240827	15.911 ug/L m
13) T 2,4-Dinitrotoluene	20.338	698447	85.886 ug/L
14) T 2-Nitrotoluene	24.808	306466	43.205 ug/L m
15) T 4-Nitrotoluene	26.657	258964	46.296 ug/L m
16) T 3-Nitrotoluene	28.898	282943	42.677 ug/L m

(f)=RT Delta > 1/2 Window

(m)=manual int.

Data Path : Z:\HPCHEM1\HPLC_L\Data\HL061813\
Data File : HL005427.D
Signal(s) : VWD1A.CH
Acq On : 18-Jun-2013, 20:15
Operator : QM
Sample : HSTDICC0020
Misc :
ALS Vial : 6 Sample Multiplier: 1

Instrument :
HPLC_L
LabSampleId :
HSTDICC0020

Integration File: autoint1.e
Quant Time: Jun 20 07:28:50 2013
Quant Method : Z:\HPCHEM1\HPLC_L\METHODS\HL061813.M
Quant Title :
QLast Update : Wed Jun 19 14:54:43 2013
Response via : Initial Calibration
Integrator: ChemStation

Volume Inj. : 50 ul
Signal Phase : ZORBAX Extend-C18
Signal Info : 4.6 X 100mm 3.5-Micron 80A

