

Method : H:\VOL1\HPLC\HPLC L\METHODS\HL091912.M (Chemstation Integrator)
Title : SW8330A/B Primary (Extend C18 3.5u 4.6x10mm)
Last Update : Thu Sep 20 09:32:22 2012
Response via : Initial Calibration

Instrument :
HPLC_L
ClientSampleId :
8330-MRL

Continuing Calibration File: HL003472.D

Min. RRF : 0.000 Min. Rel. Area : 50%
Max. RRF Dev : 20% Max. Rel. Area : 150%

	Compound	AvgRF	CCRF		%Dev	Area%
1 S	4-Nitroaniline	3.716	3.814	E3	-2.6	138
2 T	HMX	3.912	3.415	E3	12.7	112
3 T	RDX	5.200	4.500	E3	13.5	90
4 T	1,3,5-Trinitrobenzene	12.140	13.563	E3	-11.7	93
5 T	Tetryl	8.535	9.394	E3	-10.1	103
6 T	2,4,6-Trinitrotoluene	11.268	13.190	E3	-17.1	106
7 T	2-Amino-4,6-dinitrotoluene	10.886	9.111	E3	16.3	86
8 T	4-Amino-2,6-dinitrotoluene	6.911	6.378	E3	7.7	100
9 T	2,4-Dinitrotoluene	15.094	14.881	E3	1.4	95
10 T	2,6-Dinitrotoluene	8.704	7.974	E3	8.4	79
11 T	1,3-Dinitrobenzene	15.843	15.640	E3	1.3	96
12 T	Nitrobenzene	10.281	11.698	E3	-13.8	110
13 T	2-Nitrotoluene	6.595	6.457	E3	2.1	92
14 T	3-Nitrotoluene	7.072	8.470	E3	-19.8	111
15 T	4-Nitrotoluene	5.520	6.059	E3	-9.8	106
16 T	3,5-Dinitroaniline	15.294	16.123	E3	-5.4	82
17 T	Nitroglycerin	56.791	0.000		100.0#	0#
18 T	PETN	9.675	9.851	E3	-1.8	89